Robotic-assisted (RPL-4) compared with Video-assisted (VATS) lobectomy for early-stage NSCLC resulted in a greater health benefit that **favored robotic-assisted surgery over VATS surgery in utility scores.**



European Quality of Life 5-Dimension 5-Level Instrument (EQ-5D-5L) Responses to 12 months

Robotic-assisted lobectomy (RPL-4)

Video-assisted lobectomy (VATS-Lobectomy)

Note: * P value of <.05 was considered statistically significant. RPL-4: Robotic-assisted lobectomy VATS: Video-assisted lobectomy. n's are included to show the number of patients who completed EQ-5D-5L at a specific timepoint. HRQOL: Health related quality of life utilities were derived from EQ-5D-5L. Seemingly Unrelated Regression is applied to estimate the effect, adjusting for baseline characteristics (age, sex) and stratification factors (surgeon).

Purpose

To determine if robotic-assisted lobectomy (RPL-4) is costeffective and offers improved patient-reported health utility for patients with early-stage NSCLC (Non-small cell lung cancer) when compared to video-assisted lobectomy (VATS-Lobectomy)

Study design

RAVAL (Robotic-Assisted vs. Video-Assisted Lobectomy) is an international, multi-center, prospective, blinded RCT with a 1:1 allocation ratio to two parallel arms: RPL-4 (intervention) and VATS-Lobectomy (control)

Outcomes measured

Primary Endpoint for Phase A was the difference in Health Utility scores between the treatment groups at 12-weeks. **Secondary Outcomes** that were measured and reported include: differences in Health Utility scores between the treatment groups at 3- and 7-weeks, and 6- and 12-months, short-term clinical and oncological outcomes (total number of LN stations sampled/dissected), and the incremental cost per QALY (Quality Adjusted Life Years) of RPL-4 relative to VATS-Lobectomy

Key results

Primary Outcome: The mean 12-week health utility score was 0.85 (0.10) for RPL-4 and 0.80 (0.19) for VATS-Lobectomy (*p*=0.02)

Secondary Outcomes Results: At the 12-month time horizon, RPL-4 was associated with an additional cost of \$14,925.62/QALY gained (95% CI \$6,843.69, \$23,007.56) relative to VATS-Lobectomy There were no significant differences in HU (Health Utility) scores between treatment groups at 3-weeks. The mean HU score at 7-weeks was 0.84 (0.14) for the RPL-4 arm and 0.78 (0.18) for the VATS-Lobectomy arm (MD 0.04, 95% CI 0.0001, 0.08; p=0.04). The mean HU scores at 6- months and 12-months were no different between cohorts.

Significantly more Lymph nodes were sampled [10 (8-13) vs 8 (5-10); p=0.003] in the RPL-4 arm. Conversion to thoracotomy occurred in 7.41% (6/81) and 15.66% (13/83) of cases in the RPL-4 and VATS-Lobectomy arms,

respectively (p=0.10)

Study information

Robotic-assisted (RPL-4) compared with Video-assisted (VATS) lobectomy for early-stage NSCLC is **cost effective**



Note: **QALY**: Quality-Adjusted Life Year is the metric used to measure how well different medical treatments lengthen and/or improve patients' lives. : ICER: Incremental cost-effectives ratio is the difference in cost between two possible interventions, divided by the difference in effect. ICER of <\$50,000 is considered favorable when introducing new technology.

Incremental Cost Per QALY of RPL-4 Relative to VATS-Lobectomy

Economic outcomes measured

Cost-effectiveness: Incremental cost-effectiveness ratio (ICER) at a threshold of \$50,000 per QALY gained.

Health utility outcome: Quality-adjusted life-years (QALY) using Health related quality of life (HRQOL) utilities derived from EQ-5D-5L.

Resource use and unit cost: Micro-costing¹ **methods using all** direct medical and nonmedical costs from admission to discharge, capital cost and home care-related costs within 30 day.

Results:

Incremental Cost: **\$179.37** (\$154.20, \$204.54) Incremental QALY: **0.0120** (-0.1243, 0.1483) Incremental cost per QALY (ICER): **\$14925.62** (\$6843.685, \$23007.557)*

Note: Seemingly Unrelated Regression is applied to estimate the effect, adjusting for baseline characteristics (age, sex) and stratification factors (surgeon). The regression model for QALY also adjusts for baseline heath utility. Multivariate Imputation by Chained Equations is used to impute the missing data in QALY.

*CI of ICER is generated by 10,000 bootstrap samples using biascorrected and accelerated method.

¹Kaur M, Dogra , Xie F, et al. Robotic Versus Video-Assisted Thoracoscopic Lung Resection During Early Program Development. *Ann Thorac Surg.* 2018;105:1050-7. doi: https://doi.org/10.1016/j.athoracsur.2017.11.13

Study information

Robotic-assisted (RPL-4) compared with Video-assisted (VATS) lobectomy for early-stage NSCLC resulted in statistically significant increase in lymph node stations sampled and statistically significant more lymph nodes examined



Lymph Node Dissection

Total number of Lymph Node Stations, median (IQR) for RRL-4 was 6 (5-7) station compared to VATS 5 (4-6) stations and was shown to be a statistically significant, p-value of .02.

Total number of Lymph Nodes Examined, median (IQR) for RRL-4 was 10 (8-13) lymph nodes compared to VATS 8 (5-10) lymph nodes and was shown to be a statistically significant, p-value of .003.

Study information

INTUITIVE © 2023 Intuitive Surgical Operations, Inc. All rights reserved

Robotic-assisted (RPL-4) compared with Video-assisted (VATS) lobectomy for early-stage NSCLC resulted in **statistically significant reduction** in **estimated blood loss** as well as a **lower conversion rate**



Note: P value of <.05 was considered statistically significant

Short-Term Perioperative Outcomes

Amount of Blood Loss, mL, median (IQR)

Units of blood, mean (SD) for RPL-4 was 50 (30-150) units of blood transfused compared to VATS was 150 (50-200) units of blood transfused was shown to be a statistically significant, p-value of .002.

Converted to Thoracotomy, n (%)

For RPL-4 6 (7.41%) of patients converted to thoracotomy compared to VATS was 13 (15.66%) , p-value of .10.

Total procedure time, minutes median (IQR)

For RPL-4 median procedure time was 143 (118-170) and was shown to be compared to VATS 140 (119-168), p-value of .84.

Total Operating room, minutes median (IQR)

For RPL-4 median operating room (OR) time was 203 (165-234) and was shown to be comparable to VATS 193 (171-225), p-value of .62.



Robotic-assisted (RPL-4) compared with Video-assisted (VATS) lobectomy for early-stage NSCLC resulted in **reduction** in **intraoperative complications**



Note: P value of <.05 was considered statistically significant.

Short-Term Perioperative Outcomes

Length of stay, day, median (IQR)

For RPL-4 median length of stay was 3 (2-5) and was shown to be compared to VATS 3 (2-5). p-value of .84

Intraoperative complication, n (%)

For RPL-4 7 (8.64%) of patients had a intraoperative complication compared to VATS 11 (13.25%), p-value of .35.

Prolonged air leak, n (%)

For RPL-4 14 (17.28%) of patients had a prolonged air leak compared to VATS 11 (13.25%), p-value of .47.



MAT02411 v1 EU 11/2023 5 of 7

Citation

Patel YS, Baste JM, Shargall Y, Waddell TK, Yasufuku K, Machuca TN, Xie F, Thabane L, Hanna WC. Robotic Lobectomy is Cost-Effective and Provides Comparable Health Utility Scores to Video-Assisted Lobectomy: Early Results of the RAVAL Trial. Ann Surg. 2023 Aug 8. doi: 10.1097/SLA.00000000006073. Epub ahead of print. PMID: 37551615

Study design

Type: RAVAL (Robotic-Assisted vs. Video-Assisted Lobectomy) is an international, multicenter, prospective, blinded RCT with a 1:1 allocation ratio to two parallel arms: RPL-4 (intervention) and VATS-Lobectomy (control) **Timeframe:**

Started recruiting patients in 2016, and is expected to run until 2029, with multiple predefined reporting periods for various outcomes.

Patient population: 164 patients were analyzed after final eligibility review ; RPL-4 (Robotic Portal Lobectomy with 4 Arms):n=81 and VATS-Lobectomy:n=83

Eligible patients who signed consent for a minimally invasive lobectomy between 01/2016 and 07/2020 were screened for eligibility criteria. Inclusion criteria were age >18 years; clinical stagel, II, or IIIa NSCLC; and candidate for minimally invasive pulmonary lobectomy, as determined by the operating surgeon. Patients received neoadjuvant induction therapy when clinically indicated.

Exclusion criteria were clinical stage IIIb or IV NSCLC, because it was expected that factors other than surgery, namely advanced disease symptoms and systemic treatment toxicity, would affect their quality of life; or not a candidate for minimally invasive lobectomy, as deemed by the operating surgeon, due to factors such as body habitus, previous chest surgery, anatomy, or medical condition.

Outcomes measured

Primary outcome for Phase A was the difference in Health Utility scores between the treatment groups at 12-weeks.

Secondary outcomes that were measured and reported include: differences in Health Utility scores between the treatment groups at 3- and 7-weeks, and 6- and 12-months, short-term clinical and oncological outcomes (total number of LN stations sampled/dissected), and the incremental cost per QALY (Quality Adjusted Life Years) of RPL-4 relative to VATS-Lobectomy.

Results / conclusions Primary Outcome

The mean 12-week health utility score was 0.85 (0.10) for RPL-4 and 0.80 (0.19) for VATS-Lobectomy (p=0.02).

Key Secondary Outcomes

At the 12-month time horizon, RPL-4 was associated with an additional cost of \$14,925.62/QALY gained (95% CI \$6,843.69, \$23,007.56) relative to VATS-Lobectomy. Costs calculated using micro-costing methods.¹ There were no significant differences in HU (Health Utility) scores between treatment groups at 3-weeks, where both cohorts experienced a reduction in HU, when compared to preoperative baseline. The mean HU score at 7-weeks was 0.84 (0.14) for the RPL-4 arm and 0.78 (0.18) for the VATS-Lobectomy arm (MD 0.04, 95% CI 0.0001, 0.08; p=0.04), with the VATS-Lobectomy cohort continuing to record a reduction in HU when compared to preoperative baseline, but not the RPL-4 cohort. The mean HU scores at 6- months and 12-months were no different between cohorts.

Significantly more lymph nodes were sampled [10 (8-13) vs 8 (5-10); p=0.003] in the RPL-4 arm. Conversion to thoracotomy occurred in 7.41% (6/81) and 15.66% (13/83) of cases in the RPL-4 and VATS-Lobectomy arms, respectively (p=0.10). There was significantly less blood loss in the RPL-4 arm [50 (30-150) mL] than in the VATS-Lobectomy arm

[150 (50-200) mL; p=0.002]. There were no statistically significant differences between the RPL-4 and VATS-Lobectomy arms in patientreported postoperative pain during admission (p=0.88); length of stay in hospital (p=0.85); chest tube duration (p=0.98);operative time (p=0.84);incidence of prolonged air leak (p=0.47); or duration of postoperative intravenous analgesia (p=0.41) or epidural catheter placement (p=0.50). There were also no statistically significant differences between the RPL-4 and VATS-Lobectomy arms in the incidence of intraoperative adverse events (p=0.35)

¹Kaur M, Dogra , Xie F, et al. Robotic Versus Video-Assisted Thoracoscopic Lung Resection During Early Program Development. *Ann Thorac Surg.* 2018;105:1050-7. doi: https://doi.org/10.1016/j.athoracsur.2017.11.13

Study strengths

One of the specific aims of RAVAL was to generate high-quality data that could be used to calculate a reliable ICER (Incremental cost per QALY) for RPL-4, relative to VATS-Lobectomy. Prospective measurement of HU and blinded patients for 12- months postsurgery, allowed for the calculation of reliable HU data that is subject to minimal bias. The resulting calculated ICER for RPL-4, relative to VATS-Lobectomy, was less than \$15,000/QALY gained. RPL-4 is still more expensive than VATS-Lobectomy, but an ICER of <\$50,000 is considered favorable when introducing new technology. The increasing availability and diversity of the robotic platform is likely to drive costs down over time, favorably influencing the ICER.

Study limitations

Health Utility, and the resulting ICER were measured at the short time horizon of 12months. Longer term follow-up, and future data on mortality, may influence these results in either direction. A limitation of the multicentre design is the difference in postoperative care between centres. Another limitation was the high proportion of patients who reported full health at baseline, which makes small changes in HU scores difficult to measure. Due to the multi-national nature of the trial, measuring opioid consumption or delineating de novo opioid-dependence was not done.



Product Information

The Intuitive Surgical Endoscopic Instrument Control Systems (da Vinci X and da Vinci Xi Surgical Systems) are intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments during urologic surgical procedures, general laparoscopic surgical procedures, gynecologic laparoscopic surgical procedures, general thoracoscopic surgical procedures, and trans-oral otolaryngology surgical procedures restricted to benign tumors and malignant tumors classified as T1 and T2, and for benign base of tongue resection procedures. The systems are indicated for adult and pediatric use (except for trans-oral otolaryngology surgical procedures). They are intended to be used by trained physicians in an operating room environment.

The da Vinci X and da Vinci Xi Surgical Systems are class IIb medical devices CE marked (CE 2460) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

Legal Notices

In order to provide benefit and risk information, Intuitive reviews the highest available level of evidence on representative procedures. Intuitive strives to provide a complete, fair and balanced view of the clinical literature. However, our materials should not be seen as a substitute for a comprehensive scientific review. We encourage patients and physicians to review the original publications and all available literature in order to make an informed decision.

Individuals' outcomes may depend on a number of factors, including but not limited to patient characteristics, disease characteristics and/or surgeon experience.

Some products, features or technologies may not be available in all countries. Please contact your local Intuitive representative for product availability in your region. Refer to the product specific User Manual for indications, contraindications, warnings and other product information.

The information contained in this presentation has been checked and compiled with the greatest care. However, no responsibility is taken for its correctness, completeness and topicality. It is the sole responsibility of the recipient to check all information before using it in the individual case.

Privacy Notice

Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2023 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos are trademarks or registered trademarks of Intuitive Surgical or their respective owner.