

Evidence Navigator: Inguinal Hernia Repair

Systematic literature review & meta-analysis
as of March 1, 2024

Purpose

The Evidence Navigator is a slide presentation representing a summary of the meta-analysis of the highest level of evidence available specific to a given procedure and published as of a particular date. It is created by the Global Evidence Management team within Global Access, Value and Economics (GAVE). It includes information that is available in the public domain. It is a systematic review and meta-analysis of the peer-reviewed literature based on a timeframe within which a literature search has been conducted according to a set of concise inclusion and exclusion criteria. The results of the meta-analysis are presented in the form of forest plots summarized for each outcome according to a comparator and surgical approach of interest. The summary results are reflective of a specific period in time and are subject to change with increasing literature. All of the robotic-assisted surgery procedures mentioned within the Evidence Navigator were performed using a da Vinci® surgical system.

Statistical analysis

All summary measures are shown as odds ratios, risk ratios or risk differences when describing binary outcomes, or as weighted mean differences or standardized mean differences when describing continuous outcomes. Weighting is based on the study sample size and variability of the outcome. A random effect model is used if heterogeneity is statistically significant, otherwise a fixed effect model is used. The Mantel Haenszel summary statistic is used for the overall results. The meta-analysis is performed with RevMan 5.4 (Review Manager, Version 5.4. Copenhagen: The Nordic Cochrane Centre, The Cochrane Collaboration, 2014) or R software (R Foundation for Statistical Computing, Vienna, Austria. URL <https://www.R-project.org/>).

Interpretation notes

When the effect size is measured as a standardized mean difference (SMD), or a risk difference (RD), it is not possible to provide a quantitative conclusion. In such cases, a qualitative conclusion is given with reference to its statistical significance. In some instances, studies may contain some overlapping patient populations. A redundancy check is performed in order to minimize this overlap and bias due to over-reporting.

Glossary

RAS	robotic-assisted surgery
Lap	laparoscopic surgery
LOE	level of evidence
HTA	health technology assessment
RCT	randomized controlled trial
OR	odds ratio
MD	mean difference

WMD	weighted mean difference
RD	risk difference
SMD	standardized mean difference
95% CI	95% confidence interval
I²	test statistic for heterogeneity
EBL	estimated blood loss
LOS	length of hospital stay

Evidence Navigator: Inguinal Hernia Repair Summary Slides

Systematic literature review & meta-analysis
as of March 1, 2024



WHAT DOES THE LITERATURE SHOW?

Systematic literature review:

Da Vinci robotic-assisted inguinal hernia repair

Inclusion criteria

Robotic-assisted inguinal hernia repair performed with a da Vinci surgical system

January 1, 2010 – March 1, 2024

Level of Evidence 1b, 2b, 2c, 3b

RCT, large database, prospective and retrospective cohort studies (with $n \geq 20$ in each cohort)

Exclusion criteria

Not in English

Paper reports on a pediatric population

Publication is an HTA that was not published in a peer-reviewed journal

Alternate technique/approach (e.g. single-port)

No stratified analysis by study arm

Inguinal hernia repair data mixed with other procedures

Original research study does not provide quantitative results for outcomes of interest

Original research publication includes redundant patient population and similar conclusions

Concomitant inguinal hernia repair during prostatectomy

40 publications including:



Robotic-assisted patients: **23,888**



Laparoscopic patients: **91,278**



Open patients: **306,727**

Level of evidence



- 1b- RCTs
- 2b - Prospective cohort studies
- 2c - Database studies
- 3b - Retrospective cohort studies



WHAT DOES THE LITERATURE SHOW?

Systematic literature review key points:

Robotic-assisted vs. laparoscopic inguinal hernia repair



Favors robotic-assisted

- ↓ Conversions by **53%**
- ↓ VAS pain scores during hospitalization (24 hours after surgery) by an average of **1.02 points**
- ↓ 1-year hernia recurrence by **49%**
- ↓ ≥2-year hernia recurrence by **51%**



Comparable outcomes

- ≈ Estimated blood loss
- ≈ Blood transfusions
- ≈ Surgical site infection
- ≈ Inpatient length of hospital stay
- ≈ Outpatient length of hospital stay
- ≈ 30-day postoperative complications
- ≈ 30-day readmissions
- ≈ 30-day reoperations
- ≈ 30-day mortality
- ≈ 30-day hernia recurrence
- ≈ 30-day acute postoperative pain
- ≈ 1-year VAS pain score
- ≈ 2-year chronic pain



Favors laparoscopic

- ↓ Operative time for unilateral repair is on average **19.08 min shorter**
- ↓ Operative time for a bilateral repair is on average **21.42 min shorter**

Data collected: March 1, 2024



WHAT DOES THE LITERATURE SHOW?

Systematic literature review key points: Robotic-assisted vs. open inguinal hernia repair



Favors robotic-assisted

- ↓ VAS pain scores during hospitalization (24 hours after surgery) by an average of **3.37 points**
- ↓ 30-day acute postoperative pain by **68%**



Comparable outcomes

- ≈ Blood transfusions
- ≈ Surgical site infection
- ≈ Inpatient length of hospital stay
- ≈ Outpatient length of hospital stay
- ≈ 30-day postoperative complications
- ≈ 30-day readmissions
- ≈ 30-day reoperations
- ≈ 30-day mortality
- ≈ 30-day hernia recurrence
- ≈ 1-year hernia recurrence
- ≈ ≥ 2-year hernia recurrence



Favors open

- ↓ Operative time for unilateral repair is on average **22.96 min shorter**
- ↓ Operative time for a bilateral repair is on average **26.69 min shorter**

Data collected March 1, 2024

Evidence Navigator: Inguinal Hernia Repair Technical Slides

Systematic literature review & meta-analysis
as of March 1, 2024

Inguinal Hernia Repair: Literature search methods as of March 1, 2024

Monthly searches were conducted in PubMed, Scopus and Embase.

All citations were exported into a reference management system. Duplications were removed. Titles, abstracts and keywords were reviewed for literature review inclusion by Global Evidence Management team.

All robotic-assisted right colectomies performed with da Vinci® surgical systems. Publications were identified according to inclusion and exclusion criteria described.

Meta-analysis was performed using R software.


40 publications

23,888 patients who underwent RAS

91,278 patients who underwent laparoscopic surgery

306,727 patients who underwent open surgery

Level of evidence



Criteria phase	Details
Identification phase	All robotics publications (library generated from monthly search process) N=39,985 library size at the time of search March 1, 2024
Inclusion criteria	
1. Robotic-assisted Inguinal hernia repair procedure	Robotic primary inguinal hernia repair N = 317 (excluded N = 39,668)
2. Year ≥ 2010	Articles published ≥ 2010 N = 314 (excluded N = 3)
3. LOE = 1b, 2b, 2c, 3b	Articles with LOE= 1b, 2b, 2c N = 78 (excluded N = 236)
4. Study is an RCT, prospective or retrospective study or large database study with comparative cohorts (robotic-assisted vs lap and/or open surgery) and sample size N≥20	Comparator cohorts N = 75 (excluded N = 3)
Exclusion criteria	N = 35 excluded publications:
1. Not in English	N = 0 (EC#1)
2. Paper reports on a pediatric population	N = 1 (EC#2)
3. Publication is an HTA that was not published in a peer-reviewed journal	N = 0 (EC#3)
4. Alternate technique/approach (e.g., single-port)	N = 0 (EC#4)
5. No stratified analysis by study arm (e.g., combines results from robotic, laparoscopic, and/or open cohorts)	N = 13 (EC#5)
6. Inguinal hernia repair data mixed with another procedure/s	N = 7 (EC#6)
7. Original research study does not provide quantitative results for the outcomes of interest	N = 12 (EC#7)
8. Original research publication includes redundant patient population and similar conclusions	N = 0 (EC#8)
9. Inguinal hernia repair after a prostatectomy	N = 2 (EC#9)
Inguinal Hernia Repair Publications: N = 40	


Robotic-assisted vs. laparoscopic inguinal hernia repair

Summary as of March 1, 2024

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic inguinal hernia repair, the evidence for **robotic-assisted inguinal hernia repair using the da Vinci surgical system** demonstrates:

- Significantly lower pain scores (VAS) during hospitalization 24-hour by an average of 1.02 points
- Comparable estimated blood loss
- Comparable outpatient length of hospital stay
- Comparable pain scores (VAS) within 1-year of surgery
- Comparable inpatient length of hospital stay
- Significantly longer operative time for a unilateral repair by an average of 19.08 min
- Significantly longer operative time for a bilateral repair by an average of 21.42 min



Outcome	Robotic-assisted, n	Laparoscopic, n	Effect Size WMD, 95%CI	P-value
Inguinal Hernia Repair Continuous Variables (to March 1, 2024)				
Pain score (VAS) during hospitalization 24-hour ^{2,8,39}				
Subtotal	132	117	-1.02 [-1.85 ; -0.19]	p=0.02
Random, Heterogeneity: p<0.01; I²=80.78%				
EBL, mL ^{4,5,21,28}				
Subtotal	788	974	-1.05 [-3.15 ; 1.05]	p=0.33
Random, Heterogeneity: p<0.01; I² = 81.13%				
Outpatient length of stay, days ^{1,22}				
Subtotal	705	1548	-0.33 [-1.55 ; 0.88]	p=0.59
Random, Heterogeneity: p=0.01; I² = 83.24%				
Pain score (VAS) during hospitalization 1-year ^{4,25}				
Subtotal	43	70	0.09 [-1.33 ; 1.53]	p=0.89
Fixed, Heterogeneity: p=0.19; I²= 41.87%				
Inpatient length of stay, days ^{2,4,5,6,8,9,12,14,15,17,22,28,29,30,31,36,39,40}				
Subtotal	9675	58927	0.14 [-0.17 ; 0.46]	p=0.37
Random, Heterogeneity: p<0.01; I² = 97.58%				
Operative time for unilateral repair, min ^{1,4,7,9,11,13,17,20,22,26,27,28,31,35}				
Subtotal	3123	10488	19.08 [10.84 ; 27.31]	p<0.01
Random, Heterogeneity: p<0.01; I² = 95.52%				
Operative time for bilateral repair, min ^{4,9,20,22,27,29}				
Subtotal	403	533	21.42 [2.14 ; 40.69]	p=0.03
Random, Heterogeneity: p<0.01; I² = 92.90%				


Robotic-assisted vs. laparoscopic inguinal hernia repair

Summary as of March 1, 2024

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic inguinal hernia repair, the evidence for **robotic-assisted inguinal hernia repair using the da Vinci surgical system** demonstrates:

- 53% less likely to have a conversion to open surgery
- 51% less likely to have hernia recurrence ≥2-years of surgery
- 49% less likely to have hernia recurrence within 1-year of surgery
- Comparable acute postoperative pain within 30-days of surgery
- Comparable readmissions rate within 30-days of surgery
- Comparable postoperative complications rate within 30-days of surgery
- Comparable surgical site infection rate
- Comparable reoperations rate within 30-days of surgery




Robotic-assisted vs. laparoscopic inguinal hernia repair

Summary as of March 1, 2024

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic inguinal hernia repair, the evidence for **robotic-assisted inguinal hernia repair using the da Vinci surgical system** demonstrates:

- Comparable blood transfusion rate
- Comparable hernia recurrence within 30-days of surgery
- Comparable mortality rate within 30-days of surgery
- Comparable chronic pain within 2-years of surgery



Outcome	Robotic-assisted, n	Laparoscopic, n	Effect size RD 95% CI	P-value
Inguinal Hernia Repair Binary Variables (to March 1, 2024)				
Blood transfusions, n ^{6,7}				
Subtotal	348	1004	-0.0009 [-0.0071 ; 0.0054]	p=0.78
Fixed, Heterogeneity: p=0.66; I ² = 0%				
Hernia recurrence 30-day, n ^{7,8,23,32}				
Subtotal	1737	3633	0.0004 [-0.0027 ; 0.0036]	p=0.79
Fixed, Heterogeneity: p=0.88; I ² = 0%				
Mortality 30-days, n ^{2,6,7,9,11,12,30}				
Subtotal	7157	19759	0.0005 [-0.0006 ; 0.0015]	p=0.37
Fixed, Heterogeneity: p=1.00; I ² = 0%				
Chronic pain 2-years, n ^{2,28}				
Subtotal	123	123	0.0081 [-0.0231 ; 0.0392]	p=0.61
Fixed, Heterogeneity p=0.45; I ² = 0%				


Robotic-assisted vs. open inguinal hernia repair

Summary as of March 1, 2024

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open inguinal hernia repair, the evidence for **robotic-assisted inguinal hernia repair using the da Vinci surgical system** demonstrates:

- Significantly lower pain scores (VAS) during hospitalization by an average of **3.37 points**
- Comparable inpatient length of hospital stay
- Comparable outpatient length of hospital stay
- Significantly longer operative time for a unilateral repair by an average of **22.96 min**
- Significantly longer operative time for a bilateral repair by an average of **26.69 min**



Outcome	Robotic-assisted, n	Open, n	Effect Size WMD, 95%CI	P-value
Inguinal Hernia Repair Continuous Variables (to March 1, 2024)				
Pain score (VAS) during hospitalization 24-hour ^{16,19}				
Subtotal	124	197	-3.37 [-6.30 ; -0.44]	p=0.02
Random, Heterogeneity: p=0.07; I ² = 68.62%				
Inpatient length of stay, days ^{4,10,12,14,15,18,22,30,36}				
Subtotal	8831	225115	-0.63 [-1.65 ; 0.39]	p=0.23
Random, Heterogeneity: p<0.01; I ² = 99.24%				
Outpatient length of stay, days ^{10,18,22}				
Subtotal	607	609	0.79 [-0.45 ; 2.03]	p=0.21
Random, Heterogeneity: p<0.01; I ² = 94.20%				
Operative time for unilateral repair, min ^{4,7,13,17,19,22,35}				
Subtotal	1969	37721	22.96 [9.14 ; 36.78]	p<0.01
Random, Heterogeneity: p<0.01; I ² = 99.34%				
Operative time for bilateral repair, min ^{4,16,19,22,35}				
Subtotal	159	287	26.69 [10.04 ; 43.34]	p<0.01
Random, Heterogeneity: p<0.01; I ² = 70.01%				


Robotic-assisted vs. open inguinal hernia repair

Summary as of March 1, 2024

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open inguinal hernia repair, the evidence for **robotic-assisted inguinal hernia repair using the da Vinci surgical system** demonstrates:

- 68% less likely to experience early postoperative acute pain acute within 30-days of surgery
- Comparable hernia recurrence within 1-year of surgery
- Comparable readmissions within 30-days of surgery
- Comparable reoperations within 30-days of surgery
- Comparable mortality within 30-days of surgery
- Comparable surgical site infection rate
- Comparable hernia recurrence in 2-years and beyond post-surgery
- Comparable postoperative complications rate within 30-days of surgery




Robotic-assisted vs. open inguinal hernia repair

Summary as of March 1, 2024

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open inguinal hernia repair, the evidence for **robotic-assisted inguinal hernia repair using the da Vinci surgical system** demonstrates:

- Comparable blood transfusion rate
- Comparable hernia recurrence within 30-days of surgery



Outcome	Robotic-assisted, n	Open, n	Effect size RD 95% CI	P-value
Inguinal Hernia Repair Continuous Variables (to March 1, 2024)				
Blood transfusions, n ^{7,10,18}				
Subtotal	608	728	0.0000 [-0.0055 ; 0.0055]	p=1.00
Fixed, Heterogeneity: p=1.00; I ² = 0%				
Hernia recurrence 30-days, n ^{7,23}				
Subtotal	1667	1789	0.0018 [-0.0011 ; 0.0046]	p=0.23
Fixed, Heterogeneity: p=0.85; I ² = 0%				

Appendix

Weighted estimates & bibliography of included studies

Robotic-assisted vs. laparoscopic Inguinal Hernia Repair

Weighted estimates based on 37 studies
Meta-analysis covering period January 1, 2010 – March 1, 2024

This study analyzed continuous variables using weighted means and categorical variables using weighted rates with fixed or random effects models. This method gives more influence to studies with higher weights, providing a more accurate estimate of central tendency when combining results from multiple studies.

Outcomes that favor RAS	
Conversions	2.1% vs 2.8%
≥2-year hernia recurrence	1.1% vs 2.0%
1-year hernia recurrence	0.7% vs 2.3%
24-hour VAS pain score	2.4 score vs 3.4 score

Comparable outcomes	
Estimated blood loss	4.7 ml vs 5.8 ml
Blood transfusions	0% vs 0.1%
Surgical site infections	2.1% vs 2.0%
Inpatient Length of hospital stay	1.6 days vs 1.5 days
Outpatient Length of hospital stay	6.1 hours vs 6.5 hours
30-day postoperative complications	10.9% vs 7.7%
30-day readmissions	2.6% vs 1.6%
30-day reoperations	2.2% vs 0.9%
30-day mortality	0.2% vs 0.1%
30-day hernia recurrence	0.2% vs 0.2%
30-day acute postoperative pain	2% vs 6.2%
1-year VAS pain score	1.9 score vs 1.8 score
2-year chronic pain	1.6% vs 0.8%

Outcomes that favor Laparoscopic	
Operative time Unilateral	80.6 min vs 61.5 min
Operative time Bilateral	93 min vs 71.6 min

Disclaimer: The number of studies used to calculate the weighted estimates for each outcome varies

Robotic-assisted vs. open Inguinal Hernia Repair

Weighted estimates based on 19 studies

Meta-analysis covering period January 1, 2010 – March 1, 2024

This study analyzed continuous variables using weighted means and categorical variables using weighted rates with fixed or random effects models. This method gives more influence to studies with higher weights, providing a more accurate estimate of central tendency when combining results from multiple studies.

Outcomes that favor RAS

30-day Acute postoperative pain	1.4% vs 4.3%
24-hour VAS pain score	0 score vs 3.4 score

Comparable outcomes

Blood Transfusions	0% vs 0%
Surgical Site Infections	0.5% vs 0.7%
Inpatient Length of stay	1.7 days vs 2.4 days
Outpatient Length of stay	6.8 hours vs 6.0 hours
30-day postoperative complications	8.7% vs 8.5%
30-day readmissions	2.5% vs 3.3%
30-day reoperations	1.2% vs 1.1%
30-day mortality	0.1% vs 0.1%
30-day hernia recurrence	0.2% vs 0.1%
1-year hernia recurrence	1.4% vs 3.4%
≥2-year hernia recurrence	2.0% vs 1.6%

Outcomes that favor Open

Operative time Unilateral	84.4 min vs 61.4 min
Operative time Bilateral	115.2 min vs 88.5 min

Disclaimer: The number of studies used to calculate the weighted estimates for each outcome varies

Inguinal hernia repair: bibliography (1/2)

March 1, 2024

1. Abdelmoaty, W. F.D., C. M.Neighom, C.Swanstrom, L. L.Hammill, C. W. (2018). Robotic-assisted versus laparoscopic unilateral inguinal hernia repair: a comprehensive cost analysis. [3b]. Surgical Endoscopy. doi:10.1007/s00464-018-06606-9
2. Aghayeva, A.B., C.Bilgin, I. A.Bengur, F. B.Bas, M.Kirbiyik, E.Aytac, E.Baca, B. (2020). Laparoscopic totally extraperitoneal vs robotic transabdominal preperitoneal inguinal hernia repair: Assessment of short-and long-term outcomes. [3b]. International Journal of Medical Robotics and Computer Assisted Surgery, e2111. doi:10.1002/rcs.2111
3. AlMarzooqi, R.T., S.Huang, L. C.Prabhu, A.Rosen, M. (2019). Review of inguinal hernia repair techniques within the Americas Hernia Society Quality Collaborative. [2c]. Hernia. doi:10.1007/s10029-019-01968-y
4. Amundson, J. R.A., M.Forester, B.Francis, S.Kuchta, K.Denham, W.Linn, J.Ujiki, M. B.Haggerty, S. P. (2022). Laparoscopic and robotic inguinal hernia repair are safe and effective after prior pelvic or low abdominal surgery. [3b]. American Journal of Surgery. doi:10.1016/j.amjsurg.2022.08.011
5. Ayuso, S. A.M., M. N.Katzen, M. M.Aladegbami, B. G.Augenstein, V. A. (2022). Laparoscopic versus robotic inguinal hernia repair: a single-center case-matched study. [3b]. Surg Endosc. doi:10.1007/s00464-022-09368-7
6. Chao, T. C.T., H. Y.Tsai, C. H.Pen, C. M.Wu, C. C.Liao, C. H.Ou, Y. C.Tsai, C. C.Yang, S. D.Tsai, Y. C. (2023). Laparoscopic versus robotic TAPP/TEP inguinal hernia repair: a multicenter, propensity score weighted study. [3b]. Hernia. doi:10.1007/s10029-023-02916-7
7. Charles, E. J.M., J. H.Tache-Leon, C. A.Hallowell, P. T.Sawyer, R. G.Yang, Z. (2017). Inguinal hernia repair: is there a benefit to using the robot? [2c]. Surgical Endoscopy and Other Interventional Techniques, 1-6. doi:10.1007/s00464-017-5911-4
8. Choi, Y. S.K., K. D.Choi, M. S.Heo, Y. S.Yi, J. W.Choe, Y. M. (2023). Initial Experience of Robot-Assisted Transabdominal Preperitoneal (TAPP) Inguinal Hernia Repair by a Single Surgeon in South Korea. [3b]. Medicina (Kaunas), 59(3). doi:10.3390/medicina59030582
9. Ephraim, K.H., B.Mohammad, A.Dan, A.Yehonatan, N.Lior, S.Dina, O.David, H. (2022). Learning curve of robotic inguinal hernia repair in the hands of an experienced laparoscopic surgeon: a comparative study. [3b]. J Robot Surg. doi:10.1007/s11701-021-01362-w
10. Gamagami, R.D., E.Gonzalez, A.D'Amico, L.Richardson, C.Rabaza, J.Kolachalam, R. (2018). Open versus robotic-assisted transabdominal preperitoneal (R-TAPP) inguinal hernia repair: a multicenter matched analysis of clinical outcomes. [3b]. Hernia, 1-10. doi:10.1007/s10029-018-1769-1
11. Gerdes, S.B., R.Liesch, G.Freitag, B.Serra, M.Vonlanthen, R.Bueter, M.Thalheimer, A. (2022). Results of robotic TAPP and conventional laparoscopic TAPP in an outpatient setting: a cohort study in Switzerland. [3b]. Langenbecks Arch Surg. doi:10.1007/s00423-022-02552-2
12. Holleran, T. J.N., M. A.Sparks, A. D.Duncan, J. E.Garrett, M.Brody, F. J. (2021). Trends and outcomes of open, laparoscopic, and robotic inguinal hernia repair in the veterans affairs system. [2c]. Hernia. doi:10.1007/s10029-021-02419-3
13. Hsu, J. L.K., J. K.Zambito, G. M.Korakavi, N.Phillips, A.Phillips, M.Scheeres, D. E.Banks-Venegoni, A. L. (2022). Comparison of Robotic, Laparoscopic, and Open Unilateral Repair of Non-recurrent Inguinal Hernia. [3b]. American Surgeon. doi:10.1177/00031348221136572
14. Huerta, S.T., C.Argo, M.Favela, J.Pharm, T.Kukreja, S.Yan, J.Zhu, H. (2019). Open, Laparoscopic, and Robotic Inguinal Hernia Repair: Outcomes and Predictors of Complications. [3b]. Journal of Surgical Research, 241, 119-127. doi:10.1016/j.jss.2019.03.046
15. Janjua, H.C.P., E.Barry, T. M.Kuo, M. C.Baker, M. S.Kuo, P. C. (2019). The paradox of the robotic approach to inguinal hernia repair in the inpatient setting. [2c]. American Journal of Surgery. doi:10.1016/j.amjsurg.2019.09.012
16. Kakiashvili, E.B., M.Abu Shakra, I.Ganam, S.Bickel, A.Merei, F.Drobot, A.Bogouslavski, G.Kassis, W.Khatib, K.Badran, M.Kluger, Y.Almog, R. (2020). Robotic inguinal hernia repair: Is it a new era in the management of inguinal hernias? [3b]. Asian Journal of Surgery. doi:10.1016/j.asjsur.2020.03.015
17. Khoraki, J.G., P. P.Mazzini, G. S.Pessoa, B. M.Browning, M. G.Aquilina, G. R.Salluzzo, J. L.Wolfe, L. G.Campos, G. M. (2019). Perioperative outcomes and cost of robotic-assisted versus laparoscopic inguinal hernia repair. [3b]. Surgical endoscopy. doi:10.1007/s00464-019-07128-8
18. Kolachalam, R.D., E.D'Amico, L.Richardson, C.Rabaza, J.Gamagami, R.Gonzalez, A. (2017). Early outcomes of robotic-assisted inguinal hernia repair in obese patients: a multi-institutional, retrospective study. [3b]. Surgical Endoscopy and Other Interventional Techniques, 1-7. doi:10.1007/s00464-017-5665-z
19. Kosturakis, A. K.L., K. E.Carroll, N. D.Nicholl, M. B. (2018). First 100 consecutive robotic inguinal hernia repairs at a Veterans Affairs hospital. [3b]. J Robot Surg. doi:10.1007/s11701-018-0812-7
20. Kudsí, O. Y.M., J. C.Paluvui, N.Mabardy, A. S. (2017). Transition from Laparoscopic Totally Extraperitoneal Inguinal Hernia Repair to Robotic Transabdominal Preperitoneal Inguinal Hernia Repair: A Retrospective Review of a Single Surgeon's Experience. [3b]. World Journal of Surgery, 1-7. doi:10.1007/s00268-017-3998-3

Inguinal hernia repair: bibliography (2/2)

March 1, 2024

21. Kudsi, O. Y.B.-A., N.Kaoukabani, G.Gokcal, F. (2022). Comparison of perioperative and mid-term outcomes between laparoscopic and robotic inguinal hernia repair. [3b]. *Surgical Endoscopy*. doi:10.1007/s00464-022-09433-1
22. LeBlanc, K.D., E.Gonzalez, A.Gamagami, R.Pierce, R.Balentine, C.Voeller, G.The Prospective Hernia Study, Group. (2020). Prospective, multicenter, pairwise analysis of robotic-assisted inguinal hernia repair with open and laparoscopic inguinal hernia repair: early results from the Prospective Hernia Study. [2b from 2c]. *Hernia*. doi:10.1007/s10029-020-02224-4
23. Lima, D. L.N., R.Dominguez Profeta, R.Huang, L. C.Cavazzola, L. T.Malcher, F.Sreeramaju, P. (2024). Current trends and outcomes for unilateral groin hernia repairs in the United States using the Abdominal Core Health Quality Collaborative database: A multicenter propensity score matching analysis of 30-day and 1-year outcomes. [2c]. *Surgery*. doi:10.1016/j.surg.2023.11.033
24. Mehaffey, J. H.M., A. D.Mullen, M. G.Yount, K. W.Meneveau, M. O.Smith, P. W.Friel, C. M.Schirmer, B. D. (2017). Adoption of robotics in a general surgery residency program: at what cost? [3b]. *Journal of Surgical Research*, 213, 269-273. doi:10.1016/j.jss.2017.02.052
25. Miller, B. T.P., A. S.Petro, C. C.Beffa, L. R. A.Carbonell, A. M.Hope, W.Warren, J.Higgins, R. M.Jacob, B.Blatnik, J.Krpata, D. M.Tu, C.Costanzo, A.Rosen, M. J. (2022). Laparoscopic versus robotic inguinal hernia repair: 1- and 2-year outcomes from the RIVAL trial. [1b]. *Surgical Endoscopy*. doi:10.1007/s00464-022-09320-9
26. Muysoms, F.V. C., S.Kyle-Leinhase, I.Ballecer, C.Ramaswamy, A. (2018). Robotic-assisted laparoscopic groin hernia repair: observational case-control study on the operative time during the learning curve. [3b]. *Surg Endosc*. doi:10.1007/s00464-018-6236-7
27. Muysoms, F.V., M.Nachtergaele, F.Van Garsse, S.Pletinckx, P.Ramaswamy, A. (2021). Economic assessment of starting robot-assisted laparoscopic inguinal hernia repair in a single-centre retrospective comparative study: the EASTER study. [3b]. *BJS open*, 5(1). doi:10.1093/bjsopen/zraa046
28. Okamoto, N.M., S.Mishima, K.Fujiyama, Y.Wakabayashi, T.Fujita, S.Sakamoto, J.Wakabayashi, G. (2023). Comparison of short-term outcomes of robotic and laparoscopic transabdominal peritoneal repair for unilateral inguinal hernia: a propensity-score matched analysis. [3b]. *Hernia*. doi:10.1007/s10029-022-02730-7
29. Peltrini, R.C., F.Pacella, D.Castiglioni, S.Lionetti, R.Andreuccetti, J.Pignata, G.De Nisco, C.Ferraro, L.Salaj, A.Formisano, G.Bianchi, P. P.Bracale, U. (2022). Robotic versus laparoscopic transabdominal preperitoneal (TAPP) approaches to bilateral hernia repair: a multicenter retrospective study using propensity score matching analysis. [3b]. *Surgical Endoscopy*. doi:10.1007/s00464-022-09614-y
30. Pokala, B.A., P. R.Flores, L.Hennings, D.Oleynikov, D. (2019). Minimally invasive inguinal hernia repair is superior to open: a national database review. [2c]. *Hernia*. doi:10.1007/s10029-019-01934-8
31. Prabhu, A. S.C., A.Hope, W.Warren, J.Higgins, R.Jacob, B.Blatnik, J.Haskins, I.Alkhatib, H.Tastaldi, L.Fafaj, A.Tu, C.Rosen, M. J. (2020). Robotic Inguinal vs Transabdominal Laparoscopic Inguinal Hernia Repair: The RIVAL Randomized Clinical Trial. [1b]. *JAMA Surgery*. doi:10.1001/jamasurg.2020.0034
32. Saito, T.F., Y.Uchino, T.Kurahashi, S.Matsumura, T.Osawa, T.Arikawa, T.Komatsu, S.Kaneko, K.Sano, T. (2020). Preliminary results of robotic inguinal hernia repair following its introduction in a single-center trial. [3b]. *Ann Gastroenterol Surg*, 4(4), 441-447. doi:10.1002/ags3.12341
33. Shah, P. C.d. G., A.Cerfolio, R.Huang, W. C.Huang, K.Song, C.Li, Y.Kreaden, U.Oh, D. S. (2022). Impact of type of minimally invasive approach on open conversions across ten common procedures in different specialties. [2c]. *Surg Endosc*. doi:10.1007/s00464-022-09073-5
34. Sheldon, R. R.D., W. S.Weiss, J. B.Forte, D. M.Sohn, V. Y. (2019). Sage wisdom or anecdotal dictum? Equivalent opioid use after open, laparoscopic, and robotic inguinal hernia repair. [3b]. *American Journal of Surgery*. doi:10.1016/j.amjsurg.2019.02.022
35. Shenoy, R.M., M. A.Jacob, R. L.Kondo, K. K.DeVirgilio, M.Ward, R.Kansagara, D.Shekelle, P. G.Maggard-Gibbons, M.Girgis, M. D.Hynes, D. M. (2022). Robot-Assisted General Surgery Procedures at the Veterans Health Administration: A Comparison of Surgical Techniques. [2c]. *Journal of Surgical Research*, 279, 330-337. doi:10.1016/j.jss.2022.06.032
36. Tatarian, T.N., L.McPartland, C.Brown, A. M.Yang, J.Altieri, M. S.Spaniolas, K.Docimo, S.Pryor, A. D. (2021). Comparative perioperative and 5-year outcomes of robotic and laparoscopic or open inguinal hernia repair: a study of 153,727 patients in the state of New York. [2c]. *Surgical Endoscopy*. doi:10.1007/s00464-020-08211-1
37. Tonelli, C. M.L., I.Bunn, C.Kulshrestha, S.Abdelsattar, Z. M.Cohn, T.Luchette, F. A.Baker, M. S. (2022). Contemporary Matched-Cohort Comparison of Surgical Approach to Inguinal Hernia Repair: Are Minimally Invasive Approaches Associated with Higher Rates of Recurrence? [2c]. *J Am Coll Surg*, 235(1), 119-127. doi:10.1097/XCS.0000000000000235
38. Varvoglis, D. N.S.-C., M.Olson, M. A.DeAngelo, N.Garbarine, I.Lipman, J.Farrell, T. M.Overby, D. W.Perez, A.Zhou, R. (2022). Comparison of post-operative outcomes of large direct inguinal hernia repairs based on operative approach (open vs. laparoscopic vs. robotic) using the ACHQC (Abdominal Core Health Quality Collaborative) database. [2c]. *Surgical Endoscopy*. doi:10.1007/s00464-022-09805-7
39. Waite, K. E.H., M. A.Doyle, P. J. (2016). Comparison of robotic versus laparoscopic transabdominal preperitoneal (TAPP) inguinal hernia repair. [3b]. *J Robot Surg*. doi:10.1007/s11701-016-0580-1
40. Zayan, N. E.M., M. P.Schwartz, J. S.Narula, V. K. (2019). A direct comparison of robotic and laparoscopic hernia repair: patient-reported outcomes and cost analysis. [3b]. *Hernia*. doi:10.1007/s10029-019-01943-7

Disclosures

Important Safety Information

(US) Serious complications may occur in any surgery, including da Vinci surgery, up to and including death. Serious risks include, but are not limited to, injury to tissues and organs and conversion to other surgical techniques which could result in a longer operative time and/or increased complications. For summary of the risks associated with surgery refer to www.intuitive.com/safety.

Da Vinci Xi®/da Vinci X® system precaution statement

The demonstration of safety and effectiveness for the representative specific procedures did not include evaluation of outcomes related to the treatment of cancer (overall survival, disease-free survival, local recurrence), except for radical prostatectomy which was evaluated for overall survival, or treatment of the patient's underlying disease/condition. Device usage in all surgical procedures should be guided by the clinical judgment of an adequately trained surgeon.

(EU) Medical devices, CE 2460, refer to Instructions For Use for further information.

For product intended use and/or indications for use, risks, cautions, and warnings and full prescribing information, refer to the associated user manual(s) or visit <https://manuals.intuitivesurgical.com/market>.

Some products, features or technologies may not be available in all countries. Please contact your local Intuitive representative for product availability in your region.

Individual outcomes may depend on a number of factors, including but not limited to patient characteristics, disease characteristics, and/or surgeon experience.

Privacy Notice: Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2025 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos are trademarks or registered trademarks of Intuitive Surgical or their respective owner.

INTUITIVE

intuitive.com