

Evidence Navigator: Lobectomy

Systematic literature review & meta-analysis
as of December 31, 2022

Purpose

The Evidence Navigator is a slide presentation representing a summary of the meta-analysis of the highest level of evidence available specific to a given procedure and published as of a particular date. It is created by the Global Evidence Management team within Global Access, Value and Economics (GAVE). It includes information that is available in the public domain. It is a systematic review and meta-analysis of the peer-reviewed literature based on a timeframe within which a literature search has been conducted according to a set of concise inclusion and exclusion criteria. The results of the meta-analysis are presented in the form of forest plots summarized for each outcome according to a comparator and surgical approach of interest. The summary results are reflective of a specific period in time and are subject to change with increasing literature. All of the robotic-assisted surgery procedures mentioned within the Evidence Navigator were performed using a da Vinci® surgical system.

Statistical analysis

All summary measures are shown as odds ratios, risk ratios or risk differences when describing binary outcomes, or as standardized mean differences or weighted mean differences when describing continuous outcomes. Weighting is based on the study sample size and variability of the outcome. A fixed effect model is used if heterogeneity was not statistically significant or not applicable, and a random effect model is used if heterogeneity was statistically significant. Mantel Haenszel summary statistic is used for overall results. Meta-analysis is performed with RevMan 5.4 (Review Manager, Version 5.4. Copenhagen: The Nordic Cochrane Centre, The Cochrane Collaboration, 2014) or R software (R Foundation for Statistical Computing, Vienna, Austria. URL <https://www.R-project.org/>).

Interpretation notes

When the effect size is measured as a standardized mean difference (SMD), or a risk difference (RD), it is not possible to provide a quantitative conclusion. In such cases, a qualitative conclusion is given with reference to its statistical significance. In some instances, studies may contain some overlapping patient populations. A redundancy check is performed in order to minimize this overlap and bias due to over-reporting.


Glossary

RAS	robotic-assisted surgery
VATS	video-assisted thoracoscopic surgery
LOE	level of evidence
HTA	health technology assessment
RCT	randomized controlled trial
OR	odds ratio
MD	mean difference
WMD	weighted mean difference
RD	risk difference

SMD	standardized mean difference
95% CI	95% confidence interval
I²	test statistic for heterogeneity
EBL	estimated blood loss
LOS	length of hospital stay
PSM	positive surgical margins
LNy	lymph node yield
ICU	intensive care unit

Evidence Navigator: Lobectomy Summary Slides

Systematic literature review & meta-analysis
as of December 31, 2022



WHAT DOES THE LITERATURE SHOW?

Systematic literature review & meta-analysis methods: Da Vinci robotic-assisted lobectomy

Inclusion criteria

Robotic-assisted lobectomy performed with a da Vinci® surgical system

January 1, 2010 – December 31, 2022

Level of Evidence = 1b, 2b, 2c

RCT, large database, and prospective cohort studies (with $n \geq 20$ in each cohort)

Exclusion criteria

Not in English

Paper reports on a pediatric population

Publication is an HTA that was not published in a peer-reviewed journal

Alternate technique/approach (e.g. single-port)

No stratified analysis by study arm

Lobectomy data mixed with lung mediastinal resection (e.g., data from multiple surgical procedures combined)

Original research study does not provide quantitative results for outcomes of interest

Original research publication includes redundant patient population and similar conclusions

48 publications including



RAS patients: **175,043**



VATS patients: **577,880**



Open patients: **403,088**


Level of evidence



■ 1b - RCT

■ 2b - Prospective cohort studies

■ 2c - Large database studies



WHAT DOES THE LITERATURE SHOW?

Systematic literature review & meta-analysis results: Robotic-assisted vs. VATS lobectomy



Favors robotic-assisted

- ↓ Conversion to open surgery by **53%**
- ↓ Length of hospital stay by **half a day**



Comparable outcomes


- ≈ ICU admission
- ≈ ICU length of stay
- ≈ Ventilation support >48 hours
- ≈ Respiratory failure
- ≈ Blood transfusion
- ≈ Prolonged air leak
- ≈ Chest tube duration
- ≈ Atelectasis
- ≈ Positive surgical margin
- ≈ Lymph node yield
- ≈ Nodal upstaging
- ≈ Nodal stations sampled
- ≈ 30-day readmission
- ≈ 30-day post-operative complication
- ≈ 30-day reoperation
- ≈ 30-day mortality



Favors VATS

- ↓ Operative time by **12 min**

Data collected through: December 31, 2022



WHAT DOES THE LITERATURE SHOW?

Systematic literature review & meta-analysis results: Robotic-assisted vs. open lobectomy



Favors robotic-assisted

- ↓ Blood transfusion by **32%**
- ↓ Positive surgical margin by **16%**
- ↓ Lymph node yield by **1 lymph node**
- ↓ ICU length of stay by **half a day**
- ↓ ICU admission by **33%**
- ↓ Length of hospital stay by **2 days**
- ↓ 30-day postoperative complication by **27%**
- ↓ 30-day mortality by **35%**



Comparable outcomes

- ≈ Respiratory failure
- ≈ Prolonged air leak
- ≈ Nodal upstaging
- ≈ Nodal stations sampled
- ≈ 30-day reoperation
- ≈ 30-day readmission



Favors open

- ↓ Operative time by **40 min**

Data collected through: December 31, 2022

Evidence Navigator: Lobectomy Technical Slides

Systematic literature review & meta-analysis
as of December 31, 2022

Lobectomy: Literature search methods as of December 31, 2022

Monthly searches were conducted in PubMed, Scopus and Embase.

All citations were exported into a reference management system. Duplications were removed. Titles, abstracts and keywords were reviewed for literature review inclusion by the Global Evidence Management team.

All robotic-assisted lung lobectomies were performed with the da Vinci® surgical systems. Publications were identified according to inclusion and exclusion criteria described.

Meta-analysis was performed using RevMan or R software.


48 publications

175,043 patients who underwent RAS

577,880 patients who underwent video-assisted thoracoscopic surgery (VATS)

403,088 patients who underwent open surgery

Level of evidence



Criteria phase	Details
Identification phase	All unique PubMed, Scopus, and Embase references identified N=7303 December 31, 2022
Inclusion criteria	
1. Robotic-assisted lobectomy for cancer	Da Vinci® robotic-assisted lobectomy N=1217 (excluded N=2891 duplicates & N=3195 not DV lobectomy for lung cancer)
2. Year ≥ 2010	Articles published ≥ 2010 N=1213 (excluded N=4 year<2010)
3. LOE = 1b, 2b, 2c	Articles with LOE = 1b, 2b, 2c N=117 (excluded N=1096 not LOE 1b/2b/2c))
4. Comparative cohort studies n≥20 (robotic-assisted vs. VATS and/or open surgery)	Comparator cohorts N=101 (excluded N=16 sample<20)
Exclusion criteria	N=52 excluded publications:
1. Not in English	N=0 (EC#1)
2. Paper reports on a pediatric population	N=0 (EC#2)
3. Publication is an HTA that was not published in a peer-reviewed journal	N=0 (EC#3)
4. Alternate technique/approach (e.g. single-port)	N=0 (EC#4)
5. No stratified analysis by study arm (e.g., combines results from robotic, VATS and/or open cohorts)	N=33 (EC#5)
6. Lobectomy data mixed with lung mediastinal resection (e.g., data from multiple surgical procedures combined)	N=2 (EC#6)
7. Original research study does not provide quantitative results for outcomes of interest (i.e., operative time, conversions, estimated blood loss and/or transfusions, complications, length of hospital stay, mortality, etc.)	N=9 (EC#7)
8. Original research publication includes redundant patient population and similar conclusions	N=8 (EC#8)

Lobectomy publications included in review: N=49

Lobectomy publications included in meta-analysis: N=48


Robotic-assisted vs. VATS Lobectomy

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring VATS surgery

Compared to VATS lobectomy, the evidence for **robotic-assisted lobectomy** demonstrates:

- 53% less likely to convert to open surgery.
- Comparable ventilation support > 48hrs rate
- Comparable respiratory failure rate
- Comparable positive surgical margin rate
- Comparable mortality rate within 30-days of surgery
- Comparable blood transfusions rate
- Comparable readmissions rate within 30-days of surgery
- Comparable atelectasis rate
- Comparable nodal upstaging rate
- Comparable postoperative complications rate within 30-days of surgery
- Comparable ICU admission rate
- Comparable reoperations rate within 30-days of surgery
- Comparable prolonged air leak rate



Outcome	Robotic-assisted, n	VATS, n	Effect size 95% CI	P-value
Lobectomy binary variables (to December 31, 2022)				
Conversions, n	1, 3, 4, 5, 6, 9, 10, 12, 13, 17, 18, 26, 27, 28, 29, 30, 35, 36, 37, 39, 40, 44, 46, 48			
Subtotal	108724	315079	0.47 [0.39, 0.56]	p<0.01
Random, Heterogeneity: p<0.00001; I ² =96%				
Ventilation Support 48 hr, n	23, 30, 35, 38			
Subtotal	12575	28997	0.92 [0.82, 1.04]	p=0.19
Fixed, Heterogeneity: p=0.35; I ² =9%				
Respiratory Failure, n	1, 2, 23, 30, 36, 41, 42, 47			
Subtotal	23727	59475	0.93 [0.69, 1.24]	p=0.61
Random, Heterogeneity: p<0.00001; I ² =83%				
Positive Surgical Margins, n	3, 13, 25, 29, 33, 34, 39			
Subtotal	9498	18834	0.93 [0.83, 1.04]	p=0.21
Fixed, Heterogeneity: p=0.29; I ² =19%				
Mortality, n	1, 2, 3, 4, 5, 6, 7, 8, 9, 12, 13, 17, 19, 21, 23, 25, 26, 28, 30, 31, 32, 33, 36, 37, 38, 39, 41, 45, 47			
Subtotal	80824	247187	0.96 [0.83, 1.10]	p=0.52
Random, Heterogeneity: p<0.01; I ² =55%				
Blood Transfusions, n	1, 2, 17, 23, 28, 30, 37, 44, 45, 47			
Subtotal	23254	55042	0.96 [0.89, 1.04]	p=0.35
Fixed, Heterogeneity: p=0.61; I ² =0%				
Readmissions, n	3, 4, 5, 9, 11, 13, 17, 25, 33, 35, 36, 39, 41, 45, 46, 47			
Subtotal	30511	63690	0.96 [0.91, 1.02]	p=0.24
Fixed, Heterogeneity: p=0.94; I ² =0%				
Atelectasis, n	23, 30, 38, 42, 44, 46			
Subtotal	4950	16111	0.98 [0.86, 1.12]	p=0.79
Fixed, Heterogeneity: p=0.62; I ² =0%				
Nodal Upstaging, n	15, 19, 23, 25			
Subtotal	8011	28981	0.99 [0.92, 1.08]	p=0.9
Fixed, Heterogeneity: p=0.12; I ² =43%				
Postoperative Complications, n	1, 6, 8, 9, 17, 19, 28, 29, 30, 32, 36, 45, 46, 47, 48			
Subtotal	24008	82736	1.00 [0.91, 1.10]	p=0.99
Random, Heterogeneity: p<0.00001; I ² =77%				
Intensive Care (ICU) Admissions, n	9, 23, 28, 29, 35, 45, 46, 47			
Subtotal	12638	30010	1.04 [0.87, 1.25]	p=0.66
Fixed, Heterogeneity: p<0.0001; I ² =77%				
Reoperations, n	23, 29, 35, 44, 45, 47			
Subtotal	9715	27090	1.08 [0.94, 1.24]	p=0.29
Fixed, Heterogeneity: p=0.66; I ² =0%				
Prolonged Air Leak, n	2, 14, 17, 23, 30, 36, 41, 42, 44, 46, 47			
Subtotal	19680	55437	1.18 [0.99, 1.40]	p=0.06
Random, Heterogeneity: p<0.0001; I ² =75%				


Robotic-assisted vs. VATS Lobectomy

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring VATS surgery

Compared to VATS lobectomy, the evidence for **robotic-assisted lobectomy** demonstrates:

- Significantly shorter length of hospital stay by an average of 0.47 days (11 hours)
- Comparable number of nodal stations sampled
- Comparable ICU length of stay
- Comparable chest tube duration
- Comparable lymph node yield
- Significantly longer operative time by an average of 12.35 minutes



Outcome	Robotic-assisted, n	VATS, n	Effect size 95% CI	P-value
Lobectomy continuous variables (to December 31, 2022)				
LOS, days 1, 2, 4, 6, 7, 9, 13, 17, 19, 22, 23, 25, 26, 28, 29, 30, 32, 33, 35, 36, 38, 39, 41, 42, 44, 46, 47, 48				
Subtotal	89544	252845	-0.47 [-0.65, -0.28]	p<0.01
Random, Heterogeneity: p<0.00001; I ² =99%				
Nodal Stations Sampled, n 29, 46				
Subtotal	83	84	-0.15 [-2.40, 2.10]	p=0.9
Random, Heterogeneity: p<0.00001; I ² =97%				
ICU LOS, days 9, 28, 44, 46				
Subtotal	2997	2998	-0.06 [-0.21, 0.10]	p=0.49
Fixed, Heterogeneity: p=0.54; I ² =0%				
Chest Tube Duration, days 6, 9, 17, 44, 46, 48				
Subtotal	450	458	-0.0001 [-0.61, 0.61]	p=1.00
Random, Heterogeneity: p<0.00001; I ² =87%				
LNY, n 12, 13, 17, 25, 29, 33, 38, 39, 46, 47, 48				
Subtotal	13615	29827	0.18 [-0.24, 0.60]	p=0.39
Random, Heterogeneity: p<0.0001; I ² =87%				
Operative Time, min 6, 9, 17, 23, 28, 29, 30, 36, 42, 44, 46, 47, 48				
Subtotal	10039	25476	12.35 [4.62, 20.08]	p=0.002
Random, Heterogeneity: p<0.00001; I ² =91%				


Robotic-assisted vs. open lobectomy

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open lobectomy, the evidence for **robotic-assisted lobectomy** demonstrates:

- 35% lower likelihood of mortality within 30-days of surgery
- 33% less likely to be admitted to ICU
- 32% less likely to receive a blood transfusion
- 27% less likely to experience a postoperative complication within 30-days of surgery
- 16% less likely to have a positive surgical margin
- Comparable readmissions rate within 30-days of surgery
- Comparable nodal upstaging rate
- Comparable respiratory failure rate
- Comparable reoperations rate within 30-days of surgery
- Comparable prolonged air leak rate



Outcome	Robotic-assisted, n	Open, n	Effect size 95% CI	P-value
Lobectomy binary variables (to December 31, 2022)				
Mortality, n				
Subtotal	54743	245228	0.65 [0.54, 0.78]	p<0.01
Random, Heterogeneity: p<0.01; I ² =71%				
Intensive Care (ICU) Admissions, n				
Subtotal	2694	3554	0.67 [0.60, 0.75]	p<0.01
Fixed, Heterogeneity: p=0.59; I ² =0%				
Blood Transfusions, n				
Subtotal	5439	6299	0.68 [0.51, 0.90]	p=0.007
Random, Heterogeneity: p=0.02; I ² =71%				
Postoperative Complications, n				
Subtotal	5938	7617	0.73 [0.68, 0.78]	p<0.01
Fixed, Heterogeneity: p=0.20; I ² =30%				
Positive Surgical Margins, n				
Subtotal	16184	44139	0.84 [0.75, 0.92]	p=0.0005
Fixed, Heterogeneity: p=0.84; I ² =0%				
Readmissions, n				
Subtotal	18849	56241	0.92 [0.85, 1.00]	p=0.05
Fixed, Heterogeneity: p=0.83; I ² =0%				
Nodal Upstaging, n				
Subtotal	15486	64187	0.94 [0.89, 1.00]	p=0.05
Fixed, Heterogeneity: p=0.06; I ² =59%				
Respiratory Failure, n				
Subtotal	4704	11276	0.41 [0.10, 1.66]	p=0.21
Random, Heterogeneity: p<0.0001; I ² =98%				
Reoperations, n				
Subtotal	426	1283	0.60 [0.25, 1.42]	p=0.24
Fixed, Heterogeneity: p=0.71; I ² =0%				
Prolonged Air Leak, n				
Subtotal	4780	11348	1.47 [0.81, 2.65]	p=0.20
Random, Heterogeneity: p<0.00001; I ² =92%				


Robotic-assisted vs. open lobectomy

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open lobectomy, the evidence for **robotic-assisted lobectomy** demonstrates:

- Significantly shorter hospital stay by an average of 1.98 days
- Significant difference in lymph node yield by an average of 1.2 nodes
- Significantly shorter ICU length of stay by an average of 0.49 days
- Comparable number of nodal stations sampled
- Significantly longer operative time by an average of 39.61 minutes



Outcome	Robotic-assisted, n	Open, n	Effect size 95% CI	P-value
Lobectomy continuous variables (to December 31, 2022)				
LOS, days	7, 13, 16, 19, 22, 25, 26, 28, 29, 30, 33, 39, 41, 43			
Subtotal	64241	242272	-1.98 [-2.86, -1.11]	p<0.01
Random, Heterogeneity: p<0.00001; I²=100%				
LNY, n	12, 13, 15, 25, 29, 33, 39, 43			
Subtotal	85802	21655	-1.20 [-1.64, -0.75]	p<0.01
Random, Heterogeneity: p<0.00001; I²=95%				
ICU LOS, days	1,9,11,18,19			
Subtotal	2326	2326	-0.49 [-0.74, -0.23]	p=0.002
Fixed, Heterogeneity: p=0.44; I²=0%				
Nodal Stations Sampled, n	15, 29			
Subtotal	85	88	-0.58 [-1.36, 0.21]	p=0.15
Random, Heterogeneity: p=0.04; I²=77%				
Operative Time, min	16, 20, 28, 29, 30			
Subtotal	5248	5248	39.61 [22.57, 56.65]	p<0.01
Random, Heterogeneity: p<0.00001; I²=95%				

Lobectomy bibliography (1 of 3)

December 31, 2022

1. Alvarado, C. E., S. G. Worrell, A. L. Sarode, B. Jiang, S. J. Halloran, L. M. Argote-Greene, P. A. Linden and C. W. Towe (2022). "Comparing Thoracoscopic and Robotic Lobectomy Using a Nationally Representative Database." *Am Surg*: 31348221148347.
2. Alwatari, Y., J. Khoraki, L. G. Wolfe, B. Ramamoorthy, N. Wall, C. Liu, W. Julliard, C. A. Puig and R. D. Shah (2022). "Trends of utilization and perioperative outcomes of robotic and video-assisted thoracoscopic surgery in patients with lung cancer undergoing minimally invasive resection in the United States." *JTCVS Open* 12: 385-398.
3. Baig, M. Z., S. S. Razi, K. Agyabeng-Dadzie, S. Stroeve, Z. Muslim, J. Weber, L. J. Herrera and F. Y. Bhora (2022). "Robotic-assisted thoracoscopic Surgery (RATS) Demonstrates a Lower Rate of Conversion to Thoracotomy Than Video-assisted thoracoscopic Surgery (VATS) for Complex Lobectomies." *Eur J Cardiothorac Surg* 62(3).
4. Bailey, K. L., N. Merchant, Y. J. Seo, D. Elashoff, P. Benharash and J. Yanagawa (2019). "Short-Term Readmissions After Open, Thoracoscopic, and Robotic Lobectomy for Lung Cancer Based on the Nationwide Readmissions Database." *World J Surg* 43(5): 1377-1384.
5. Cui, Y., E. L. Grogan, S. A. Deppen, F. Wang, P. P. Massion, C. E. Bailey, W. Zheng, H. Cai and X. O. Shu (2020). "Mortality for Robotic- vs Video-Assisted Lobectomy-Treated Stage I Non-Small Cell Lung Cancer Patients." *JNCI Cancer Spectr* 4(5): pkaa028.
6. Duclos, G., A. Charvet, N. Resseguier, D. Trousse, X. B. D'Journo, L. Zieleskiewicz, P. A. Thomas and M. Leone (2018). "Postoperative morphine consumption and anaesthetic management of patients undergoing video-assisted or robotic-assisted lung resection: A prospective, propensity score-matched study." *Journal of Thoracic Disease* 10(6): 3558-3567.
7. Erhunmwunsee, L., P. Bhandari, E. Sosa, M. Sur, P. H. G. Ituarte and N. S. Lui (2020). "Socioeconomic, rural, and insurance-based inequities in robotic lung cancer resections." *Video-Assisted Thoracic Surgery* 5(june): 1-13.
8. Glenn, Z. F., M. Zubair, L. Hussain and K. Grannan (2019). "Comparison of pulmonary lobectomies using robotic and video-assisted thoracoscopic approaches: results from 2010-2013 National Inpatient Sample." *J Cardiovasc Surg (Torino)* 60(4): 526-531.
9. Gonde, H., M. Laurent, A. Gillibert, O. M. Sarsam, R. Varin, G. Grimandi, C. Peillon and J. M. Baste (2017). "The affordability of minimally invasive procedures in major lung resection: a prospective study." *Interact Cardiovasc Thorac Surg* 25(3): 469-475.
10. Hendriksen, B. S., C. S. Hollenbeak, M. D. Taylor and M. F. Reed (2019). "Minimally Invasive Lobectomy Modality and Other Predictors of Conversion to Thoracotomy." *Innovations (Phila)* 14(4): 342-352.
11. Hendriksen, B. S., M. F. Reed, M. D. Taylor and C. S. Hollenbeak (2019). "Readmissions After Lobectomy in an Era of Increasing Minimally Invasive Surgery: A Statewide Analysis." *Innovations (Phila)* 14(5): 453-462.
12. Hennon, M. W., L. H. DeGraaff, A. Groman, T. L. Demmy and S. Yendamuri (2020). "The association of nodal upstaging with surgical approach and its impact on long-term survival after resection of non-small-cell lung cancer." *Eur J Cardiothorac Surg* 57(5): 888-895.
13. Herb, J. N., D. G. Kindell, P. D. Strassle, K. B. Stitzenberg, B. E. Haithcock, G. N. Mody and J. M. Long (2021). "Trends and Outcomes in Minimally Invasive Surgery for Locally Advanced Non-Small Cell Lung Cancer with N2 Disease." *Semin Thorac Cardiovasc Surg* 33(2): 547-555.
14. Hoeijmakers, F., K. J. Hartemink, A. F. Verhagen, W. H. Steup, E. Marra, W. F. B. Röell, D. J. Heineman, W. H. Schreurs, R. A. E. M. Tollenaar and M. W. J. M. Wouters (2021). "Variation in incidence, prevention and treatment of persistent air leak after lung cancer surgery." *Eur J Cardiothorac Surg* 61(1): 110-117.
15. Huang, J., C. Li, H. Li, F. Lv, L. Jiang, H. Lin, P. Lu, Q. Luo and W. Xu (2019). "Robot-assisted thoracoscopic surgery versus thoracotomy for c-N2 stage NSCLC: Short-term outcomes of a randomized trial." *Translational Lung Cancer Research* 8(6): 951-958.
16. Huang, J., Y. Tian, C. Li, Y. Shen, H. Li, F. Lv, H. Lin, P. Lu, J. Lin, C. Lau, R. M. Terra, L. Jiang and Q. Luo (2021). "Robotic-assisted thoracic surgery reduces perioperative complications and achieves a similar long-term survival profile as posterolateral thoracotomy in clinical N2 stage non-small cell lung cancer patients: A multicenter, randomized, controlled trial." *Translational Lung Cancer Research* 10(11): 4281-4292.
17. Jin, R., Y. Zheng, Y. Yuan, D. Han, Y. Cao, Y. Zhang, C. Li, J. Xiang, Z. Zhang, Z. Niu, T. Lerut, J. Lin, A. E. Abbas, A. Pardolesi, T. Suda, D. Amore, S. Schraag, C. Aigner, J. Li, J. Che, J. Hang, J. Ren, L. Zhu and H. Li (2022). "Robotic-assisted Versus Video-assisted Thoracoscopic Lobectomy: Short-term Results of a Randomized Clinical Trial (RVlob Trial)." *Ann Surg* 275(2): 295-302.
18. Kamel, M. K., A. N. Sholi, S. W. Harrison, B. Lee, J. L. Port, N. K. Altorki and B. M. Stiles (2022). "Minimally Invasive Surgery for Lung Cancer Following Neoadjuvant Therapy in the United States." *J Laparoendosc Adv Surg Tech A* 32(8): 860-865.
19. Kent, M., T. Wang, R. Whyte, T. Curran, R. Flores and S. Gangadharan (2014). "Open, video-assisted thoracic surgery, and robotic lobectomy: review of a national database." *Ann Thorac Surg* 97(1): 236-244.

Lobectomy bibliography (2 of 3)

December 31, 2022

20. Lacroix, V., D. Kahn, P. Matte, T. Pieters, P. Noirhomme, A. Poncelet and A. Steyaert (2021). "Robotic-assisted lobectomy favors early lung recovery versus limited thoracotomy." *Thoracic and Cardiovascular Surgeon* 69(6): 557-563.
21. Liberman, M., E. Goudie, C. Morse, W. Hanna, N. Evans, K. Yasufuku, J. Sampalis and V. P. S. W. Group (2019). "Prospective, multicenter, international phase 2 trial evaluating ultrasonic energy for pulmonary artery branch sealing in video-assisted thoracoscopic surgery lobectomy." *J Thorac Cardiovasc Surg*.
22. Lo, T., R. Schiller, K. Raghunathan, V. Krishnamoorthy, O. K. Jawitz, S. Pyati, T. Van De Ven, R. R. Bartz, A. Thompson and T. Ohnuma (2021). "Changes in analgesic strategies for lobectomy from 2009 to 2018." *JTCVS Open* 6: 224-236.
23. Louie, B. E., J. L. Wilson, S. Kim, R. J. Cerfolio, B. J. Park, A. S. Farivar, E. Vallieres, R. W. Aye, W. R. Burfeind, Jr. and M. I. Block (2016). "Comparison of Video-Assisted Thoracoscopic Surgery and Robotic Approaches for Clinical Stage I and Stage II Non-Small Cell Lung Cancer Using The Society of Thoracic Surgeons Database." *Ann Thorac Surg* 102(3): 917-924.
24. Martin, J. T., E. B. Durbin, L. Chen, T. Gal, A. Mahan, V. Ferraris and J. Zwischenberger (2016). "Nodal Upstaging During Lung Cancer Resection Is Associated With Surgical Approach." *Annals of Thoracic Surgery* 101(1): 238-244.
25. Merritt, R. E., M. Abdel-Rasoul, D. M. D'Souza and P. J. Kneuert (2023). "Lymph Node Upstaging for Robotic, Thoracoscopic, and Open Lobectomy for Stage T2-3N0 Lung Cancer." *Ann Thorac Surg* 115(1): 175-182.
26. Mitzman, B., X. Wang, B. Haaland and T. K. Varghese, Jr. (2022). "Trends and factors affecting approach choice to pulmonary resection." *Journal of Surgical Oncology* 126(3): 599-608.
27. Muslim, Z., S. Stroeve, K. Poulikidis, J. F. Weber, C. P. Connery, L. J. Herrera and F. Y. Bhora (2022). "Conversion to Thoracotomy in Non-Small Cell Lung Cancer: Risk Factors and Perioperative Outcomes." *Innovations: Technology and Techniques in Cardiothoracic and Vascular Surgery* 17(2): 148-155.
28. Nguyen, D. M., I. S. Sarkaria, C. Song, R. M. Reddy, N. Villamizar, L. J. Herrera, L. Shi, E. Liu, D. Rice and D. S. Oh (2020). "Clinical and economic comparative effectiveness of robotic-assisted, video-assisted thoracoscopic, and open lobectomy." *Journal of Thoracic Disease* 12(3): 296-306.
29. Novellis, P., P. Maisonneuve, E. Dieci, E. Voulaz, E. Bottoni, S. Di Stefano, M. Solinas, A. Testori, U. Cariboni, M. Alloisio and G. Veronesi (2021). "Quality of life, postoperative pain, and lymph node dissection in a robotic approach compared to vats and open for early stage lung cancer." *Journal of Clinical Medicine* 10(8).
30. Oh, D. S., R. M. Reddy, M. L. Gorrepati, S. Mehendale and M. F. Reed (2017). "Robotic-Assisted, Video-Assisted Thoracoscopic and Open Lobectomy: Propensity-Matched Analysis of Recent Premier Data." *Ann Thorac Surg* 104(5): 1733-1740.
31. Park, B. J., J. M. Snider, N. R. Bates, S. D. Cassivi, G. K. Jett, J. R. Sonett and E. M. Toloza (2016). "Prospective evaluation of biodegradable polymeric sealant for intraoperative air leaks." *Journal of cardiothoracic surgery* 11(1): 168.
32. Paul, S., J. Jalbert, A. J. Isaacs, N. K. Altorki, O. W. Isom and A. Sedrakyan (2014). "Comparative effectiveness of robotic-assisted vs thoracoscopic lobectomy." *Chest* 146(6): 1505-1512.
33. Rajaram, R., S. Mohanty, D. J. Bentrem, E. S. Pavey, D. D. Odell, A. Bharat, K. Y. Bilimoria and M. M. DeCamp (2017). "Nationwide Assessment of Robotic Lobectomy for Non-Small Cell Lung Cancer." *Ann Thorac Surg* 103(4): 1092-1100.
34. Ray, M. A., N. R. Faris, M. P. Smeltzer, C. Fehnel, C. Houston-Harris, P. Levy, L. Wiggins, V. Sachdev, T. Robbins, D. Spencer and R. U. Osarogiagbon (2018). "Effectiveness of Implemented Interventions on Pathologic Nodal Staging of Non-Small Cell Lung Cancer." *Annals of Thoracic Surgery* 106(1): 228-234.
35. Ricciardi, R., R. N. Goldstone, T. Francone, M. Wszolek, H. Auchincloss, A. de Groot, I. F. Shih and Y. Li (2022). "Healthcare Resource Utilization After Surgical Treatment of Cancer: Value of Minimally Invasive Surgery." *Surgical Endoscopy* 36(10): 7549-7560.
36. Seder, C. W., F. Farrokhyar, R. Nayak, J. M. Baste, Y. Patel, J. Agzarian, C. J. Finley, Y. Shargall, P. A. Thomas, M. Dahan, J. P. Verhoye, F. Mbadinga and W. C. Hanna (2022). "Robotic vs. Thoracoscopic Anatomic Lung Resection in Obese Patients: A Propensity Adjusted Analysis." *Ann Thorac Surg* 114(5): 1879-1885.
37. Servais, E. L., D. L. Miller, D. Thibault, M. G. Hartwig, A. S. Kosinski, C. T. Stock, T. Price, S. M. Quadri, R. S. D'Agostino and W. R. Burfeind (2022). "Conversion to Thoracotomy During Thoracoscopic versus Robotic Lobectomy: Predictors and Outcomes." *Ann Thorac Surg* 114(2): 409-417.
38. Sesti, J., R. C. Langan, J. Bell, A. Nguyen, A. L. Turner, P. Hilden, K. Leshchuk, M. Dabrowski and S. Paul (2020). "A Comparative Analysis of Long-term Survival of Robotic vs. Thoracoscopic Lobectomy." *Ann Thorac Surg* 110(4): 1139-1146.
39. Shagabayeva, L., B. Fu, N. Panda, A. L. Potter, H. G. Auchincloss, A. Mansur, C. F. Jeffrey Yang and L. Schumacher (2023). "Open, Video- and Robot-Assisted Thoracoscopic Lobectomy for Stage II-IIIa Non-Small Cell Lung Cancer." *Ann Thorac Surg* 115(1): 184-190.

Lobectomy bibliography (3 of 3)

December 31, 2022

40. Shah, P. C., A. de Groot, R. Cerfolio, W. C. Huang, K. Huang, C. Song, Y. Li, U. Kreaden and D. S. Oh (2022). "Impact of type of minimally invasive approach on open conversions across ten common procedures in different specialties." *Surg Endosc* 36(8): 6067-6075.
41. Subramanian, M. P., J. Liu, W. C. Chapman, M. A. Olsen, Y. Yan, Y. Liu, T. R. Semenkovich, B. F. Meyers, V. Puri and B. D. Kozower (2019). "Utilization Trends, Outcomes, and Cost in Minimally Invasive Lobectomy." *The Annals of thoracic surgery* 108(6): 1648-1655.
42. Swanson, S. J., D. L. Miller, R. J. McKenna, Jr., J. Howington, M. B. Marshall, A. C. Yoo, M. Moore, C. L. Gunnarsson and B. F. Meyers (2014). "Comparing robot-assisted thoracic surgical lobectomy with conventional video-assisted thoracic surgical lobectomy and wedge resection: results from a multihospital database (Premier)." *J Thorac Cardiovasc Surg* 147(3): 929-937.
43. Tang, A., S. Raja, A. C. Bribresco, D. P. Raymond, M. Sudarshan, S. C. Murthy and U. Ahmad (2020). "Robotic Approach Offers Similar Nodal Upstaging to Open Lobectomy for Clinical Stage I Non-small Cell Lung Cancer." *Ann Thorac Surg* 110(2): 424-433.
44. Terra, R. M., P. H. X. N. Araujo, L. L. Lauricella, J. R. M. Campos, J. R. M. Trindade and P. M. Pêgo-Fernandes (2022). "A Brazilian randomized study: Robotic-Assisted vs. Video-assisted lung lobectomy Outcomes (BRAVO trial)." *Jornal brasileiro de pneumologia : publicacao oficial da Sociedade Brasileira de Pneumologia e Tisilogia* 48(4): e20210464.
45. Veluswamy, R. R., S. A. Whittaker Brown, G. Mhango, K. Sigel, D. G. Nicastrì, C. B. Smith, M. Bonomi, M. D. Galsky, E. Taioli, A. I. Neugut and J. P. Wisnivesky (2020). "Comparative Effectiveness of Robotic-Assisted Surgery for Resectable Lung Cancer in Older Patients." *Chest* 157(5): 1313-1321.
46. Veronesi, G., A. E. S. Abbas, P. Muriana, R. Lembo, E. Bottoni, G. Perroni, A. Testori, E. Dieci, C. T. Bakhos, S. Car, L. Luzzi, M. Alloisio and P. Novellis (2021). "Perioperative Outcome of Robotic Approach Versus Manual Videothoroscopic Major Resection in Patients Affected by Early Lung Cancer: Results of a Randomized Multicentric Study (ROMAN Study)." *Frontiers in Oncology* 11.
47. Williams, A. M., L. Zhao, T. R. Grenda, R. G. Kathawate, B. E. Biesterveld, U. F. Bhatti, P. W. Carrott, K. H. Lagisetty, A. C. Chang, W. Lynch, J. Lin and R. M. Reddy (2022). "Higher Long-term Quality of Life Metrics After Video-Assisted Thoracoscopic Surgery Lobectomy Compared With Robotic-Assisted Lobectomy." *Ann Thorac Surg* 113(5): 1591-1597.
48. Zhou, J. C., W. P. Wang, S. Q. Wu, J. L. Wang and W. H. Li (2022). "Clinical Efficacy of Thoracoscopic Surgery with the da Vinci Surgical System versus Video-Assisted Thoracoscopic Surgery for Lung Cancer." *Journal of Oncology* 2022.

Disclosures

Important Safety Information

(US) Serious complications may occur in any surgery, including da Vinci surgery, up to and including death. Serious risks include, but are not limited to, injury to tissues and organs and conversion to other surgical techniques which could result in a longer operative time and/or increased complications. For summary of the risks associated with surgery refer to www.intuitive.com/safety.

Da Vinci Xi®/da Vinci X® system precaution statement

The demonstration of safety and effectiveness for the representative specific procedures did not include evaluation of outcomes related to the treatment of cancer (overall survival, disease-free survival, local recurrence), except for radical prostatectomy which was evaluated for overall survival, or treatment of the patient's underlying disease/condition. Device usage in all surgical procedures should be guided by the clinical judgment of an adequately trained surgeon.

(EU) Da Vinci X & Xi Surgical Systems

The Intuitive Surgical Endoscopic Instrument Control Systems (da Vinci X and da Vinci Xi Surgical Systems) are intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments during urologic surgical procedures, general laparoscopic surgical procedures, gynecologic laparoscopic surgical procedures, general thoracoscopic surgical procedures, and trans-oral otolaryngology surgical procedures restricted to benign tumors and malignant tumors classified as T1 and T2, and for benign base of tongue resection procedures. The systems are indicated for adult and pediatric use (except for trans-oral otolaryngology surgical procedures). They are intended to be used by trained physicians in an operating room environment.

The da Vinci X and da Vinci Xi Surgical Systems are class IIb medical devices CE marked (CE 2460) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

For product intended use and/or indications for use, risks, cautions, and warnings and full prescribing information, refer to the associated user manual(s) or visit <https://manuals.intuitivesurgical.com/market>. Some products, features or technologies may not be available in all countries. Product availability is subject to regulatory approval in the specific market. Please contact your local Intuitive representative for product availability in your region.

Individual outcomes may depend on a number of factors—including but not limited to—patient characteristics, disease characteristics, and/or surgeon experience.

Privacy Notice: Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2025 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos, including Intuitive, Da Vinci, and Ion, are trademarks or registered trademarks of Intuitive Surgical or their respective owner.

INTUITIVE

intuitive.com