

Evidence Navigator: Hysterectomy For Endometrial Cancer

Systematic literature review summary
as of December 31, 2022

Purpose

The Evidence Navigator is a slide presentation representing a summary of the meta-analysis of the highest level of evidence available specific to a given procedure and published as of a particular date. It is created by the Global Evidence Management team within Global Access, Value and Economics (GAVE). It includes information that is available in the public domain. It is a systematic review and meta-analysis of the peer-reviewed literature based on a timeframe within which a literature search has been conducted according to a set of concise inclusion and exclusion criteria. The results of the meta-analysis are presented in the form of forest plots summarized for each outcome according to a comparator and surgical approach of interest. The summary results are reflective of a specific period in time and are subject to change with increasing literature. All of the robotic-assisted surgery procedures mentioned within the Evidence Navigator were performed using a da Vinci® surgical system.

Statistical analysis

All summary measures are shown as odds ratios, risk ratios or risk differences when describing binary outcomes, or as weighted mean differences or standardized mean differences when describing continuous outcomes. Weighting is based on the study sample size and variability of the outcome. A random effect model is used if heterogeneity is statistically significant, otherwise a fixed effect model is used. The Mantel Haenszel summary statistic is used for the overall results. The meta-analysis is performed with RevMan 5.4 (Review Manager, Version 5.4. Copenhagen: The Nordic Cochrane Centre, The Cochrane Collaboration, 2014) or R software (R Foundation for Statistical Computing, Vienna, Austria. URL <https://www.R-project.org/>).

Interpretation notes

When the effect size is measured as a standardized mean difference (SMD), or a risk difference (RD), it is not possible to provide a quantitative conclusion. In such cases, a qualitative conclusion is given with reference to its statistical significance. In some instances, studies may contain some overlapping patient populations. A redundancy check is performed in order to minimize this overlap and bias due to over-reporting.

Glossary

RAS	robotic-assisted surgery
Lap	laparoscopic surgery
LOE	level of evidence
HTA	health technology assessment
RCT	randomized controlled trial
OR	odds ratio
MD	mean difference
LNY	lymph node yield
LVSI	lymphovascular space invasion

WMD	weighted mean difference
RD	risk difference
SMD	standardized mean difference
95% CI	95% confidence interval
I²	test statistic for heterogeneity
EBL	estimated blood loss
LOS	length of hospital stay
PSM	positive surgical margins

Evidence Navigator: Hysterectomy for endometrial cancer Summary Slides

Systematic literature review summary
as of December 31, 2022



WHAT DOES THE LITERATURE SHOW?

Systematic literature review: Hysterectomy for endometrial cancer — clinical outcomes

Inclusion criteria

Robotic-assisted hysterectomy for endometrial cancer performed with a da Vinci surgical system

January 1, 2010 – December 31, 2022

Level of Evidence = 1b, 2b, 2c

RCT, prospective cohort studies, or large database study (with $n \geq 20$ in each cohort)

Exclusion criteria

Not in English

Paper on a pediatric population

Publication is a HTA not published in a peer-reviewed journal

Alternate technique/approach

No stratified analysis by study arm

Hysterectomy data mixed with other procedure(s)

Original research study does not provide quantitative results

Original research publication includes redundant patient population and similar conclusions

34 publications including



Robotic-assisted patients: **350,460**



Laparoscopic patients: **113,246**



Open patients: **219,787**

Level of evidence



■ 1b - RCTs

■ 2b - Prospective cohort studies

■ 2c - Database studies

Data collected through: December 31, 2022



WHAT DOES THE LITERATURE SHOW?

Systematic literature review key points:

Robotic-assisted vs. laparoscopic hysterectomy for endometrial cancer



Favors robotic-assisted

- ↓ Conversions by **55%**
- ↓ Estimated blood loss by **24 mL**
- ↓ Length of stay by **0.5 days**
- ↓ 30-day mortality by **37%**



Comparable outcomes

- ≈ Operative time
- ≈ Blood transfusions
- ≈ Lymph node yield
- ≈ Pelvic lymph node yield
- ≈ Para-aortic lymph node yield
- ≈ Positive lymphovascular space invasion
- ≈ Intraoperative complications
- ≈ 30-day postoperative complications
- ≈ 30-day reoperations
- ≈ 30-day readmissions



Favors laparoscopic

None

Data collected through: December 31, 2022

■ Significant difference favoring
robotic-assisted surgery

■ No significant difference;
comparable outcomes

■ Significant difference favoring
laparoscopic surgery



WHAT DOES THE LITERATURE SHOW?

Systematic literature review key points:

Robotic-assisted vs. open hysterectomy for endometrial cancer



Favors robotic-assisted

- ↓ Estimated blood loss by **152 mL**
- ↓ Length of stay by **2.5 days**
- ↓ Blood transfusions by **71%**
- ↓ 30-day postoperative complications by **62%**
- ↓ 30-day reoperations by **92%**
- ↓ 30-day readmissions by **45%**
- ↓ 30-day mortality by **64%**



Comparable outcomes

- ≈ Operative time
- ≈ Lymph node yield
- ≈ Para-aortic lymph node yield
- ≈ Positive lymphovascular space invasion
- ≈ Positive surgical margins
- ≈ Intraoperative complications



Favors open

- ↑ Pelvic lymph yield by **6 nodes**

Data collected through: December 31, 2022

■ Significant difference favoring robotic-assisted surgery

■ No significant difference; comparable outcomes

■ Significant difference favoring laparoscopic surgery

Evidence Navigator: Hysterectomy for endometrial cancer Technical Slides

**Systematic literature review summary
as of December 31, 2022**

Hysterectomy for Endometrial Cancer: literature search methods

as of December 31, 2022

Monthly searches were conducted in PubMed, Scopus and Embase.

All citations were exported into a reference management system. Duplications were removed. Titles, abstracts and keywords were reviewed for literature review inclusion by Global Evidence Management team.

All robotic-assisted hysterectomies performed with da Vinci® surgical systems. Publications were identified according to inclusion and exclusion criteria described.

Meta-analysis was performed using RevMan or R software.

34 publications

350,460 patients who underwent RAS

113,246 patients who underwent laparoscopic surgery

219,787 patients who underwent open surgery

Level of evidence



- 1b - RCTs
- 2b - Prospective cohort studies
- 2c - Database studies

Criteria phase	Details
Identification phase	Unique records identified from PubMed, Scopus, Embase search N=9,440 to December 31, 2022
Inclusion criteria	
1. Robotic-assisted hysterectomy (radical, total, simple) with or without salpingectomy, oophorectomy, and lymphadenectomy for cancer or other gynecologic oncology procedure	N=2,626 (excluded N=6,814)
2. Year≥2010	N=2,623 (excluded N=3)
3. LOE=1b, 2b, 2c	N=253 (excluded N=2,370)
4. Study is an RCT, prospective study or large database study with comparative cohorts (robotic-assisted vs lap and/or open surgery) and sample size N≥20	N=231 (excluded N=22)
Exclusion criteria	N=179 excluded publications:
1. Not in English	N=1 (EC#1)
2. Paper reports on a pediatric population	N=0 (EC#2)
3. Publication is an HTA that was not published in a peer-reviewed journal	N=0 (EC#3)
4. Alternate technique/approach (e.g. single-port)	N=4 (EC#4)
5. No stratified analysis by study arm (e.g. combines results from robotic, laparoscopic and/or open cohorts)	N=113 (EC#5)
6. Hysterectomy cancer data mixed with other procedures (e.g. data from multiple surgical procedures combined)	N=32 (EC#6)
7. Original research study does not provide quantitative results for at least one of the findings relative to the outcomes of interest	N=21 (EC#7)
8. Original research publication includes redundant patient population and similar conclusions	N=8 (EC#8)

Gyn Onc publications: N=52 (34 endometrial cancer)


Robotic-assisted vs. laparoscopic hysterectomy for endometrial cancer

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic hysterectomy, the evidence for **robotic-assisted hysterectomy for endometrial cancer** demonstrates:

- Significantly less estimated blood loss by an average of 24 mL
- Significantly shorter hospital length of stay by an average of 0.5 days
- Comparable lymph node yield
- Comparable pelvic lymph node yield
- Comparable para-aortic lymph node yield
- Comparable operative time



Outcome	Robotic-assisted, n	Laparoscopic, n	Effect size 95% CI	P-value
Endometrial cancer continuous variables (to December 31, 2022)				
EBL, mL ^{6, 16, 17, 21-23, 27}				
Subtotal	826	689	WMD: -23.67 [-45.44, -1.91]	p=0.03
Random, Heterogeneity: p = 0.03; I ² = 53%				
LOS, days ^{6, 8, 9, 11, 13, 16, 19, 21, 27, 32}				
Subtotal	26570	13786	WMD: -0.53 [-0.90, -0.17]	p<0.01
Random, Heterogeneity: p<0.01; I ² = 88%				
LNY, n ^{21, 23, 27}				
Subtotal	128	128	WMD: -1.98 [-7.08, 3.13]	p=0.45
L-R Random, Heterogeneity: p=0.07; I ² = 62%				
Pelvic LNY, n ^{16, 19, 27}				
Subtotal	203	169	WMD: 0.82 [-1.99, 3.63]	p=0.57
L-R Random, Heterogeneity: p=0.03; I ² = 64%				
Para-aortic LNY, n ^{16, 19, 27}				
Subtotal	144	147	WMD: 1.01 [-3.56, 5.58]	p=0.66
L-R Random, Heterogeneity: p<0.01; I ² = 91%				
Operative time, min ^{6, 9, 16, 17, 19, 21, 22, 27}				
Subtotal	2355	2215	WMD: 13.82 [-7.00, 34.63]	p=0.19
Random, Heterogeneity: p<0.01; I ² = 93%				


Robotic-assisted vs. laparoscopic hysterectomy for endometrial cancer

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery
 ■ No significant difference; comparable outcomes
 ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic hysterectomy, the evidence for **robotic-assisted hysterectomy for endometrial cancer** demonstrates:

- 55% less likely to be converted to open surgery
- 37% less likely to experience mortality within 30-days of surgery
- Comparable reoperation rate within 30-days of surgery
- Comparable blood transfusion rate
- Comparable intraoperative complication rate
- Comparable rate of positive lymphovascular space invasion
- Comparable postoperative complication rate within 30-days of surgery
- Comparable readmission rate within 30-days of surgery



Outcome	Robotic-assisted, n	Laparoscopic, n	Effect size 95% CI	P-value
Endometrial cancer binary variables (to December 31, 2022)				
Conversions, n	5-7, 9, 15, 16, 19, 21, 24, 27, 33, 34			
Subtotal	87685	33816	OR: 0.45 [0.29, 0.67]	p<0.01
Random, Heterogeneity: p<0.01; I ² = 86%				
Mortality, n	8, 11, 21, 22, 27, 31, 33, 34			
Subtotal	31824	13897	OR: 0.63 [0.49, 0.82]	p<0.01
Fixed, Heterogeneity: p=0.68; I ² = 0%				
Reoperations, n	16, 21, 31			
Subtotal	1620	1177	OR: 0.78 [0.39, 1.57]	p=0.49
Fixed, Heterogeneity: p=0.25; I ² =25%				
Transfusions, n	6, 9, 10, 16, 19, 21, 30, 31, 33, 34			
Subtotal	12055	8331	OR: 0.90 [0.80, 1.02]	p=0.10
Fixed, Heterogeneity: p=0.37; I ² = 8%				
Intraop complications, n	6, 16, 19, 21, 30, 31, 33, 34			
Subtotal	10096	6696	OR: 0.91 [0.77, 1.07]	p=0.25
Fixed, Heterogeneity: p=0.54; I ² = 0%				
Positive LVSI, n	1, 5, 8			
Subtotal	92379	32590	OR: 0.92 [0.83, 1.02]	p=0.11
Fixed, Heterogeneity: p=0.30; I ² =19%				
Postop complications, n	4, 6, 9, 10, 13, 16, 19, 30, 31, 33, 34			
Subtotal	14007	11930	OR: 0.98 [0.91, 1.06]	p=0.58
Random, Heterogeneity: p<0.01; I ² = 81%				
Readmissions, n	3, 4, 8, 9, 31			
Subtotal	141005	48010	OR: 1.03 [0.95, 1.11]	p=0.47
Fixed, Heterogeneity: p=0.16; I ² =39%				


Robotic-assisted vs. open hysterectomy for endometrial cancer

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery
 ■ No significant difference; comparable outcomes
 ■ Significant difference favoring open surgery

Compared to open hysterectomy, the evidence for **robotic-assisted hysterectomy for endometrial cancer** demonstrates:

- Significantly less estimated blood loss by an average of 152 mL
- Significantly shorter hospital length of stay by an average of 2.5 days
- Comparable para-aortic lymph node yield
- Comparable lymph node yield
- Comparable operative time
- Significantly lower pelvic lymph node yield by an average of 6 nodes



Outcome	Robotic-assisted, n	Open, n	Effect size 95% CI	P-value
Endometrial cancer continuous variables (to December 31, 2022)				
EBL, mL ^{6, 12, 20, 23, 26, 28}				
Subtotal	580	2843	WMD: -151.99 [-225.83, -78.14]	p<0.01
Random, Heterogeneity: p<0.01; I ² = 96%				
LOS, days ^{6, 8, 9, 11, 13, 18, 12, 19, 20, 26, 28, 32}				
Subtotal	30284	35697	WMD: -2.50 [-3.06, -1.93]	p<0.01
Random, Heterogeneity: p<0.01; I ² = 97%				
Para-aortic LNY, n ^{6, 19, 26}				
Subtotal	274	751	WMD: -0.22 [-2.33, 1.89]	p=0.84
O-R Random, Heterogeneity: p = 0.12; I ² = 53%				
LNY, n ^{6, 12, 23, 26, 28}				
Subtotal	347	805	WMD: 0.74 [-5.85, 7.32]	p=0.83
O-R Random, Heterogeneity: p<0.01; I ² = 93%				
Operative time, min ^{6, 9, 12, 19, 20, 26}				
Subtotal	3083	5373	WMD: 12.97 [-14.63, 40.56]	p=0.36
Random, Heterogeneity: p<0.01; I ² = 94%				
Pelvic LNY, n ^{6, 19, 26}				
Subtotal	274	751	WMD: 6.36 [-3.61, 9.11]	p<0.01
O-R Random, Heterogeneity: p = 0.09; I ² = 59%				


Robotic-assisted vs. open hysterectomy for endometrial cancer

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery
 ■ No significant difference; comparable outcomes
 ■ Significant difference favoring open surgery

Compared to open hysterectomy, the evidence for **robotic-assisted hysterectomy for endometrial cancer** demonstrates:

- 92% less likely to experience a reoperation within 30-days of surgery
- 71% less likely to experience a blood transfusion
- 64% less likely to experience mortality within 30-days of surgery
- 62% less likely to experience a postoperative complications within 30-days of surgery
- 45% less likely to experience a readmission within 30-days of surgery
- Comparable positive surgical margins
- Comparable rate of positive lymphovascular space invasion
- Comparable intraoperative complication rate



Outcome	Robotic-assisted, n	Open, n	Effect size 95% CI	P-value
Endometrial cancer binary variables (to December 31, 2022)				
Reoperations, n 18, 26				
Subtotal	2850	14226	OR: 0.08 [0.02, 0.35]	p<0.01
Fixed, Heterogeneity: p=0.18; I ² = 41%				
Transfusions, n 6, 9, 10, 12, 16, 19, 20, 26				
Subtotal	6331	20198	OR: 0.29 [0.22, 0.38]	p<0.01
Random, Heterogeneity: p=0.04; I ² = 49%				
Mortality, n 8, 11, 18, 25, 26				
Subtotal	50057	47634	OR: 0.36 [0.29, 0.43]	p<0.01
Fixed, Heterogeneity: p=0.35; I ² = 10%				
Postop complications, n 4, 6, 9, 10, 12, 13, 18, 19, 20, 26, 28				
Subtotal	8398	23399	OR: 0.38 [0.34, 0.42]	p<0.01
Fixed, Heterogeneity: p=0.06; I ² = 41%				
Readmissions, n 3, 4, 8, 9, 12, 25, 26				
Subtotal	164527	101938	OR: 0.55 [0.46, 0.66]	p<0.01
Random, Heterogeneity: p<0.01; I ² = 85%				
Positive surgical margins, n 1, 25				
Subtotal	25106	21418	OR: 0.86 [0.45, 1.62]	p=0.63
Random, Heterogeneity: p=0.01; I ² = 84%				
Positive LVSI, n 1, 5, 8				
Subtotal	92379	60497	OR: 0.87 [0.60, 1.26]	p=0.46
Random, Heterogeneity: p<0.01; I ² = 99%				
Intraop complications, n 6, 12, 18, 19, 26, 28, 30				
Subtotal	4530	21858	OR: 1.04 [0.66, 1.64]	p=0.85
Random, Heterogeneity: p<0.01; I ² = 78%				

Hysterectomy for endometrial cancer: bibliography (1 of 2)

December 31, 2022

1. Abel, M. K., J. K. Chan, S. Chow, K. Darcy, C. Tian, D. S. Kapp, A. K. Mann and C. I. Liao (2020). "Trends and survival outcomes of robotic, laparoscopic, and open surgery for stage II uterine cancer." *Int J Gynecol Cancer* 30(9): 1347-1355.
2. Albright, B. B., D. Nasioudis, M. E. Byrne, N. A. Latif, E. M. Ko and A. F. Haggerty (2020). "Perioperative outcomes and disparities in utilization of sentinel lymph node biopsy in minimally invasive staging of endometrial cancer." *Gynecol Oncol* 159(3): 758-766. [reports risk of conversion to open with no raw data]
3. Barrington, D. A., C. E. Meade, C. M. Cosgrove, D. E. Cohn and A. S. Felix (2022). "Racial and ethnic disparities in readmission risk following the surgical management of endometrial cancer." *Gynecol Oncol* 166(3): 543-551.
4. Beck, T. L., M. A. Schiff, B. A. Goff and R. R. Urban (2018). "Robotic, Laparoscopic, or Open Hysterectomy: Surgical Outcomes by Approach in Endometrial Cancer." *J Minim Invasive Gynecol* 25(6): 986-993.
5. Bixel, K., D. A. Barrington, M. H. Vetter, A. A. Suarez and A. S. Felix (2022). "Determinants of surgical approach and survival among women with endometrial carcinoma." *J Minim Invasive Gynecol* 29(2): 219-230.
6. Borgfeldt, C., G. Kalapotharakos, K. C. Asciutto, M. Lofgren and T. Hogberg (2016). "A population-based registry study evaluating surgery in newly diagnosed uterine cancer." *Acta Obstet Gynecol Scand* 95(8): 901-911.
7. Bregar, A. J., A. Melamed, E. Diver, J. T. Clemmer, S. Uppal, J. O. Schorge, L. W. Rice, M. G. Del Carmen and J. A. Rauh-Hain (2017). "Minimally Invasive Staging Surgery in Women with Early-Stage Endometrial Cancer: Analysis of the National Cancer Data Base." *Ann Surg Oncol* 24(6): 1677-1687.
8. Cardenas-Goicoechea, J., Y. U. Wang, J. H. Lee, M. Shoraka, S. L. Carbajal-Mamani, D. Fishman, A. N. Riner and J. G. Trevino (2022). "Survival After Minimally Invasive Surgery in Older Women With Endometrial Carcinoma." *Anticancer Res* 42(1): 75-85.
9. Casarin, J., C. Song, F. Multinu, S. Cappuccio, E. Liu, K. A. Butler, G. E. Glaser, W. A. Cliby, C. L. Langstraat, F. Ghezzi, A. Z. Fu and A. Mariani (2020). "Implementing robotic surgery for uterine cancer in the United States: Better outcomes without increased costs." *Gynecol Oncol* 156(2): 451-458.
10. Chan, J. K., A. B. Gardner, K. Taylor, C. A. Thompson, K. Blansit, X. Yu and D. S. Kapp (2015). "Robotic versus laparoscopic versus open surgery in morbidly obese endometrial cancer patients - a comparative analysis of total charges and complication rates." *Gynecol Oncol* 139(2): 300-305.
11. Dubeshter, B., C. Angel, E. Toy, S. Thomas and J. C. Glantz (2013). "Current role of robotic hysterectomy." *J Gynecol Surg* 29(4): 174-178.
12. Eklin, S., A. Lindfors, P. Sjöli and P. Dahm-Kahler (2015). "A prospective, comparative study on robotic versus open-surgery hysterectomy and pelvic lymphadenectomy for endometrial carcinoma." *Int J Gynecol Cancer* 25(2): 250-256.
13. Eoh, K. J., E. J. Nam, S. W. Kim, M. Shin, S. J. H. Kim, J. A. Kim and Y. T. Kim (2021). "Nationwide Comparison of Surgical and Oncologic Outcomes in Endometrial Cancer Patients Undergoing Robotic, Laparoscopic, and Open Surgery: A Population-Based Cohort Study." *Cancer Res Treat* 53(2): 549-557.
14. Esselen, K. M., A. Vitonis, J. Einarsson, M. G. Muto and S. Cohen (2015). "Health Care Disparities in Hysterectomy for Gynecologic Cancers: Data From the 2012 National Inpatient Sample." *Obstet Gynecol* 126(5): 1029-1039. [LOS is categorical]
15. Ferguson, S. E., T. Panzarella, S. Lau, L. T. Gien, V. Samouelian, C. Giede, H. Steed, T. Le, B. Renkosinski and M. Q. Bernardini (2018). "Prospective cohort study comparing quality of life and sexual health outcomes between women undergoing robotic, laparoscopic and open surgery for endometrial cancer." *Gynecol Oncol* 149(3): 476-483.
16. Gracia, M., J. García-Santos, M. Ramirez, M. Bellón, M. A. Herraiz and P. J. Coronado (2020). "Value of robotic surgery in endometrial cancer by body mass index." *Int J Gynaecol Obstet* 150(3): 398-405.
17. Gueli Alletti, S., E. Perrone, C. Fedele, S. Cianci, T. Pasciuto, V. Chiantera, S. Uccella, A. Ercoli, G. Vizzielli, A. Fagotti, V. Gallotta, F. Cosentino, B. Costantini, S. Restaino, G. Monterossi, A. Rosati, L. C. Turco, V. A. Capozzi, F. Fanfani and G. Scambia (2021). "A Multicentric Randomized Trial to Evaluate the ROle of Uterine MANipulator on Laparoscopic/Robotic HYsterectomy for the Treatment of Early-Stage Endometrial Cancer: The ROMANHY Trial." *Frontiers in Oncology* 11: 720894.
18. Guy, M. S., J. Sheeder, K. Behbakht, J. D. Wright and S. R. Guntupalli (2016). "Comparative outcomes in older and younger women undergoing laparotomy or robotic surgical staging for endometrial cancer." *Am J Obstet Gynecol* 214(3): 350 e351-350 e310.

Hysterectomy for endometrial cancer: bibliography (2 of 2)

December 31, 2022

19. Jung, Y. W., D. W. Lee, S. W. Kim, E. J. Nam, J. H. Kim, J. W. Kim and Y. T. Kim (2010). "Robot-assisted staging using three robotic arms for endometrial cancer: comparison to laparoscopy and laparotomy at a single institution." *J Surg Oncol* 101(2): 116-121.
20. Lundin, E. S., N. B. Wodlin, L. Nilsson and P. Kjolhede (2019). "A prospective randomized assessment of quality of life between open and robotic hysterectomy in early endometrial cancer." *Int J Gynecol Cancer*.
21. Maenpaa, M. M., K. Nieminen, E. I. Tomas, M. Laurila, T. H. Luukkaala and J. U. Maenpaa (2016). "Robotic-assisted vs traditional laparoscopic surgery for endometrial cancer: a randomized controlled trial." *Am J Obstet Gynecol* 215(5): 588 e581-588 e587.
22. Narducci, F., E. Bogart, T. Hebert, T. Gauthier, P. Collinet, J. M. Classe, F. Lecuru, A. Delest, S. Motton, V. Conri, C. Ferrer, F. Marchal, G. Ferron, A. Probst, J. Thery, M. C. Le Deley, D. Lefebvre, D. Francon, E. Leblanc and E. Lambaudie (2020). "Severe perioperative morbidity after robot-assisted versus conventional laparoscopy in gynecologic oncology: Results of the randomized ROBOGYN-1004 trial." *Gynecol Oncol* 158(2): 382-389.
23. Pilka, R., R. Marek, T. Adam, M. Kudela, D. Ondrova, D. Neubert, J. Hambalek, M. Maderka, D. Solichova, L. K. Krcmova and B. Melichar (2016). "Systemic Inflammatory Response After Open, Laparoscopic and Robotic Surgery in Endometrial Cancer Patients." *Anticancer Res* 36(6): 2909-2922.
24. Rozenholc, A., V. Samouelian, T. Warkus, P. Gauthier, D. Provencher, P. Sauthier, F. Gauthier, P. Drakopoulos and B. Cormier (2019). "Green versus blue: Randomized controlled trial comparing indocyanine green with methylene blue for sentinel lymph node detection in endometrial cancer." *Gynecol Oncol* 153(3): 500-504.
25. Safdieh, J., Y. C. Lee, A. Wong, A. Lee, J. P. Weiner, D. Schwartz and D. Schreiber (2017). "A Comparison of Outcomes Between Open Hysterectomy and Robotic-Assisted Hysterectomy for Endometrial Cancer Using the National Cancer Database." *Int J Gynecol Cancer* 27(7): 1508-1516.
26. Salehi, S., E. Avall-Lundqvist, B. Legerstam, J. W. Carlson and H. Falconer (2017). "Robot-assisted laparoscopy versus laparotomy for infrarenal paraaortic lymphadenectomy in women with high-risk endometrial cancer: A randomised controlled trial." *Eur J Cancer* 79: 81-89.
27. Silva, E. S. A., J. P. M. de Carvalho, C. Anton, R. P. Fernandes, E. C. Baracat and J. P. Carvalho (2018). "Introduction of robotic surgery for endometrial cancer into a Brazilian cancer service: a randomized trial evaluating perioperative clinical outcomes and costs." *Clinics (Sao Paulo)* 73(suppl 1): e522s.
28. Somashekhar, S. P., R. C. Jaka and S. S. Zaveri (2014). "Prospective Randomized Study Comparing Robotic-Assisted Hysterectomy and Regional Lymphadenectomy with Traditional Laparotomy for Staging of Endometrial Carcinoma -Initial Indian Experience." *Indian J Surg Oncol* 5(3): 217-223.
29. Wright, J. D., L. J. Havrilesky, D. E. Cohn, Y. Huang, J. Rathbun, L. W. Rice, C. L. Brown, R. D. Alvarez and E. M. Ko (2018). "Estimating potential for savings for low risk endometrial cancer using the Endometrial Cancer Alternative Payment Model (ECAP): A companion paper to the Society of Gynecologic Oncology Report on the Endometrial Cancer Alternative Payment Model." *Gynecol Oncol* 149(2): 241-247. [LOS is categorical]
30. Wright, J. D., W. M. Burke, A. I. Tergas, J. Y. Hou, Y. Huang, J. C. Hu, G. C. Hillyer, C. V. Ananth, A. I. Neugut and D. L. Hershman (2016). "Comparative Effectiveness of Minimally Invasive Hysterectomy for Endometrial Cancer." *J Clin Oncol* 34(10): 1087-1096.
31. Wright, J. D., W. M. Burke, E. T. Wilde, S. N. Lewin, A. S. Charles, J. H. Kim, N. Goldman, A. I. Neugut, T. J. Herzog and D. L. Hershman (2012). "Comparative effectiveness of robotic versus laparoscopic hysterectomy for endometrial cancer." *J Clin Oncol* 30(8): 783-791.
32. Yu, X., D. Lum, T. K. Kiet, K. C. Fuh, J. Orr, Jr., R. A. Brooks, S. M. Ueda, L. M. Chen, D. S. Kapp and J. K. Chan (2013). "Utilization of and charges for robotic versus laparoscopic versus open surgery for endometrial cancer." *J Surg Oncol* 107(6): 653-658.
33. Zakhari, A., N. Czuzoj-Shulman, A. R. Spence, W. H. Gotlieb and H. A. Abenhaim (2015). "Laparoscopic and robot-assisted hysterectomy for uterine cancer: a comparison of costs and complications." *Am J Obstet Gynecol* 213(5): 665 e661-667.
34. Zakhari, A., N. Czuzoj-Shulman, A. R. Spence, W. H. Gotlieb and H. A. Abenhaim (2016). "Hysterectomy for Uterine Cancer in the Elderly: A Comparison Between Laparoscopic and Robot-Assisted Techniques." *Int J Gynecol Cancer* 26(7): 1222-1227.

Disclosures

Important Safety Information

(US) Serious complications may occur in any surgery, including da Vinci surgery, up to and including death. Serious risks include, but are not limited to, injury to tissues and organs and conversion to other surgical techniques which could result in a longer operative time and/or increased complications. For summary of the risks associated with surgery refer to www.intuitive.com/safety.

Da Vinci Xi®/da Vinci X® system precaution statement

The demonstration of safety and effectiveness for the representative specific procedures did not include evaluation of outcomes related to the treatment of cancer (overall survival, disease-free survival, local recurrence), except for radical prostatectomy which was evaluated for overall survival, or treatment of the patient's underlying disease/condition. Device usage in all surgical procedures should be guided by the clinical judgment of an adequately trained surgeon.

(EU) Da Vinci X & Xi Surgical Systems

The Intuitive Surgical Endoscopic Instrument Control Systems (da Vinci X and da Vinci Xi Surgical Systems) are intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments during urologic surgical procedures, general laparoscopic surgical procedures, gynecologic laparoscopic surgical procedures, general thoracoscopic surgical procedures, and trans-oral otolaryngology surgical procedures restricted to benign tumors and malignant tumors classified as T1 and T2, and for benign base of tongue resection procedures. The systems are indicated for adult and pediatric use (except for trans-oral otolaryngology surgical procedures). They are intended to be used by trained physicians in an operating room environment.

The da Vinci X and da Vinci Xi Surgical Systems are class IIb medical devices CE marked (CE 2460) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

For product intended use and/or indications for use, risks, cautions, and warnings and full prescribing information, refer to the associated user manual(s) or visit <https://manuals.intuitivesurgical.com/market>. Some products, features or technologies may not be available in all countries. Product availability is subject to regulatory approval in the specific market. Please contact your local Intuitive representative for product availability in your region.

Individual outcomes may depend on a number of factors—including but not limited to—patient characteristics, disease characteristics, and/or surgeon experience.

Privacy Notice: Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2025 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos, including Intuitive, Da Vinci, and Ion, are trademarks or registered trademarks of Intuitive Surgical or their respective owner.

INTUITIVE

intuitive.com