

Evidence Navigator: Radical Prostatectomy

Systematic literature review summary
as of December 31, 2022

Purpose

The Evidence Navigator is a slide presentation representing a summary of the meta-analysis of the highest level of evidence available specific to a given procedure and published as of a particular date. It is created by the Global Evidence Management team within Global Access, Value and Economics (GAVE). It includes information that is available in the public domain. It is a systematic review and meta-analysis of the peer-reviewed literature based on a timeframe within which a literature search has been conducted according to a set of concise inclusion and exclusion criteria. The results of the meta-analysis are presented in the form of forest plots summarized for each outcome according to a comparator and surgical approach of interest. The summary results are reflective of a specific period in time and are subject to change with increasing literature. All of the robotic-assisted surgery procedures mentioned within the Evidence Navigator were performed using a da Vinci® surgical system.

Statistical analysis

All summary measures are shown as odds ratios, risk ratios or risk differences when describing binary outcomes, or as weighted mean differences or standardized mean differences when describing continuous outcomes. Weighting is based on the study sample size and variability of the outcome. A random effect model is used if heterogeneity is statistically significant, otherwise a fixed effect model is used. The Mantel Haenszel summary statistic is used for the overall results. The meta-analysis is performed with RevMan 5.4 (Review Manager, Version 5.4. Copenhagen: The Nordic Cochrane Centre, The Cochrane Collaboration, 2014) or R software (R Foundation for Statistical Computing, Vienna, Austria. URL <https://www.R-project.org/>).

Interpretation notes

When the effect size is measured as a standardized mean difference (SMD), or a risk difference (RD), it is not possible to provide a quantitative conclusion. In such cases, a qualitative conclusion is given with reference to its statistical significance. In some instances, studies may contain some overlapping patient populations. A redundancy check is performed in order to minimize this overlap and bias due to over-reporting.

Glossary

RAS	robotic-assisted surgery
Lap	laparoscopic surgery
LOE	level of evidence
HTA	health technology assessment
RCT	randomized controlled trial
OR	odds ratio
MD	mean difference
LNY	lymph node yield

WMD	weighted mean difference
RD	risk difference
SMD	standardized mean difference
95% CI	95% confidence interval
I²	test statistic for heterogeneity
EBL	estimated blood loss
LOS	length of hospital stay
PSM	positive surgical margins

Evidence Navigator: Radical Prostatectomy Summary Slides

Systematic literature review summary
as of December 31, 2022



WHAT DOES THE LITERATURE SHOW?

Systematic literature review & meta-analysis: Radical prostatectomy— clinical outcomes

Inclusion criteria

Robotic-assisted radical prostatectomy performed with a da Vinci® surgical system

January 1, 2010 – December 31, 2022

Level of Evidence = 1b, 2b, 2c

RCT, prospective cohort studies, or large database study (with $n \geq 20$ in each cohort)

Exclusion criteria

Not in English

Paper reports on a pediatric population

Publication is an HTA not published in a peer-reviewed journal

Alternate technique/approach (e.g. single-port)

No stratified analysis by study arm

Radical prostatectomy data mixed with other procedures

Original research study does not provide quantitative results

Original research publication includes redundant patient population and similar conclusions

95 publications including



Robotic-assisted patients: **669,924**




Laparoscopic patients: **59,442**



Open patients: **573,660**

Level of evidence





WHAT DOES THE LITERATURE SHOW?

Systematic literature review & meta-analysis key points: Robotic-assisted vs. laparoscopic radical prostatectomy



Favors robotic-assisted

- ↓ Blood transfusions by **44%**
- ↓ Length of stay by **0.5 days**
- ↓ 30-day postoperative complications by **28%**
- ↓ 30-day readmissions by **37%**
- ↓ 30-day mortality by **86%**
- ↓ Urinary incontinence at 6 months by **26%**
- ↓ Urinary incontinence at 12 months by **22%**
- ↓ Erectile dysfunction at 12 months by **45%**



Comparable outcomes

- ≈ Operative time
- ≈ Conversions
- ≈ Estimated blood loss
- ≈ Lymph node yield
- ≈ Positive surgical margins
- ≈ Intraoperative complications
- ≈ 30-day reoperations
- ≈ Erectile dysfunction at 6 months



Favors laparoscopic

None

Data collected through: December 31, 2022

■ Significant difference favoring
robotic-assisted surgery

■ No significant difference;
comparable outcomes

■ Significant difference favoring
laparoscopic surgery



WHAT DOES THE LITERATURE SHOW?

Systematic literature review & meta-analysis key points: Robotic-assisted vs. open radical prostatectomy



Favors robotic-assisted

- ↓ Estimated blood loss by **448 mL**
- ↓ Blood transfusions by **82%**
- ↑ Lymph node yield by **2 nodes**
- ↓ Positive surgical margins by **19%**
- ↓ Intraoperative complications by **42%**
- ↓ Length of stay by **1.4 days**
- ↓ 30-day postoperative complications by **47%**
- ↓ 30-day readmissions by **39%**
- ↓ 30-day mortality by **60%**
- ↓ Time to continence by **1.9 months**
- ↓ Urinary incontinence at 6 months by **22%**
- ↓ Erectile dysfunction at 6 months by **59%**
- ↓ Erectile dysfunction at 12 months by **34%**



Comparable outcomes

- ≈ 30-day reoperations
- ≈ Urinary incontinence at 12 months



Favors open

- ↓ Operative time by **49 minutes**

Data collected through: December 31, 2022

■ Significant difference favoring
robotic-assisted surgery

■ No significant difference;
comparable outcomes

■ Significant difference favoring
open surgery

Evidence Navigator: Radical Prostatectomy Technical Slides

Systematic literature review summary
as of December 31, 2022

radical prostatectomy: Literature search methods as of December 31, 2022

Monthly searches were conducted in PubMed, Scopus and Embase.

All citations were exported into a reference management system. Duplications were removed. Titles, abstracts and keywords were reviewed for literature review inclusion by Global Evidence Management team.

All robotic-assisted radical prostatectomies were performed with da Vinci® surgical systems. Publications were identified according to inclusion and exclusion criteria described.

Meta-analysis was performed using RevMan or R software.


95 publications

669,924 patients who underwent RAS

59,442 patients who underwent laparoscopic surgery

573,660 patients who underwent open surgery

Level of evidence



Criteria phase	Details
Identification phase	All unique PubMed, Scopus, and Embase references identified N=24,042 December 31, 2022
Inclusion criteria	
1. Robotic-assisted radical prostatectomy for cancer	Da Vinci® robotic-assisted radical prostatectomy N=6,767 (excluded N=17,275)
2. Year ≥ 2010	Articles published ≥ 2010 N=6,748 (excluded N=19)
3. LOE = 1b, 2b, 2c	Articles with LOE = 1b, 2b, 2c N=418 (excluded N=6,330)
4. Study is an RCT, prospective study or large database with comparative cohorts (robotic-assisted vs. laparoscopic and/or open surgery), with n≥20 in each cohort	Comparator cohorts N=394 (excluded N=24)
Exclusion criteria	N=299 excluded publications:
1. Not in English	N=1 (EC#1)
2. Paper reports on a pediatric population	N=0 (EC#2)
3. Publication is an HTA that was not published in a peer-reviewed journal	N=0 (EC#3)
4. Alternate technique/approach (e.g. single-port)	N=1 (EC#4)
5. No stratified analysis by study arm (e.g., combines results from robotic, laparoscopic and/or open cohorts)	N=169 (EC#5)
6. Radical prostatectomy data mixed with other procedures (e.g., data from multiple surgical procedures combined)	N=3 (EC#6)
7. Original research study does not provide quantitative results for at least one of the findings relative to the outcomes of interest (i.e., operative time, conversions, estimated blood loss and/or transfusions, complications, length of hospital stay, mortality, etc.)	N=106 (EC#7)
8. Original research publication includes redundant patient population and similar conclusions	N=19 (EC#8)

radical prostatectomy publications: N=95


Robotic-assisted vs. laparoscopic radical prostatectomy (1 of 3)

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic radical prostatectomy, the evidence for **robotic-assisted radical prostatectomy using the da Vinci surgical system** demonstrates:

- Significantly shorter hospital length of stay by an average of 0.5 days
- Comparable estimated blood loss
- Comparable operative time
- Comparable lymph node yield



Outcome	Robotic-assisted, n	Laparoscopic, n	Effect size 95% CI	P-value
Radical prostatectomy continuous variables (to December 31, 2022)				
LOS, days 1, 2, 6, 24, 27, 37, 42, 63, 79, 91, 92, 94, 95				
Subtotal	76129	12128	MD: -0.49 [-0.86, -0.12]	p<0.01
Random, Heterogeneity: p<0.01; I² = 99%				
EBL, mL 9, 63, 75, 77, 91, 94				
Subtotal	991	716	MD: -5.46 [-30.52, 19.60]	p=0.67
Random, Heterogeneity: p<0.01; I² = 79%				
Operative time, min 9, 24, 63, 75, 77, 79, 91, 94				
Subtotal	3507	3351	MD: -0.68 [-12.03, 10.68]	p=0.91
Random, Heterogeneity: p<0.01; I² = 94%				
LNY (Lap-RAS), n 2, 62				
Subtotal	3206	1182	MD: 0.70 [-1.80, 3.20]	p=0.58
L-R Random, Heterogeneity: p<0.01; I² = 98%				


Robotic-assisted vs. laparoscopic radical prostatectomy (2 of 3)

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic radical prostatectomy, the evidence for **robotic-assisted radical prostatectomy using the da Vinci surgical system** demonstrates:

- 86% lower likelihood of mortality within 30 days of surgery
- 45% less likely to experience erectile dysfunction at 12 months
- 44% less likely to receive a blood transfusion
- 37% less likely to experience a readmission within 30 days of surgery
- 28% less likely to experience a postoperative complication within 30 days of surgery
- 26% less likely to experience urinary incontinence at 6 months
- 22% less likely to experience urinary incontinence at 12 months



Outcome	Robotic-assisted, n	Laparoscopic, n	Effect size 95% CI	P-value
Radical prostatectomy binary variables (to December 31, 2022)				
Mortality, n(%)	1, 4, 37, 63, 76, 77, 95	5900	OR: 0.14 [0.04, 0.51]	p<0.01
Subtotal 45952				
Fixed, Heterogeneity: p=0.45; I²=0%				
Erectile dysfunction at 12 mo., n(%)	3, 4, 9, 24, 62, 78, 91, 92	827	OR: 0.55 [0.34, 0.90]	p=0.02
Subtotal 1945				
Random, Heterogeneity: p<0.01; I²=79%				
Transfusions, n(%)	2, 4, 6, 24, 28, 37, 39, 42, 76, 79, 92, 94	39975	OR: 0.56 [0.44, 0.71]	p<0.01
Subtotal 61493				
Random, Heterogeneity: p=0.03; I²=49%				
Readmissions, n(%)	37, 47, 60	8332	OR: 0.63 [0.47, 0.83]	p<0.01
Subtotal 39911				
Random, Heterogeneity: p<0.01; I²=84%				
Postoperative complications, n(%)	4, 5, 24, 26, 37, 42, 63, 76, 77, 79, 94, 95	8914	OR: 0.72 [0.53, 0.98]	p=0.04
Subtotal 34708				
Random, Heterogeneity: p<0.01; I²=80%				
Urinary incontinence at 6 mo., n(%)	3, 4, 62, 78, 91, 94	524	OR: 0.74 [0.58, 0.94]	p=0.02
Subtotal 1109				
Fixed, Heterogeneity: p=0.34; I²=12%				
Urinary incontinence at 12 mo., n(%)	3, 4, 9, 24, 62, 78, 91, 92	1122	OR: 0.78 [0.66, 0.92]	p<0.01
Subtotal 2428				
Random, Heterogeneity: p<0.01; I²=79%				


Robotic-assisted vs. laparoscopic radical prostatectomy (3 of 3)

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring laparoscopic surgery

Compared to laparoscopic radical prostatectomy, the evidence for **robotic-assisted radical prostatectomy using the da Vinci surgical system** demonstrates:

- Comparable erectile dysfunction rate at 6 months
- Comparable positive surgical margin rate
- Comparable reoperation rate within 30 days of surgery
- Comparable risk of conversion to open surgery
- Comparable risk of intraoperative complications



Outcome	Robotic-assisted, n	Laparoscopic, n	Effect size 95% CI	P-value
Radical prostatectomy binary variables (to December 31, 2022)				
Erectile dysfunction at 6 mo., n(%) 3, 4, 62, 78, 94				
Subtotal	695	288	OR: 0.47 [0.21, 1.05]	p=0.07
Random, Heterogeneity: p<0.01; I ² =82%				
Positive surgical margins, n(%) 4, 5, 9, 16, 24, 37, 39, 42, 63, 70, 73, 77, 84, 91, 94				
Subtotal	36514	7454	OR: 0.97 [0.81, 1.16]	p=0.72
Random, Heterogeneity: p<0.01; I ² =76%				
Reoperations, n(%) 24, 77				
Subtotal	952	351	OR: 1.17 [0.38, 3.62]	p=0.79
Fixed, Heterogeneity: p=0.4; I ² =0%				
Conversions, n(%) 4, 11, 26, 42, 77, 89, 94				
Subtotal	83886	6766	RD: -0.0112 [-0.0452, 0.022]	p=0.52
Random, Heterogeneity: p<0.01; I ² =98%				
Intraoperative complications, n(%) 42, 63				
Subtotal	382	918	RD: -0.0260 [-0.1061, 0.0541]	p=0.52
Random, Heterogeneity: p<0.01; I ² =88%				


Robotic-assisted vs. open radical prostatectomy (1 of 3)

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open radical prostatectomy, the evidence for **robotic-assisted radical prostatectomy using the da Vinci surgical system** demonstrates:

- Significantly less estimated blood loss by an average of 448 mL
- Significantly higher lymph node yield by an average of 2 nodes
- Significantly shorter time to continence by an average of 1.9 months
- Significantly shorter hospital length of stay by an average of 1.4 days
- Significantly longer operative time by an average of 49 minutes



Outcome	Robotic-assisted, n	Open, n	Effect size 95% CI	P-value
Radical prostatectomy continuous variables (to December 31, 2022)				
EBL, mL	7, 10, 15-17, 22, 29, 32, 34, 38, 41, 45, 46, 51, 55, 58, 61, 67, 86, 94			
Subtotal	10063	4649	MD: -448.12 [-550.93, -345.30]	p<0.01
Random, Heterogeneity: p<0.01; I² = 97%				
LNY (Open-RAS), n	2, 31, 34, 71, 82, 93			
Subtotal	4885	2650	MD: -2.01 [-3.05, -0.98]	p<0.01
O-R Random, Heterogeneity: p<0.01; I² = 88%				
Time to continence, months	45, 67			
Subtotal	207	107	MD: -1.90 [-3.66, -0.13]	p=0.03
Random, Heterogeneity: p<0.01; I² = 97%				
LOS, days	2, 6, 16, 19, 21, 23, 24, 27, 29, 34, 35, 37, 38, 42, 45, 46, 48, 52, 55, 58, 61, 64, 67, 79, 82, 86, 92, 93, 95			
Subtotal	168748	98289	MD: -1.41 [-1.81, -1.02]	p<0.01
Random, Heterogeneity: p=0; I² = 100%				
Operative time, min	7, 10, 15, 17, 19-21, 24, 29, 32, 34, 38, 45, 51, 58, 59, 61, 66, 71, 74, 79, 82, 83, 86, 93			
Subtotal	44109	46870	MD: 48.79 [31.76, 65.82]	p<0.01
Random, Heterogeneity: p=0; I² = 100%				


Robotic-assisted vs. open radical prostatectomy (2 of 3)

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open radical prostatectomy, the evidence for **robotic-assisted radical prostatectomy using the da Vinci surgical system** demonstrates:

- 82% less likely to receive a transfusion
- 60% lower likelihood of mortality within 30 days of surgery
- 59% less likely to experience erectile dysfunction at 6 months
- 47% less likely to experience a postoperative complication within 30 days of surgery
- 42% less likely to experience an intraoperative complication
- 39% less likely to experience a readmission within 30 days of surgery
- 34% less likely to experience erectile dysfunction at 12 months



Outcome	Robotic-assisted, n	Open, n	Effect size 95% CI	P-value
Radical prostatectomy binary variables (to December 31, 2022)				
Transfusions, n(%) 2, 6-8, 12, 13, 16, 17, 19-21, 23, 24, 28, 32, 34, 37-39, 41, 42, 58, 59, 64, 65, 67, 76, 79, 81, 85, 92, 93, 95				
Subtotal	194691	284388	OR: 0.18 [0.16, 0.21]	p<0.01
Random, Heterogeneity: p<0.01; I²=86%				
Mortality, n(%) 13, 19, 20, 28, 34, 37, 56-59, 61, 64, 65, 76, 81, 86, 88, 93, 95				
Subtotal	198152	293386	OR: 0.40 [0.31, 0.52]	p<0.01
Fixed, Heterogeneity: p=0.31; I²=13%				
Erectile dysfunction at 6 mo., n(%) 3, 18, 22, 40, 53, 80				
Subtotal	2817	2388	OR: 0.41 [0.28, 0.61]	p<0.01
Random, Heterogeneity: p<0.01; I²=79%				
Postoperative complications, n(%) 19, 20, 23, 24, 34, 37, 42, 45, 49, 55, 58, 59, 61, 64, 65, 67, 69, 71, 76, 79, 81, 93, 95				
Subtotal	209460	206754	OR: 0.53 [0.45, 0.62]	p<0.01
Random, Heterogeneity: p<0.01; I²=97%				
Intraoperative complications, n(%) 34, 42, 64, 65, 81, 93				
Subtotal	80389	60984	OR: 0.58 [0.41, 0.82]	p<0.01
Random, Heterogeneity: p<0.01; I²=80%				
Readmissions, n(%) 16, 23, 34, 37, 47, 57, 59-61, 67, 88, 93				
Subtotal	74199	36793	OR: 0.61 [0.50, 0.75]	p<0.01
Random, Heterogeneity: p<0.01; I²=85%				
Erectile dysfunction at 12 mo., n(%) 3, 18, 20, 22, 24, 29, 40, 45, 53, 54, 68, 92				
Subtotal	8613	5469	OR: 0.66 [0.49, 0.89]	p<0.01
Random, Heterogeneity: p<0.01; I²=90%				


Robotic-assisted vs. open radical prostatectomy (3 of 3)

Summary as of December 31, 2022

■ Significant difference favoring robotic-assisted surgery ■ No significant difference; comparable outcomes ■ Significant difference favoring open surgery

Compared to open radical prostatectomy, the evidence for **robotic-assisted radical prostatectomy using the da Vinci surgical system** demonstrates:

- 22% less likely to experience urinary incontinence at 6 months
- 19% lower likelihood of a positive surgical margin
- Comparable urinary incontinence rate at 12 months
- Comparable reoperation rate within 30 days of surgery



Outcome	Robotic-assisted, n	Open, n	Effect size 95% CI	P-value
Radical prostatectomy binary variables (to December 31, 2022)				
Urinary incontinence at 6 mo., n(%)	3, 18, 22, 25, 51, 53, 67, 80			
Subtotal	2871	2668	OR: 0.78 [0.62, 0.98]	p=0.04
Random, Heterogeneity: p=0.03; I²=54%				
Positive surgical margins, n(%)	2, 7, 12, 15, 16, 20, 22, 24, 25, 29-31, 33, 37, 39-45, 50, 51, 53, 54, 57, 58, 61, 67-71, 84, 88, 93			
Subtotal	88319	52921	OR: 0.81 [0.73, 0.90]	p<0.01
Random, Heterogeneity: p<0.01; I²=90%				
Urinary incontinence at 12 mo., n(%)	3, 18, 20, 22, 24, 25, 29, 30, 45, 51, 53, 54, 67, 68, 92			
Subtotal	10252	6411	OR: 0.97 [0.90, 1.05]	p=0.43
Fixed, Heterogeneity: p=0.06; I²=39%				
Reoperations, n(%)	13, 20, 24, 28, 58, 59, 61, 67, 86			
Subtotal	33210	151636	OR: 0.99 [0.87, 1.13]	p=0.9
Fixed, Heterogeneity: p=0.39; I²=6%				

Radical prostatectomy bibliography (1 of 5)

December 31, 2022

1. Anderson, J. E., D. C. Chang, J. K. Parsons and M. A. Talamini (2012). "The first national examination of outcomes and trends in robotic surgery in the United States." *J Am Coll Surg* 215(1): 107-114; discussion 114-106.
2. Aning, J. J., G. S. Reilly, S. Fowler, B. Challacombe, J. S. McGrath, P. Sooriakumaran and B. S. o. Oncology (2019). "Perioperative and oncological outcomes of radical prostatectomy for high-risk prostate cancer in the UK: an analysis of surgeon-reported data." *BJU Int* 124(3): 441-448.
3. Antonelli, A., C. Palumbo, M. Noale, A. Porreca, S. Maggi, C. Simeone, P. Bassi, F. Bertoni, S. Bracarda, M. Buglione, G. N. Conti, R. Corvo, M. Gacci, V. Mirone, R. Montironi, L. Triggiani, A. Tubaro, W. Artibani and I. T. C. N. R. s. g. Pros (2019). "Impact of Surgical Approach on Patient-Reported Outcomes after Radical prostatectomy: A Propensity Score-Weighted Analysis from a Multicenter, Prospective, Observational Study (The Pros-IT CNR Study)." *Urol Int* 103(1): 8-18.
4. Asimakopoulos, A. D., C. T. Pereira Fraga, F. Annino, P. Pasqualetti, A. A. Calado and C. Mugnier (2011). "Randomized comparison between laparoscopic and robot-assisted nerve-sparing radical prostatectomy." *J Sex Med* 8(5): 1503-1512.
5. Asimakopoulos, A. D., R. Miano, N. Di Lorenzo, E. Spera, G. Vespasiani and C. Mugnier (2013). "Laparoscopic versus robot-assisted bilateral nerve-sparing radical prostatectomy: comparison of pentalecta rates for a single surgeon." *Surg Endosc* 27(11): 4297-4304.
6. Basto, M., N. Sathianathan, L. Te Marvelde, S. Ryan, J. Goad, N. Lawrentschuk, A. J. Costello, D. A. Moon, A. G. Heriot, J. Butler and D. G. Murphy (2016). "Patterns-of-care and health economic analysis of robot-assisted radical prostatectomy in the Australian public health system." *BJU Int* 117(6): 930-939.
7. Beck, S., L. Zins, C. Holthausen, C. Rademacher, F. von Breunig, P. Tennstedt, A. Haese, M. Graefen, C. Zöllner and M. Fischer (2020). "Comparison of Cognitive Function After Robot-Assisted radical prostatectomy and Open Retropubic Radical prostatectomy: A Prospective Observational Single-Center Study." *Urology* 139: 110-117.
8. Beilstein, C. M., M. Huber, M. A. Furrer, L. M. Löffel, P. Y. Wuethrich and D. Engel (2022). "Impact of analgesic techniques on early quality of recovery after radical prostatectomy: A 3-arm, randomized trial." *European Journal of Pain (United Kingdom)* 26(9): 1990-2002.
9. Berge, V., R. E. Berg, J. R. Hoff, N. Wessel, L. M. Diep, S. J. Karlsen and L. M. Eri (2013). "A prospective study of transition from laparoscopic to robot-assisted radical prostatectomy: quality of life outcomes after 36-month follow-up." *Urology* 81(4): 781-786.
10. Bier, S., J. Hennenlotter, S. Rausch, S. Aufderklamm, J. C. Martzog, A. Stenzl, C. Schwentner and T. Todenhofer (2016). "Return to Work and Normal Daily Life Activity after Open and Robot-Assisted Radical prostatectomy--A Single Surgeon Analysis." *Urol Int* 96(3): 280-286.
11. Binhas, M., L. Salomon, F. Roudot-Thoraval, C. Armand, B. Plaud and J. Marty (2012). "Radical prostatectomy with robot-assisted radical prostatectomy and laparoscopic radical prostatectomy under low-dose aspirin does not significantly increase blood loss." *Urology* 79(3): 591-595.
12. Breyer, B. N., C. B. Davis, J. E. Cowan, C. J. Kane and P. R. Carroll (2010). "Incidence of bladder neck contracture after robot-assisted laparoscopic and open radical prostatectomy." *BJU Int* 106(11): 1734-1738.
13. Carlsson, S., A. E. Nilsson, M. C. Schumacher, M. N. Jonsson, D. S. Volz, G. Steineck and P. N. Wiklund (2010). "Surgery-related complications in 1253 robot-assisted and 485 open retropubic radical prostatectomies at the Karolinska University Hospital, Sweden." *Urology* 75(5): 1092-1097.
14. Carter, S. C., S. Lipsitz, Y. C. Shih, P. L. Nguyen, Q. D. Trinh and J. C. Hu (2014). "Population-Based Determinants of Radical prostatectomy Operative Time." *BJU International* 113(5b): E112-118.
15. Cazzaniga, W., R. A. Godtman, S. Carlsson, G. Ahlgren, E. Johansson, D. Robinson, J. Hugosson and P. Stattin (2019). "Population-based, nationwide registration of prostatectomies in Sweden." *Journal of Surgical Oncology* 120(4): 803-812.
16. Chang, P., A. A. Wagner, M. M. Regan, J. A. Smith, C. S. Saigal, M. S. Litwin, J. C. Hu, M. R. Cooperberg, P. R. Carroll, E. A. Klein, A. S. Kibel, G. L. Andriole, M. Han, A. W. Partin, D. P. Wood, C. M. Crociani, T. K. Greenfield, D. Patil, L. A. Hembroff, K. Davis, L. Stork, D. E. Spratt, J. T. Wei and M. G. Sanda (2022). "Prospective Multicenter Comparison of Open and Robotic Radical prostatectomy: The PROST-QA/RP2 Consortium." *The Journal of urology* 207(1): 127-136.
17. Cheetham, P. J., D. J. Lee, A. Rose-Morris, S. F. Brewster and K. K. Badani (2010). "Does the presence of robotic surgery affect demographics in patients choosing to undergo radical prostatectomy? A multi-center contemporary analysis." *Journal of Robotic Surgery* 4: 155-160.
18. Coughlin, G. D., J. W. Yaxley, S. K. Chambers, S. Occhipinti, H. Samaratunga, L. Zajdlewicz, P. Teloken, N. Dunlison, S. Williams, M. F. Lavin and R. A. Gardiner (2018). "Robot-assisted laparoscopic radical prostatectomy versus open radical retropubic radical prostatectomy: 24-month outcomes from a randomised controlled study." *Lancet Oncol* 19(8): 1051-1060.
19. Davis, J. W., U. S. Kreaden, J. Gabbert and R. Thomas (2014). "Learning curve assessment of robot-assisted radical prostatectomy compared with open-surgery controls from the premier perspective database." *J Endourol* 28(5): 560-566.
20. Di Pierro, G. B., P. Baumeister, P. Stucki, J. Beatrice, H. Danuser and A. Mattei (2011). "A prospective trial comparing consecutive series of open retropubic and robot-assisted laparoscopic radical prostatectomy in a centre with a limited caseload." *Eur Urol* 59(1): 1-6.

Radical prostatectomy bibliography (2 of 5)

December 31, 2022

21. Doumerc, N., C. Yuen, R. Savdie, M. B. Rahman, K. K. Rasiah, R. Pe Benito, W. Delprado, J. Matthews, A. M. Haynes and P. D. Stricker (2010). "Should experienced open prostatic surgeons convert to robotic surgery? The real learning curve for one surgeon over 3 years." *BJU Int* 106(3): 378-384.
22. Fode, M., J. Sonksen and H. Jakobsen (2014). "Radical radical prostatectomy: initial experience with robot-assisted laparoscopic procedures at a large university hospital." *Scand J Urol* 48(3): 252-258.
23. Gandaglia, G., J. D. Sammon, S. L. Chang, T. K. Choueiri, J. C. Hu, P. I. Karakiewicz, A. S. Kibel, S. P. Kim, R. Konijeti, F. Montorsi, P. L. Nguyen, S. Sukumar, M. Menon, M. Sun and Q. D. Trinh (2014). "Comparative effectiveness of robot-assisted and open radical radical prostatectomy in the postdissemination era." *J Clin Oncol* 32(14): 1419-1426.
24. García Cortés, Á., J. Colombás Vives, C. Gutiérrez Castañé, S. Chiva San Román, P. Doménech López, F. J. Ancizu Marckert, M. Hevia Suárez, I. Merino Narro, J. M. Velis Campillo, F. Guillén Grima, M. Torres Roca, F. Díez-Caballero y Alonso, D. Rosell Costa, F. Villacampa Aubá, F. R. de Fata Chillón, G. Andrés Boville, G. Barbas Bernardos, B. Miñana López, J. E. Robles García and J. I. Pascual Piédrola (2022). "Comparison of surgical approaches to radical radical prostatectomy in our series beyond oncological and functional outcomes." *Actas Urológicas Españolas* 46(5): 275-284.
25. Geraerts, I., H. Van Poppel, N. Devoogdt, B. Van Cleynenbreugel, S. Joniau and M. Van Kampen (2013). "Prospective evaluation of urinary incontinence, voiding symptoms and quality of life after open and robot-assisted radical radical prostatectomy." *BJU Int* 112(7): 936-943.
26. Gozen, A. S., Y. Akin, M. Akgul, C. Yazici, J. Klein and J. Rassweiler (2014). "A Novel Practical Trocar Placement Technique for Extraperitoneal Laparoscopic and Robotic-Assisted Laparoscopic Radical radical prostatectomy in Patients with Lower Midline Abdominal Incisions." *J Laparoendosc Adv Surg Tech A* 24(6): 417-421.
27. Gray, W. K., J. Day, T. W. R. Briggs and S. Harrison (2022). "An observational study of volume-outcome effects for robot-assisted radical radical prostatectomy in England." *BJU international* 129(1): 93-103.
28. Groeben, C., R. Koch, M. Baunacke, M. P. Wirth and J. Huber (2017). "High volume is the key for improving in-hospital outcomes after radical radical prostatectomy: a total population analysis in Germany from 2006 to 2013." *World J Urol* 35(7): 1045-1053.
29. Haglind, E., S. Carlsson, J. Stranne, A. Wallerstedt, U. Wilderang, T. Thorsteinsdottir, M. Lagerkvist, J. E. Damber, A. Bjartell, J. Hugosson, P. Wiklund, G. Steineck and L. s. committee (2015). "Urinary Incontinence and Erectile Dysfunction After Robotic Versus Open Radical radical prostatectomy: A Prospective, Controlled, Nonrandomised Trial." *Eur Urol* 68(2): 216-225.
30. Hagman, A., A. Lantz, S. Carlsson, J. Höijer, J. Stranne, S. I. Tyritzis, E. Haglind, A. Bjartell, J. Hugosson, O. Akre, G. Steineck and P. Wiklund (2021). "Urinary continence recovery and oncological outcomes after surgery for prostate cancer analysed by risk category: results from the LAParoscopic radical prostatectomy robot and open trial." *World Journal of Urology* 39(9): 3239-3249.
31. Herlemann, A., J. E. Cowan, P. R. Carroll and M. R. Cooperberg (2018). "Community-based Outcomes of Open versus Robot-assisted Radical radical prostatectomy." *Eur Urol* 73(2): 215-223.
32. Hong, J. Y., J. Y. Kim, Y. D. Choi, K. H. Rha, S. J. Yoon and H. K. Kil (2010). "Incidence of venous gas embolism during robotic-assisted laparoscopic radical radical prostatectomy is lower than that during radical retropubic radical prostatectomy." *Br J Anaesth* 105(6): 777-781.
33. Hu, J. C., G. Gandaglia, P. I. Karakiewicz, P. L. Nguyen, Q.-D. Trinh, Y.-C. T. Shih, F. Abdollah, K. Chamie, J. L. Wright and P. A. Ganz (2014). "Comparative Effectiveness of Robot-assisted Versus Open Radical radical prostatectomy Cancer Control." *European Urology* 66: 666-672.
34. Huang, M. M., Z. T. Su, R. E. N. Becker, C. P. Pavlovich, A. W. Partin, M. E. Allaf and H. D. Patel (2021). "Complications After Open and Robotic-Assisted Radical radical prostatectomy and Association with Post-Operative Opioid Use: An Analysis of Data from the PREVENTER Trial." *BJU international* 127(2): 190-197.
35. Hyams, E. S., J. K. Mullins, P. M. Pierorazio, A. W. Partin, M. E. Allaf and B. R. Matlaga (2013). "Impact of robotic technique and surgical volume on the cost of radical radical prostatectomy." *J Endourol* 27(3): 298-303.
36. Inoue, S., T. Hayashi, J. Teishima and A. Matsubara (2020). "Comparative analysis of robot-assisted, laparoscopic, and open radical radical prostatectomy regarding lower urinary tract symptoms: A longitudinal study." *Urological Science* 31(1): 21-27.
37. John, J. B., J. Pascoe, S. Fowler, T. Walton, M. Johnson, J. Aning, B. Challacombe, R. Bufacchi, A. J. Dickinson and J. S. McGrath (2022). "A 'real-world' standard for radical radical prostatectomy: Analysis of the British Association of Urological Surgeons Complex Operations Reports, 2016–2018." *Journal of Clinical Urology* 16(5): 432-443.
38. Karaveli, A., A. S. Kavakli, O. Cakin, G. Aykal, A. Yildiz and M. Ates (2022). "Comparison of plasma neutrophil gelatinase-associated lipocalin (NGAL) levels after robot-assisted laparoscopic and retropubic radical radical prostatectomy: an observational study." *Braz J Anesthesiol* 72(1): 21-28.
39. Khadhour, S., C. Miller, S. Fowler, L. Hounsborne, A. McNeill, J. Adshead, J. S. McGrath and B. S. o. Oncology (2018). "The British Association of Urological Surgeons (BAUS) radical radical prostatectomy audit 2014/2015 - an update on current practice and outcomes by centre and surgeon case-volume." *BJU Int* 121(6): 886-892.

Radical prostatectomy bibliography (3 of 5)

December 31, 2022

40. Kim, S. C. S., C.; Kim, W.; Kang, T.; Park, J.; Jeong, I. G.; Lee, S.; Cho, Y. M.; Ahn, H. (2011). "Factors Determining Functional Outcomes After Radical prostatectomy: Robot-Assisted Versus Retropubic." *European Urology* 60(3): 413-419.
41. Kordan, Y., D. A. Barocas, H. O. Altamar, P. E. Clark, S. S. Chang, R. Davis, S. D. Herrell, R. Baumgartner, V. Mishra, R. C. Chan, J. A. Smith, Jr. and M. S. Cookson (2010). "Comparison of transfusion requirements between open and robotic-assisted laparoscopic radical prostatectomy." *BJU Int* 106(7): 1036-1040.
42. Laird, A., S. Fowler, D. W. Good, G. D. Stewart, V. Srinivasan, D. Cahill, S. F. Brewster, S. A. McNeill and S. British Association of Urological (2015). "Contemporary practice and technique-related outcomes for radical prostatectomy in the UK: a report of national outcomes." *BJU Int* 115(5): 753-763.
43. Lantz, A., D. Bock, O. Akre, E. Angenete, A. Bjartell, S. Carlsson, K. K. Modig, M. Nyberg, K. S. Kollberg, G. Steineck, J. Stranne, P. Wiklund and E. Haglind (2021). "Functional and Oncological Outcomes After Open Versus Robot-assisted Laparoscopic Radical prostatectomy for Localised Prostate Cancer: 8-Year Follow-up." *European Urology* 80(5): 650-660.
44. Loeb, S., J. I. Epstein, A. E. Ross, L. Schultz, E. B. Humphreys and J. P. Jarow (2010). "Benign prostate glands at the bladder neck margin in robotic vs open radical prostatectomy." *BJU Int* 105(10): 1446-1449.
45. Ludovico, G. M., G. Dacheille, G. Pagliarulo, C. D'Elia, N. Mondaini, M. Gacci, B. Detti, G. Malossini, R. Bartoletti and T. Cai (2013). "Bilateral nerve sparing robotic-assisted radical prostatectomy is associated with faster continence recovery but not with erectile function recovery compared with retropubic open radical prostatectomy: the need for accurate selection of patients." *Oncol Rep* 29(6): 2445-2450.
46. Martinschek, A., L. Stumm, M. Ritter, E. Heinrich, C. Bolenz and L. Trojan (2017). "Prospective, Controlled Study of Invasiveness and Post-Aggression Metabolism in Patients Undergoing Robotic-Assisted Radical prostatectomy." *Urol Int* 99(2): 201-206.
47. Maynou, L., W. T. Mehtsun, V. Serra-Sastre and I. Papanicolas (2021). "Patterns of adoption of robotic radical prostatectomy in the United States and England." *Health Services Research* 56 Suppl 3: 1441-1461.
48. Nabi, J., D. F. Friedlander, X. Chen, A. P. Cole, J. C. Hu, A. S. Kibel, P. Dasgupta and Q. D. Trinh (2020). "Assessment of Out-of-Pocket Costs for Robotic Cancer Surgery in US Adults." *JAMA network open* 3(1): e1919185.
49. Ning, C., X. Hu, F. Liu, J. Lin, J. Zhang, Z. Wang and Y. Zhu (2019). "Post-surgical outcomes of patients with chronic kidney disease and end stage renal disease undergoing radical prostatectomy: 10-year results from the US National Inpatient Sample." *BMC Nephrology* 20(1): 278.
50. Nyberg, M., O. Akre, D. Bock, S. V. Carlsson, S. Carlsson, J. Hugosson, A. Lantz, G. Steineck, J. Stranne, S. Tyritzis, P. Wiklund, E. Haglind and A. Bjartell (2020). "Risk of Recurrent Disease 6 Years After Open or Robotic-assisted Radical prostatectomy in the Prospective Controlled Trial LAPPRO." *European Urology Open Science* 20: 54-61.
51. Ohara, E., N. Kawamorita, Y. Satake, Y. Kaiho, K. Mitsuzuka, H. Saito, S. Ishidoya, Y. Arai and A. Ito (2022). "Minimal residual membranous urethral length and membranous urethral length predict poor recovery from incontinence after robot-assisted radical prostatectomy and after open radical prostatectomy." *International Journal of Urology* 29(12): 1517-1523.
52. Okhawere, K. E., I. F. Shih, S. H. Lee, Y. Li, J. A. Wong and K. K. Badani (2021). "Comparison of 1-Year Health Care Costs and Use Associated with Open vs Robotic-Assisted Radical prostatectomy." *JAMA Network Open* 4(3): e212265.
53. O'Neil, B., T. Koyama, J. Alvarez, R. M. Conwill, P. C. Albertsen, M. R. Cooperberg, M. Goodman, S. Greenfield, A. S. Hamilton, K. E. Hoffman, R. M. Hoffman, S. H. Kaplan, J. L. Stanford, A. M. Stroup, L. E. Paddock, X. C. Wu, R. A. Stephenson, M. J. Resnick, D. A. Barocas and D. F. Penson (2016). "The Comparative Harms of Open and Robotic radical prostatectomy in Population Based Samples." *J Urol* 195(2): 321-329.
54. Ong, W. L., S. M. Evans, T. Spelman, P. A. Kearns, D. G. Murphy and J. L. Millar (2016). "Comparison of oncological and health-related quality of life outcomes between open and robot-assisted radical prostatectomy for localised prostate cancer - findings from the population-based Victorian Prostate Cancer Registry." *BJU Int* 118(4): 563-569.
55. Patel, H. D., A. Srivastava, N. D. Patel, F. A. Faisal, W. Ludwig, G. A. Joice, Z. R. Schwen, M. E. Allaf, M. Han and A. S. Herati (2019). "A Prospective Cohort Study of Postdischarge Opioid Practices After Radical prostatectomy: The ORIOLES Initiative." *Eur Urol* 75(2): 215-218.
56. Patel, H. D., F. A. Faisal, N. D. Patel, C. P. Pavlovich, M. E. Allaf, M. Han and A. S. Herati (2020). "Effect of a prospective opioid reduction intervention on opioid prescribing and use after radical prostatectomy: results of the Opioid Reduction Intervention for Open, Laparoscopic, and Endoscopic Surgery (ORIOLES) Initiative." *BJU Int* 125(3): 426-432.
57. Pearce, S. M., J. J. Pariser, T. Karrison, S. G. Patel and S. E. Eggener (2016). "Comparison of Perioperative and Early Oncologic Outcomes between Open and Robotic Assisted Laparoscopic radical prostatectomy in a Contemporary Population Based Cohort." *J Urol* 196(1): 76-81.
58. Philippou, P., E. Waane and E. Rowe (2012). "Robot-assisted laparoscopic radical prostatectomy versus open: comparison of the learning curve of a single surgeon." *J Endourol* 26(8): 1002-1008.
59. Pilecki, M. A., B. B. McGuire, U. Jain, J. Y. Kim and R. B. Nadler (2014). "National multi-institutional comparison of 30-day postoperative complication and readmission rates between open retropubic radical prostatectomy and robot-assisted laparoscopic radical prostatectomy using NSQIP." *J Endourol* 28(4): 430-436.

Radical prostatectomy bibliography (4 of 5)

December 31, 2022

60. Ploussard, G., A. Grabia, E. Barret, J. B. Beauval, L. Brureau, G. Créhange, C. Dariane, G. Fiard, G. Fromont, M. Gauthé, R. Mathieu, R. Renard-Penna, G. Roubaud, A. Ruffion, P. Sargos, M. Rouprêt and C. E. Lequeu (2022). "Annual nationwide analysis of costs and post-operative outcomes after radical radical prostatectomy according to the surgical approach (open, laparoscopic, and robotic)." *World J Urol* 40(2): 419-425.
61. Porcaro, A. B., N. De Luyk, P. Corsi, M. Sebben, A. Tafuri, D. Inverardi, D. De Marchi, I. Tamanini, M. Brunelli, M. A. Cerruto, G. L. Salvagno, G. C. Guidi and W. Artibani (2016). "Robotic assisted radical prostatectomy accelerates postoperative stress recovery: Final results of a contemporary prospective study assessing pathophysiology of cortisol peri-operative kinetics in prostate cancer surgery." *Asian Journal of Urology* 3(2): 88-95.
62. Porpiglia, F., C. Fiori, R. Bertolo, M. Manfredi, F. Mele, E. Checcucci, S. De Luca, R. Passera and R. M. Scarpa (2018). "Five-year Outcomes for a Prospective Randomised Controlled Trial Comparing Laparoscopic and Robot-assisted Radical prostatectomy." *Eur Urol Focus* 4(1): 80-86.
63. Porpiglia, F., I. Morra, M. Lucci Chiarissi, M. Manfredi, F. Mele, S. Grande, F. Ragni, M. Poggio and C. Fiori (2013). "Randomised controlled trial comparing laparoscopic and robot-assisted radical prostatectomy." *Eur Urol* 63(4): 606-614.
64. Preisser, F., S. Nazzani, E. Mazzone, S. Knipper, M. Bandini, Z. Tian, A. Haese, F. Saad, K. C. Zorn, F. Montorsi, S. F. Shariat, M. Graefen, D. Tilki and P. I. Karakiewicz (2019). "Regional differences in total hospital charges between open and robotically assisted radical prostatectomy in the United States." *World J Urol* 37(7): 1305-1313.
65. Qin, Y., H. Han, Y. Xue, C. Wu, X. Wei, Y. Liu, Y. Cao, Y. Ruan and J. He (2020). "Comparison and trend of perioperative outcomes between robot-assisted radical prostatectomy and open radical prostatectomy: nationwide inpatient sample 2009-2014." *International braz j urol : official journal of the Brazilian Society of Urology* 46(5): 754-771.
66. Quinto, D., S. T. Reis, L. J. Zampolli, R. Pimenta, V. R. Guimarães, N. I. Viana, G. A. dos Santos, M. P. Gimenez, K. R. Leite, H. Zampolli, J. A. S. da Cruz, M. Srougi and C. C. Passerotti (2021). "Robotically assisted laparoscopic radical prostatectomy induces lower tissue trauma than radical retropubic radical prostatectomy." *Journal of Robotic Surgery* 15(1): 147-151.
67. Raheem, A. A., A. Hagra, A. Ghaith, M. J. Alenzi, A. Elghiaty, T. Gameel, I. Alowidah, W. S. Ham, Y. D. Choi, A. H. El-Bahnasy, A. Omar, M. El-Bendary and K. H. Rha (2020). "Retzius-sparing robot-assisted radical prostatectomy versus open retropubic radical prostatectomy: a prospective comparative study with 19-month follow-up." *Minerva Urol Nefrol* 72(5): 586-594.
68. Rechtman, M., A. Forbes, J. L. Millar, M. Evans, L. Dodds, D. G. Murphy and S. M. Evans (2022). "Comparison of urinary and sexual patient-reported outcomes between open radical prostatectomy and robot-assisted radical prostatectomy: a propensity score matched, population-based study in Victoria." *BMC Urol* 22(1): 18.
69. Reisz, P. A., A. A. Laviana, Z. Zhao, L. C. Huang, T. Koyama, R. Conwill, K. Hoffman, M. Goodman, A. S. Hamilton, X. C. Wu, L. E. Paddock, A. Stroup, M. R. Cooperberg, M. Hashibe, B. B. O'Neil, S. H. Kaplan, S. Greenfield, D. F. Penson and D. A. Barocas (2020). "Assessing the Quality of Surgical Care for Clinically Localized Prostate Cancer: Results from the CEASAR Study." *J Urol* 204(6): 1236-1241.
70. Retel, V. P., C. Bouchardy, M. Usel, I. Neyroud-Caspar, F. Schmidlin, G. Wirth, C. Iselin, R. Miralbell and E. Rapiti (2014). "Determinants and effects of positive surgical margins after radical prostatectomy on prostate cancer mortality: a population-based study." *BMC Urol* 14(1): 86.
71. Roscigno, M., G. La Croce, R. Naspro, M. Nicolai, M. Manica, M. Scarcello, D. Chinaglia and L. F. Da Pozzo (2019). "Extended pelvic lymph node dissection during radical prostatectomy: comparison between initial robotic experience of a high volume open surgeon and his contemporary open series." *Minerva Urol Nefrol* 71(6): 597-604.
72. Salciccia, S., D. Rosati, P. Viscuso, V. Canale, E. Scarrone, M. Frisenda, R. Catuzzi, M. Moriconi, V. Asero, S. Signore, M. De Dominicis, P. Emiliozzi, A. Carbone, A. L. Pastore, A. Fuschi, G. B. Di Pierro, A. Gentilucci, S. Cattarino, G. Mariotti, G. M. Busetto, M. Ferro, E. De Berardinis, G. P. Ricciuti, V. Panebianco, F. M. Magliocca, F. Del Giudice, M. Maggi and A. Sciarra (2021). "Influence of operative time and blood loss on surgical margins and functional outcomes for laparoscopic versus robotic-assisted radical prostatectomy: a prospective analysis." *Cent European J Urol* 74(4): 503-515.
73. Sciarra, A., M. Frisenda, M. Maggi, F. M. Magliocca, A. Ciardi, V. Panebianco, E. Berardinis, S. Salciccia, G. B. Di Pierro, A. Gentilucci, F. Del Giudice, G. M. Busetto and A. Tufano (2021). "Prospective comparative trial on nerve-sparing radical prostatectomy using a robot-assisted versus laparoscopic technique: expectation versus satisfaction and impact on surgical margins." *Cent European J Urol* 74(2): 169-177.
74. Simon, R. M., L. E. Howard, D. M. Moreira, M. K. Terris, C. J. Kane, W. J. Aronson, C. L. Amling, M. R. Cooperberg and S. J. Freedland (2017). "Predictors of operative time during radical retropubic radical prostatectomy and robot-assisted laparoscopic radical prostatectomy." *Int J Urol* 24(8): 618-623.
75. Sofra, M., A. Antenucci, M. Gallucci, C. Mandoj, R. Papalia, C. Claroni, I. Monteferrante, G. Torregiani, V. Gianaroli, I. Sperduti, L. Tomao and E. Forastiere (2014). "Perioperative changes in pro and anticoagulant factors in prostate cancer patients undergoing laparoscopic and robotic radical prostatectomy with different anaesthetic techniques." *J Exp Clin Cancer Res* 33(1): 63.
76. Stolzenburg, J. U., I. Kyriazis, C. Fahlenbrach, C. Gilfrich, C. Gunster, E. Jeschke, G. Popken, L. Weissbach, C. von Zastrow and H. Leicht (2016). "National trends and differences in morbidity among surgical approaches for radical prostatectomy in Germany." *World J Urol* 34(11): 1515-1520.

Radical prostatectomy bibliography (5 of 5)

December 31, 2022

77. Stolzenburg, J. U., S. Holze, P. Neuhaus, I. Kyriazis, H. M. Do, A. Dietel, M. C. Truss, C. I. Grzella, D. Teber, M. Hohenfellner, R. Rabenalt, P. Albers and M. Mende (2021). "Robotic-assisted Versus Laparoscopic Surgery: Outcomes from the First Multicentre, Randomised, Patient-blinded Controlled Trial in Radical radical prostatectomy (LAP-01)." *Eur Urol* 79(6): 750-759.
78. Stolzenburg, J. U., S. Holze, V. K. A. Arthanareeswaran, P. Neuhaus, H. M. Do, C. M. Haney, A. Dietel, M. C. Truss, K. D. Stützel, D. Teber, M. Hohenfellner, R. Rabenalt, P. Albers and M. Mende (2022). "Robotic-assisted Versus Laparoscopic Radical radical prostatectomy: 12-month Outcomes of the Multicentre Randomised Controlled LAP-01 Trial." *European Urology Focus* 8(6): 1583-1590.
79. Sugihara, T., H. Yasunaga, H. Horiguchi, H. Matsui, T. Fujimura, H. Nishimatsu, H. Fukuhara, H. Kume, Y. Changhong, M. W. Kattan, K. Fushimi and Y. Homma (2014). "Robot-assisted versus other types of radical radical prostatectomy: population-based safety and cost comparison in Japan, 2012-2013." *Cancer Sci* 105(11): 1421-1426.
80. Thompson, J. E., S. Egger, M. Bohm, A. R. Siriwardana, A. M. Haynes, J. Matthews, M. J. Scheltema and P. D. Stricker (2018). "Superior Biochemical Recurrence and Long-term Quality-of-life Outcomes Are Achievable with Robotic Radical radical prostatectomy After a Long Learning Curve-Updated Analysis of a Prospective Single-surgeon Cohort of 2206 Consecutive Cases." *Eur Urol* 73(5): 664-671.
81. Trinh, Q. D., J. Sammon, M. Sun, P. Ravi, K. R. Ghani, M. Bianchi, W. Jeong, S. F. Shariat, J. Hansen, J. Schmitges, C. Jeldres, C. G. Rogers, J. O. Peabody, F. Montorsi, M. Menon and P. I. Karakiewicz (2012). "Perioperative outcomes of robot-assisted radical radical prostatectomy compared with open radical radical prostatectomy: results from the nationwide inpatient sample." *Eur Urol* 61(4): 679-685.
82. Tyritzis, S. I., U. Wilderang, W. Lantz Alpha, G. Steineck, J. Hugosson, A. Bjartell, J. Stranne, E. Haglind and N. P. Wiklund (2020). "Hospital readmissions after limited vs. extended lymph node dissection during open and robot-assisted radical radical prostatectomy." *Urol Oncol* 38(1): 5 e1-5 e8.
83. von Bodman, C., M. Brock, F. Roghmann, A. Byers, B. Loppenberg, K. Braun, J. Pastor, F. Sommerer, J. Noldus and R. J. Palisaar (2013). "Intraoperative frozen section of the prostate decreases positive margin rate while ensuring nerve sparing procedure during radical radical prostatectomy." *J Urol* 190(2): 515-520.
84. Wah, W., N. Papa, M. Evans, S. Ahern and A. Earnest (2021). "A multi-level spatio-temporal analysis on prostate cancer outcomes." *Cancer Epidemiology* 72: 101939.
85. Wallerstedt Lantz, A., J. Stranne, S. I. Tyritzis, D. Bock, D. Wallin, H. Nilsson, S. Carlsson, T. Thorsteinsdottir, O. Gustafsson, J. Hugosson, A. Bjartell, P. Wiklund, G. Steineck and E. Haglind (2019). "90-Day readmission after radical radical prostatectomy-a prospective comparison between robot-assisted and open surgery." *Scand J Urol* 53(1): 26-33.
86. Wallerstedt, A., S. I. Tyritzis, T. Thorsteinsdottir, S. Carlsson, J. Stranne, O. Gustafsson, J. Hugosson, A. Bjartell, U. Wilderang, N. P. Wiklund, G. Steineck, E. Haglind and L. s. committe (2015). "Short-term results after robot-assisted laparoscopic radical radical prostatectomy compared to open radical radical prostatectomy." *Eur Urol* 67(4): 660-670.
87. Wang, M., Y. F. Yang, B. D. Guo, H. M. Hou, L. F. Meng, X. Wang, M. Amankwah, B. Q. Liu, C. Y. Jin, M. Liu and J. Y. Wang (2021). "The impact of atrial fibrillation on outcomes in patients undergoing radical radical prostatectomy." *World J Urol* 39(5): 1509-1519.
88. Wang, Y., H. Gieschen, M. Greenberger, X. Yu, G. Tian, N. VanderWalde, T. Stockstill, M. Farmer, L. Rinker, E. W. Izaguirre, B. Somer and M. T. Ballo (2021). "Survival After Robotic-Assisted radical prostatectomy for Localized Prostate Cancer: An Epidemiologic Study." *Annals of surgery* 274(6): e507-e514.
89. 1Weiner, A. B., P. Murthy, K. A. Richards, S. G. Patel and S. E. Eggener (2015). "Population based analysis of incidence and predictors of open conversion during minimally invasive radical radical prostatectomy." *J Urol* 193(3): 826-831.
90. Wen, T., C. M. Deibert, F. S. Siringo and B. A. Spencer (2014). "Positioning-related complications of minimally invasive radical prostatectomies." *J Endourol* 28(6): 660-667.
91. Willis, D. L., M. L. Gonzalgo, M. Brotzman, Z. Feng, B. Trock and L. M. Su (2012). "Comparison of outcomes between pure laparoscopic vs robot-assisted laparoscopic radical radical prostatectomy: a study of comparative effectiveness based upon validated quality of life outcomes." *BJU International* 109(6): 898-905.
92. Wu, S. Y., C. L. Chang, C. I. Chen and C. C. Huang (2021). "Comparison of Acute and Chronic Surgical Complications following Robot-Assisted, Laparoscopic, and Traditional Open Radical radical prostatectomy among Men in Taiwan." *JAMA Network Open* 4(8).
93. Yaxley, J. W., G. D. Coughlin, S. K. Chambers, S. Occhipinti, H. Samaratunga, L. Zajdlewicz, N. Dunglison, R. Carter, S. Williams, D. J. Payton, J. Perry-Keene, M. F. Lavin and R. A. Gardiner (2016). "Robot-assisted laparoscopic radical prostatectomy versus open radical retropubic radical prostatectomy: early outcomes from a randomised controlled phase 3 study." *Lancet* 388(10049): 1057-1066.
94. Yilmazel, F. K., E. Sam, A. E. Cinislioglu, I. H. Tor, F. Akkas, F. Bedir, I. Karabulut, H. R. Aydin, S. Adanur and O. Polat (2022). "Comparison of Perioperative, Oncological, and Functional Outcomes of Three-Dimensional Versus Robot-Assisted Laparoscopic Radical radical prostatectomy: A Preliminary Study." *Journal of laparoendoscopic & advanced surgical techniques. Part A* 32(3): 304-309.
95. Yu, H. Y., N. D. Hevelone, S. R. Lipsitz, K. J. Kowalczyk and J. C. Hu (2012). "Use, costs and comparative effectiveness of robotic assisted, laparoscopic and open urological surgery." *J Urol* 187(4): 1392-1398.

Disclosures

Important Safety Information

(US) Serious complications may occur in any surgery, including da Vinci surgery, up to and including death. Serious risks include, but are not limited to, injury to tissues and organs and conversion to other surgical techniques which could result in a longer operative time and/or increased complications. For summary of the risks associated with surgery refer to www.intuitive.com/safety.

Da Vinci Xi®/da Vinci X® system precaution statement

The demonstration of safety and effectiveness for the representative specific procedures did not include evaluation of outcomes related to the treatment of cancer (overall survival, disease-free survival, local recurrence), except for radical prostatectomy which was evaluated for overall survival, or treatment of the patient's underlying disease/condition. Device usage in all surgical procedures should be guided by the clinical judgment of an adequately trained surgeon.

(EU) Da Vinci X & Xi Surgical Systems

The Intuitive Surgical Endoscopic Instrument Control Systems (da Vinci X and da Vinci Xi Surgical Systems) are intended to assist in the accurate control of Intuitive Surgical Endoscopic Instruments during urologic surgical procedures, general laparoscopic surgical procedures, gynecologic laparoscopic surgical procedures, general thoracoscopic surgical procedures, and trans-oral otolaryngology surgical procedures restricted to benign tumors and malignant tumors classified as T1 and T2, and for benign base of tongue resection procedures. The systems are indicated for adult and pediatric use (except for trans-oral otolaryngology surgical procedures). They are intended to be used by trained physicians in an operating room environment.

The da Vinci X and da Vinci Xi Surgical Systems are class IIb medical devices CE marked (CE 2460) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

For product intended use and/or indications for use, risks, cautions, and warnings and full prescribing information, refer to the associated user manual(s) or visit <https://manuals.intuitivesurgical.com/market>. Some products, features or technologies may not be available in all countries. Product availability is subject to regulatory approval in the specific market. Please contact your local Intuitive representative for product availability in your region.

Individual outcomes may depend on a number of factors—including but not limited to—patient characteristics, disease characteristics, and/or surgeon experience.

Privacy Notice: Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2025 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos, including Intuitive, Da Vinci, and Ion, are trademarks or registered trademarks of Intuitive Surgical or their respective owner.

INTUITIVE

intuitive.com