

Maturing clinical evidence


Recent large scale meta-analysis covering period 2010-2020

Results are based on internal, unpublished meta-analyses of peer reviewed literature for robotic-assisted procedures using the da Vinci surgical system


Maturing clinical evidence

Recent large scale meta-analysis covering period 2010+

Results are based on a meta-analysis of peer reviewed literature for robotic-assisted procedures (right colectomy, LAR/TME, prostatectomy, partial nephrectomy, lobectomy, hysterectomy for endometrial and cervical cancer) published between 2010- 2020. This work was presented at the ISPOR 2021 annual congress. The summary of clinical results are reflective of a pooled analysis of 7 systematic literature reviews, presented by outcome across different surgical procedures. While the meta-analysis results provide a single conclusion that is statistically significant or not statistically significant, these results are subject to variability. The results of this analysis may depend on several factors, including but not limited to patient characteristics, disease characteristics, the procedure of interest, and/or surgeon experience.




Hysterectomy for cervical cancer and endometrial cancer




Partial Nephrectomy




Right Colectomy



LAR/TME



Prostatectomy



Lobectomy

Outcomes measured

Conversions

Blood transfusions

30-day complications

Length of stay

30-day mortality

30-day readmissions

Operative time

Procedure specific 253 publications from 7 procedures

1,380,094 patients who underwent RAS

963,125 patients who underwent laparoscopic surgery or VATS

1,585,084 patients who underwent open surgery

Robotic-assisted procedure with the da Vinci surgical system

Year \geq 2010

Level of Evidence \leq 2a, 2c

RCT, meta-analysis/systematic review, or large database study

Maturing clinical evidence

Clinical value of robotic-assisted surgery using the da Vinci surgical system (unpublished meta-analysis)

Results are based on 253 peer-reviewed publications where the da Vinci surgical system was used. The summary of clinical results reflects a pooled analysis of 7 systematic literature reviews, presented by outcome across different surgical procedures (right colectomy, LAR/TME, prostatectomy, partial nephrectomy, lobectomy, hysterectomy for endometrial and cervical cancer) published between 2010-2020. There was a total of 138,094 patients that underwent RAS using the Da Vinci surgical system 963,125 patients who underwent laparoscopic or VATS, and 1,585,084 patients who underwent open surgery. This work was presented at the ISPOR 2021 annual congress. While the meta-analysis results provide a single conclusion that is statistically significant or not statistically significant, these results are subject to variability. The results of this analysis may depend on several factors, including but not limited to patient characteristics, disease characteristics, the procedure of interest, and/or surgeon experience.

Outcomes that favor RAS

	vs. Lap	vs. Open
Conversions	55% less likely	
Blood transfusions	28% less likely	70% less likely
30-day complications	14% less likely	39% less likely
Length of stay	0.4 days shorter	1.9 days shorter
30-day mortality	33% less likely	57% less likely
30-day readmissions	23% less likely	

Comparable outcomes

30-day readmissions		comparable
---------------------	--	------------

Outcomes that favor lap/open

Operative time	18.8 min longer	37.3 min longer
----------------	------------------------	------------------------

Maturing clinical evidence

Clinical value of robotic-assisted surgery using the da Vinci surgical system (unpublished meta-analysis)

Results are based on 253 peer-reviewed publications where the da Vinci surgical system was used. The summary of clinical results reflects a pooled analysis of 7 systematic literature reviews, presented by outcome across different surgical procedures (right colectomy, LAR/TME, prostatectomy, partial nephrectomy, lobectomy, hysterectomy for endometrial and cervical cancer) published between 2010-2020. There was a total of 138,094 patients that underwent RAS using the Da Vinci surgical system 963,125 patients who underwent laparoscopic or VATS, and 1,585,084 patients who underwent open surgery. This work was presented at the ISPOR 2021 annual congress. While the meta-analysis results provide a single conclusion that is statistically significant or not statistically significant, these results are subject to variability. The results of this analysis may depend on several factors, including but not limited to patient characteristics, disease characteristics, the procedure of interest, and/or surgeon experience.

Outcomes that favor RAS

	vs. Lap	vs. Open
Conversions	3.7% vs. 7.3%	
Blood transfusions	4.6% vs. 5.7%	4.3% vs. 12.3%
30-day complications	17.5% vs. 20.2%	17.2% vs. 23.7%
Length of stay	0.4 days shorter	1.9 days shorter
30-day mortality	0.4% vs. 0.5%	0.3% vs. 0.8%
30-day readmissions	5.9% vs. 8.0%	

Comparable outcomes

30-day readmissions		6.2% vs. 7.8%
---------------------	--	-----------------------------

Outcomes that favor lap/open

Operative time	18.8 min longer	37.3 min longer
----------------	------------------------	------------------------

Appendix

Bibliography of included studies, glossary

Radical prostatectomy: bibliography July 6, 2020

- Anderson, J. E., D. C. Chang, J. K. Parsons and M. A. Talamini (2012). "The first national examination of outcomes and trends in robotic surgery in the United States." *J Am Coll Surg* 215(1): 107-114; discussion 114-106.
- Aning, J. J., G. S. Reilly, S. Fowler, B. Challacombe, J. S. McGrath, P. Sooriakumaran and B. S. o. Oncology (2019). "Perioperative and oncological outcomes of radical prostatectomy for high-risk prostate cancer in the UK: an analysis of surgeon-reported data." *BJU Int* 124(3): 441-448.
- Antonelli, A., C. Palumbo, M. Noale, A. Porreca, S. Maggi, C. Simeone, P. Bassi, F. Bertoni, S. Bracarda, M. Buglione, G. N. Conti, R. Corvo, M. Gacci, V. Mirone, R. Montironi, L. Triggiani, A. Tubaro, W. Artibani and J. T. C. N. R. s. g. Pros (2019). "Impact of Surgical Approach on Patients-Reported Outcomes after Radical Prostatectomy: A Propensity Score-Weighted Analysis from a Multicenter, Prospective, Observational Study (The Pros-IT CNR Study)." *Urol Int* 103(1): 8-18.
- Asimakopoulos, A. D., C. T. Pereira Fraga, F. Annino, P. Pasqualetti, A. A. Calado and C. Mugnier (2011). "Randomized comparison between laparoscopic and robot-assisted nerve-sparing radical prostatectomy." *J Sex Med* 8(5): 1503-1512.
- Basiri, A., J. J. de la Rosette, S. Tabatabaei, H. H. Woo, M. P. Laguna and H. Shemshaki (2018). "Comparison of retropubic, laparoscopic and robotic radical prostatectomy: who is the winner?" *World J Urol* 36(4): 609-621.
- Basto, M., N. Sathianathan, L. Te Marvelde, S. Ryan, J. Goad, N. Lawrentschuk, A. J. Costello, D. A. Moon, A. G. Heriot, J. Butler and D. G. Murphy (2016). "Patterns-of-care and health economic analysis of robot-assisted radical prostatectomy in the Australian public health system." *BJU Int* 117(6): 930-939.
- Bijani, A., A. E. Hebert, M. Davitian, H. May, M. Speers, R. Leung, N. E. Mohamed, H. S. Sacks and A. Tewari (2016). "A Multidimensional Analysis of Prostate Surgery Costs in the United States: Robotic-Assisted versus Retropubic Radical Prostatectomy." *Value Health* 19(4): 391-403.
- Bjorklund, J., Y. Folkvaljon, A. Cole, S. Carlsson, D. Robinson, S. Loeb, P. Stattin and O. Akre (2016). "Postoperative mortality 90 days after robot-assisted laparoscopic prostatectomy and retropubic radical prostatectomy: a nationwide population-based study." *BJU Int* 118(2): 302-306.
- Cao, L., Z. Yang, L. Qi and M. Chen (2019). "Robot-assisted and laparoscopic vs open radical prostatectomy in clinically localized prostate cancer: perioperative, functional, and oncological outcomes: A systematic review and meta-analysis." *Medicine (Baltimore)* 98(22): e15770.
- Carter, S. C., S. Lipsitz, Y. C. Shih, P. L. Nguyen, Q. D. Trinh and J. C. Hu (2014). "Population-Based Determinants of Radical Prostatectomy Operative Time." *BJU International* 113(5b): E112-118.
- Cazaniga, W., R. A. Godtman, S. Carlsson, G. Ahlgren, E. Johansson, D. Robinson, J. Hugosson and P. Stattin (2019). "Population-based, nationwide registration of prostatectomies in Sweden." *Journal of Surgical Oncology* 120(4): 803-812.
- Coughlin, G. D., J. W. Yaxley, S. K. Chambers, S. Occhipinti, H. Samarantunga, L. Zajdelwicz, P. Teloken, N. Dunglison, S. Williams, M. F. Lavin and R. A. Gardiner (2018). "Robot-assisted laparoscopic prostatectomy versus open radical retropubic prostatectomy: 24-month outcomes from a randomised controlled study." *Lancet Oncol* 19(8): 1051-1060.
- Davis, J. W., U. S. Kreaden, J. Gabbert and R. Thomas (2014). "Learning curve assessment of robot-assisted radical prostatectomy compared with open-surgery controls from the premier perspective database." *J Endourol* 28(5): 560-566.
- De Carlo, F., F. Celestino, C. Verri, F. Masedu, E. Liberati and S. M. Di Stasi (2014). "Retropubic, laparoscopic, and robot-assisted radical prostatectomy: surgical, oncological, and functional outcomes: a systematic review." *Urol Int* 93(4): 373-383.
- Ficarra, V., G. Novara, T. E. Ahlering, A. Costello, J. A. Eastham, M. Graefen, G. Guazzoni, M. Menon, A. Mottrie, V. R. Patel, H. Van der Poel, R. C. Rosen, A. K. Tewari, T. G. Wilson, F. Zattoni and F. Montorsi (2012). "Systematic review and meta-analysis of studies reporting potency rates after robot-assisted radical prostatectomy." *Eur Urol* 62(3): 418-430.
- Ficarra, V., G. Novara, R. C. Rosen, W. Artibani, P. R. Carroll, A. Costello, M. Menon, F. Montorsi, V. R. Patel, J. U. Stolzenberg, H. Van der Poel, T. G. Wilson, F. Zattoni and A. Mottrie (2012). "Systematic review and meta-analysis of studies reporting urinary continence recovery after robot-assisted radical prostatectomy." *Eur Urol* 62(3): 405-417.
- Fridriksson, J. O., Y. Folkvaljon, K. J. Lundstrom, D. Robinson, S. Carlsson, P. Stattin and P. Stattin (2017). "Long-term adverse effects after retropubic and robot-assisted radical prostatectomy. Nationwide, population-based study." *J Surg Oncol* 116(4): 500-506.
- Gandaglia, G., J. D. Sammon, S. L. Chang, T. K. Choueiri, J. C. Hu, P. I. Karakiewicz, A. S. Kibel, S. P. Kim, R. Konijeti, F. Montorsi, P. L. Nguyen, S. Sukumar, M. Menon, M. Sun and Q. D. Trinh (2014). "Comparative effectiveness of robot-assisted and open radical prostatectomy in the postdissemination era." *J Clin Oncol* 32(14): 1419-1426.
- Gandaglia, G., Q. D. Trinh, J. C. Hu, J. Schifflmann, A. Becker, F. Roghmann, I. Popa, Z. Tian, P. Perrotte, F. Montorsi, A. Briganti, P. I. Karakiewicz, M. Sun and F. Abdollah (2014). "The impact of robot-assisted radical prostatectomy on the use and extent of pelvic lymph node dissection in the "post-dissemination" period." *Eur J Surg Oncol* 40(9): 1080-1086.
- Groeben, C., R. Koch, M. Baunacke, M. P. Wirth and J. Huber (2017). "High volume is the key for improving in-hospital outcomes after radical prostatectomy: a total population analysis in Germany from 2006 to 2013." *World J Urol* 35(7): 1045-1053.
- Herlemann, A., J. E. Cowan, P. R. Carroll and M. R. Cooperberg (2018). "Community-based Outcomes of Open versus Robot-assisted Radical Prostatectomy." *Eur Urol* 73(2): 215-223.
- Hu, J. C., G. Gandaglia, P. I. Karakiewicz, P. L. Nguyen, Q. D. Trinh, Y.-C. T. Shih, F. Abdollah, K. Chiamie, J. L. Wright and P. A. Ganz (2014). "Comparative Effectiveness of Robot-assisted Versus Open Radical Prostatectomy Cancer Control." *European Urology* 66: 666-672.
- Huang, X., L. Wang, X. Zheng and X. Wang (2017). "Comparison of perioperative, functional, and oncologic outcomes between standard laparoscopic and robotic-assisted radical prostatectomy: a systemic review and meta-analysis." *Surg Endosc* 31(3): 1045-1060.
- Hyams, E. S., J. K. Mullins, P. M. Pierorazio, A. W. Partin, M. E. Allaf and B. R. Matlaga (2013). "Impact of robotic technique and surgical volume on the cost of radical prostatectomy." *J Endourol* 27(3): 298-303.
- Khadhour, S., C. Miller, S. Fowler, L. Hounsme, A. McNeill, J. Adshedge, J. S. McGrath and B. S. o. Oncology (2018). "The British Association of Urological Surgeons (BAUS) radical prostatectomy audit 2014/2015 - an update on current practice and outcomes by centre and surgeon case-volume." *BJU Int* 121(6): 886-892.
- (2015). "Contemporary practice and technique-related outcomes for radical prostatectomy in the UK: a report of national outcomes." *BJU Int* 115(5): 753-763.
- Lee, S. H., H. J. Seo, N. R. Lee, S. K. Son, D. K. Kim and K. H. Rha (2017). "Robot-assisted radical prostatectomy has lower biochemical recurrence than laparoscopic radical prostatectomy: Systematic review and meta-analysis." *Investig Clin Urol* 58(3): 152-163.
- Lundstrom, K. J., Y. Folkvaljon, S. Loeb, A. B. Axelson, P. Stattin and P. Nordin (2016). "Small bowel obstruction and abdominal pain after retrobic versus open radical prostatectomy." *Scand J Urol* 50(3): 155-159.
- Moran, P. S., M. O'Neill, C. Tejelar, M. Flattery, L. A. Murphy, G. Smyth and M. Ryan (2013). "Robot-assisted radical prostatectomy compared with open and laparoscopic approaches: a systematic review and meta-analysis." *Int J Urol* 20(3): 312-321.
- Nabi, J., D. F. Friedlander, X. Chen, A. P. Cole, J. C. Hu, A. S. Kibel, P. Dasgupta and Q. D. Trinh (2020). "Assessment of Out-of-Pocket Costs for Robotic Cancer Surgery in US Adults." *JAMA network open* 3(1): e1919185.
- Ning, C., X. Hu, F. Liu, J. Lin, J. Zhang, Z. Wang and Y. Zhu (2019). "Post-surgical outcomes of patients with chronic kidney disease and end stage renal disease undergoing radical prostatectomy: 10-year results from the US National Inpatient Sample." *BMC Nephrology* 20(1): 278.
- Nossiter, J., A. Sujenthiran, S. C. Charman, P. J. Cathcart, A. Aggarwal, H. Payne, N. W. Clarke and J. van der Meulen (2018). "Robot-assisted radical prostatectomy vs laparoscopic and open retropubic radical prostatectomy: functional outcomes 18 months after diagnosis from a national cohort study in England." *Br J Cancer* 118(4): 489-494.
- Novara, G., V. Ficarra, S. Moccellini, T. E. Ahlering, P. R. Carroll, M. Graefen, G. Guazzoni, M. Menon, V. R. Patel, S. F. Shariat, A. K. Tewari, H. Van Poppel, F. Zattoni, F. Montorsi, A. Mottrie, R. C. Rosen and T. G. Wilson (2012). "Systematic review and meta-analysis of studies reporting oncologic outcome after robot-assisted radical prostatectomy." *Eur Urol* 62(3): 382-404.
- Novara, G., V. Ficarra, R. C. Rosen, W. Artibani, A. Costello, J. A. Eastham, M. Graefen, G. Guazzoni, S. F. Shariat, J. U. Stolzenberg, H. Van Poppel, F. Zattoni, F. Montorsi, A. Mottrie and T. G. Wilson (2012). "Systematic review and meta-analysis of perioperative outcomes and complications after robot-assisted radical prostatectomy." *Eur Urol* 62(3): 431-452.
- O'Neil, B., T. Koyama, J. Alvarez, R. M. Conwill, P. C. Albertsen, M. R. Cooperberg, M. Goodman, S. Greenfield, A. S. Hamilton, K. E. Hoffman, R. M. Hoffman, S. H. Kaplan, J. L. Stanford, A. M. Stroup, L. E. Paddock, X. C. Wu, R. A. Stephenson, M. J. Resnick, D. A. Barocas and D. F. Penson (2016). "The Comparative Harms of Open and Robotic Prostatectomy in Population Based Samples." *J Urol* 195(2): 321-329.
- Ong, W. L., S. M. Evans, T. Spelman, P. A. Kearns, D. G. Murphy and J. L. Millar (2016). "Comparison of oncological and health-related quality of life outcomes between open and robot-assisted radical prostatectomy for localised prostate cancer - findings from the population-based Victorian Prostate Cancer Registry." *BJU Int* 118(4): 563-569.
- Pan, X., W., X. M. Cui, J. F. Teng, D. X. Zhang, Z. J. Wang, F. J. Qu, Y. Gao, X. G. Cui and D. F. Xu (2015). "Robot-Assisted Radical Prostatectomy vs. Open Retropubic Radical Prostatectomy for Prostate Cancer: A Systematic Review and Meta-analysis." *Indian J Surg* 77(Suppl 3): 1326-1333.
- Pearce, S. M., J. J. Pariser, T. Karrison, S. G. Patel and S. E. Eggener (2016). "Comparison of Perioperative and Early Oncologic Outcomes between Open and Robotic Assisted Laparoscopic Prostatectomy in a Contemporary Population Based Cohort." *J Urol* 196(1): 76-81.
- Pielecki, M. A., B. B. McGuire, U. Jain, J. Y. Kim and R. B. Nadler (2014). "National multi-institutional comparison of 30-day postoperative complication and readmission rates between open retropubic radical prostatectomy and robot-assisted laparoscopic prostatectomy using NSQIP." *J Endourol* 28(4): 430-436.
- Porpiglia, F., C. Fiori, R. Bertolo, M. Manfredi, F. Mele, E. Checucci, S. De Luca, R. Passera and R. M. Scarpa (2018). "Five-year Outcomes for a Prospective Randomised Controlled Trial Comparing Laparoscopic and Robot-assisted Radical Prostatectomy." *Eur Urol Focus* 4(1): 80-86.
- Porpiglia, F., I. Morra, M. Lucci Chiarissi, M. Manfredi, F. Mele, S. Grande, F. Ragni, M. Poggio and C. Fiori (2013). "Randomised controlled trial comparing laparoscopic and robot-assisted radical prostatectomy." *Eur Urol* 63(4): 606-614.
- Preisser, F., S. Nazzari, E. Mazzone, S. Knipper, M. Bandini, Z. Tian, A. Haese, F. Saad, K. C. Zorn, F. Montorsi, S. F. Shariat, M. Graefen, D. Tilki and P. I. Karakiewicz (2019). "Regional differences in total hospital charges between open and robotically assisted radical prostatectomy in the United States." *World J Urol* 37(7): 1305-1313.
- Retel, V. P., C. Bouchardy, M. Usel, I. Neyroud-Caspar, F. Schmidlin, G. Wirth, C. Iselin, R. Miralbell and E. Rapiti (2014). "Determinants and effects of positive surgical margins after prostatectomy on prostate cancer mortality: a population-based study." *BMC Urol* 14(1): 86.
- Robertson, C., A. Close, C. Fraser, T. Gurung, X. Jia, P. Sharma, L. Vale, C. Ramsay and R. Pickard (2013). "Relative effectiveness of robot-assisted and standard laparoscopic prostatectomy as alternatives to open radical prostatectomy for treatment of localised prostate cancer: a systematic review and mixed treatment comparison meta-analysis." *BJU Int* 112(6): 798-812.
- Seo, H. J., N. R. Lee, S. K. Son, D. K. Kim, K. H. Rha and S. H. Lee (2016). "Comparison of Robot-Assisted Radical Prostatectomy and Open Radical Prostatectomy Outcomes: A Systematic Review and Meta-Analysis." *Yonsei Med J* 57(5): 1165-1177.
- Simon, R. M., L. E. Howard, D. M. Moreira, M. K. Terris, C. J. Kane, W. J. Aronson, C. L. Amling, M. R. Cooperberg and S. J. Freedland (2017). "Predictors of operative time during radical retropubic prostatectomy and robot-assisted laparoscopic prostatectomy." *Int J Urol* 24(8): 618-623.
- Sridharan, K. and G. Sivaramkrishnan (2018). "Prostatectomies for localized prostate cancer: a mixed comparison network and cumulative meta-analysis." *J Robot Surg* 12(4): 633-639.
- Srougi, V., J. Bessa, Jr., M. Baghdadi, I. Nunes-Silva, J. B. da Costa, S. Garcia-Barreras, E. Barret, F. Rozet, M. Galiano, R. Sanchez-Salas and X. Cathelineau (2017). "Surgical method influences specimen margins and biochemical recurrence during radical prostatectomy for high-risk prostate cancer: a systematic review and meta-analysis." *World J Urol* 35(10): 1481-1488.
- Steffens, D., R. Thanigaalam, S. Leslie, B. Maneck, J. M. Young and M. Solomon (2017). "Robotic Surgery in Uro-oncology: A Systematic Review and Meta-analysis of Randomized Controlled Trials." *Urology* 106: 9-17.
- Stolzenburg, J. U., I. Kyriazi, C. Fahlenbrach, C. Gilfrich, C. Gunster, E. Jeschke, G. Popken, L. Weissbach, C. von Zastrow and H. Leicht (2016). "National trends and differences in morbidity among surgical approaches for radical prostatectomy in Germany." *World J Urol* 34(11): 1515-1520.
- Sugihara, T., H. Yasunaga, H. Horiguchi, H. Matsui, T. Fujimura, H. Nishimatsu, H. Fukuhara, H. Kume, Y. Changhong, M. W. Kattan, K. Fukushima and Y. Homma (2014). "Robot-assisted versus other types of radical prostatectomy: population-based safety and cost comparison in Japan, 2012-2013." *Cancer Sci* 105(11): 1421-1426.
- Tang, K., K. Jiang, H. Chen, Z. Chen, H. Xu and Z. Ye (2017). "Robotic vs. Retropubic radical prostatectomy in prostate cancer: A systematic review and a meta-analysis update." *Oncotarget* 8(19): 32237-32257.
- Tewari, A., P. Sooriakumaran, D. A. Bloch, U. Seshadri-Kreadan, A. E. Hebert and P. Wiklund (2012). "Positive surgical margin and perioperative complication rates of primary surgical treatments for prostate cancer: a systematic review and meta-analysis comparing retropubic, laparoscopic, and robotic prostatectomy." *Eur Urol* 62(1): 1-15.
- Trinh, Q. D., J. Sammon, M. Sun, P. Ravi, K. R. Ghani, M. Bianchi, W. Jeong, S. F. Shariat, J. Hansen, J. Schmitges, C. Jeldres, C. G. Rogers, J. O. Peabody, F. Montorsi, M. Menon and P. I. Karakiewicz (2012). "Perioperative outcomes of robot-assisted radical prostatectomy compared with open radical prostatectomy: results from the nationwide inpatient sample." *Eur Urol* 61(4): 679-685.
- Tyritzis, S. I., U. Wilderang, W. Lantz Alpha, G. Steineck, J. Hugosson, A. Bjartell, J. Stranne, E. Haglund and N. P. Wiklund (2019). "Hospital readmissions after limited vs. extended lymph node dissection during open and robot-assisted radical prostatectomy." *Urol Oncol* 38(1): 5 e1-5 e8.
- Wallerstedt Lantz, A., J. Stranne, S. I. Tyritzis, D. Bock, D. Wallin, H. Nilsson, S. Carlsson, T. Thorstensdotir, O. Gustafsson, J. Hugosson, A. Bjartell, P. Wiklund, G. Steineck and E. Haglund (2019). "90-Day readmission after radical prostatectomy: a prospective comparison between robot-assisted and open surgery." *Scand J Urol* 53(1): 26-33.
- Wang, M., Y. F. Yang, B. D. Guo, H. M. Hou, L. F. Meng, X. Wang, M. Amankwah, B. Q. Liu, C. Y. Jin, M. Liu and J. Y. Wang (2020). "The impact of arial fibrillation on outcomes in patients undergoing radical prostatectomy." *World J Urol*.
- Wang, Y., H. Gieschen, M. Greenberger, X. Yu, G. Tian, N. VanderWalde, T. Stockstill, M. Farmer, L. Rinker, E. W. Izaguirre, B. Somer and M. T. Ballo (2019). "Survival After Robotic-Assisted Prostatectomy for Localized Prostate Cancer: An Epidemiologic Study." *Annals of surgery*.
- Weiner, A. B., P. Murthy, K. A. Richards, S. G. Patel and S. E. Eggener (2015). "Population based analysis of incidence and predictors of open versus radical minimally invasive radical prostatectomy." *J Urol* 193(3): 826-831.
- Wen, T., C. M. Delbert, F. S. Siringo and B. A. Spencer (2014). "Positioning-related complications of minimally invasive radical prostatectomies." *J Endourol* 28(6): 660-667.
- Yaxley, J. W., G. D. Coughlin, S. K. Chambers, S. Occhipinti, H. Samarantunga, L. Zajdelwicz, N. Dunglison, R. Carter, S. Williams, D. J. Payton, J. Perry-Keene, M. F. Lavin and R. A. Gardiner (2016). "Robot-assisted laparoscopic prostatectomy versus open radical retropubic prostatectomy: early outcomes from a randomised controlled phase 3 study." *Lancet* 388(10049): 1057-1066.
- Yu, H. Y., N. D. Hevelone, S. R. Lipsitz, K. J. Kowalczyk and J. C. Hu (2012). "Use, costs and comparative effectiveness of robotic assisted, laparoscopic and open urological surgery." *J Urol* 187(4): 1392-1398.
- Not Included in Forest Plots
- Barashi, N. S., S. M. Pearce, A. J. Cohen, J. J. Pariser, V. T. Packiam and S. E. Eggener (2018). "Incidence, Risk Factors, and Outcomes for Rectal Injury During Radical Prostatectomy: A Population-based Study." *European Urology* 1(6): 501-506. [reports rectal injury]
- Du, Y., Q. Long, B. Guan, L. Mu, J. Tian, Y. Jiang, X. Bai and D. Wu (2018). "Robot-Assisted Radical Prostatectomy Is More Beneficial for Prostate Cancer Patients: A System Review and Meta-Analysis." *Med Sci Monit* 24: 272-287. [errors]
- Friethriksson, J. O., E. Holmberg, J. Adolffson, M. Lambe, A. Bill-Axelsson, S. Carlsson, J. Hugosson and P. Stattin (2014). "Rehospitalization after radical prostatectomy in a nationwide, population based study." *J Urol* 192(1): 112-119. [reports perioperative hospital stays, which includes readmissions]
- Hyltdgard, V. B., K. R. Larsen, J. Poulsen and R. Sogaard (2017). "Robot-assisted surgery in a broader healthcare perspective: a difference-in-difference-based cost analysis of a national prostatectomy cohort." *BMJ Open* 7(7): e015580. [reports hospital stay for 12 months before surgery]
- Kostakis, I. D., H. Sran, R. Uwechue, P. Chandak, J. Olsburgh, N. Mamode, I. Loukopoulos and N. Kessarisi (2019). "Comparison Between Robotic and Laparoscopic or Open Anastomoses: A Systematic Review and Meta-Analysis." *Robot Surg* 6: 27-40. [reports anastomotic leaks and strictures]
- Marra, A. R., M. Puig-Asensio, M. B. Edmond, M. L. Schweizer and K. G. Nepple (2019). "Infectious Complications of Conventional Laparoscopic vs Robotic Laparoscopic Prostatectomy: A Systematic Literature Review and Meta-Analysis." *J Endourol* 33(3): 179-188. [reports infectious complication rate]
- Pucheril, D., S. A. Fletcher, X. Chen, D. F. Friedlander, A. P. Cole, M. J. Krimphove, A. C. Fields, N. Melnitchouk, A. S. Kibel, P. Dasgupta and Q. D. Trinh (2020). "Workplace absenteeism amongst patients undergoing open vs. robotic radical prostatectomy, hysterectomy, and partial colectomy." *Surgical Endoscopy*. [reports cumulative number of missed work-days from 14 days preoperatively to 352 days postoperatively]
- Shkolyar, E., I. F. Shih, Y. Li, J. Wong and J. C. Liao (2020). "Robotic-Assisted Radical Prostatectomy Associated With Decreased Persistent Postoperative Opioid Use." *J Endourol* 34(4): 475-481. [reports new persistent opioid use 90-180d postoperatively]
- Sujenthiran, A., J. Nossiter, M. Parry, S. C. Charman, A. Aggarwal, H. Payne, P. Dasgupta, N. W. Clarke, J. van der Meulen and P. Cathcart (2018). "National cohort study comparing severe medium-term urinary complications after robot-assisted vs laparoscopic vs retropubic open radical prostatectomy." *BJU Int* 121(3): 445-452. [reports severe urinary complication rate]
- Wang, T., Q. Wang and S. Wang (2019). "A Meta-analysis of Robot Assisted Laparoscopic Radical Prostatectomy Versus Laparoscopic Radical Prostatectomy." *Open Med (Wars)* 14(1): 485-490. [errors]

Partial nephrectomy: bibliography July 6, 2020

1. Aboumarzouk, O. M., R. J. Stein, R. Eyraud, G. P. Haber, P. L. Chlosta, B. K. Somani and J. H. Kaouk (2012). "Robotic versus laparoscopic partial nephrectomy: a systematic review and meta-analysis." *Eur Urol* 62(6): 1023-1033.
2. Ayangbesan, A., D. M. Golombos, R. Golan, P. O'Malley, P. Lewicki, X. Wu and D. S. Scherr (2019). "Surgical Approach Does Not Impact Margin Status After Partial Nephrectomy for Large Renal Masses." *J Endourol* 33(1): 50-60.
3. Bahler, C. D., M. F. Monn, C. K. Flack, A. R. Gramm, T. A. Gardner and C. P. Sundaram (2018). "Assessing Cost of Robotic Utilization in Partial Nephrectomy with Increasing Utilization." *J Endourol* 32(8): 710-716.
4. Cacciamani, G. E., L. G. Medina, T. Gill, A. Abreu, R. Sotelo, W. Artibani and I. S. Gill (2018). "Impact of Surgical Factors on Robotic Partial Nephrectomy Outcomes: Comprehensive Systematic Review and Meta-Analysis." *J Urol* 200(2): 258-274.
5. Camp, C., J. O'Hara, D. Hughes and J. Adshear (2018). "Short-term Outcomes and Costs Following Partial Nephrectomy in England: A Population-based Study." *Eur Urol Focus* 4(4): 579-585.
6. Choi, J. E., J. H. You, D. K. Kim, K. H. Rha and S. H. Lee (2015). "Comparison of perioperative outcomes between robotic and laparoscopic partial nephrectomy: a systematic review and meta-analysis." *Eur Urol* 67(5): 891-901.
7. Davaro, F., J. Roberts, A. May, C. McFerrin, S. Siddiqui and Z. Hamilton (2020). "Robotic surgery does not affect upstaging of T1 renal masses." *J Robot Surg* 14(3): 447-454.
8. Froghi, S., K. Ahmed, M. S. Khan, P. Dasgupta and B. Challacombe (2013). "Evaluation of robotic and laparoscopic partial nephrectomy for small renal tumours (T1a)." *BJU International* 112(4): E322-333.
9. Ghani, K. R., S. Sukumar, J. D. Sammon, C. G. Rogers, Q. D. Trinh and M. Menon (2014). "Practice patterns and outcomes of open and minimally invasive partial nephrectomy since the introduction of robotic partial nephrectomy: results from the nationwide inpatient sample." *J Urol* 191(4): 907-912.
10. Hadjipavlou, M., F. Khan, S. Fowler, A. Joyce, F. X. Keeley and S. Sriprasad (2016). "Partial vs radical nephrectomy for T1 renal tumours: an analysis from the British Association of Urological Surgeons Nephrectomy Audit." *BJU Int* 117(1): 62-71.
11. Kates, M., M. W. Ball, H. D. Patel, M. A. Gorin, P. M. Piorozio and M. E. Allaf (2015). "The financial impact of robotic technology for partial and radical nephrectomy." *J Endourol* 29(3): 317-322.
12. Khanna, A., S. Campbell, P. B. Murthy, K. Ericson, Y. A. Nyame and R. Abouassaly (2020). "Unplanned Conversion from Minimally Invasive to Open Kidney Surgery: The Impact of Robotics." *J Endourol* 34(9): 955-963.
13. Leow, J. J., N. H. Heah, S. L. Chang, Y. L. Chong and K. S. Png (2016). "Outcomes of Robotic versus Laparoscopic Partial Nephrectomy: an Updated Meta-Analysis of 4,919 Patients." *J Urol* 196(5): 1371-1377.
14. Nabi, J., D. F. Friedlander, X. Chen, A. P. Cole, J. C. Hu, A. S. Kibel, P. Dasgupta and Q. D. Trinh (2020). "Assessment of Out-of-Pocket Costs for Robotic Cancer Surgery in US Adults." *JAMA network open* 3(1): e1919185.
15. Pak, J. S., J. J. Lee, K. Bilal, M. Finkelstein and M. A. Palese (2017). "Utilization trends and outcomes up to 3 months of open, laparoscopic, and robotic partial nephrectomy." *J Robot Surg* 11(2): 223-229.
16. Shen, Z., L. Xie, W. Xie, H. Hu, T. Chen, C. Xing, X. Liu, H. Xu, Y. Zhang, Z. Wu, D. Tian and C. Wu (2016). "The comparison of perioperative outcomes of robot-assisted and open partial nephrectomy: a systematic review and meta-analysis." *World J Surg Oncol* 14(1): 220.
17. Tsai, S. H., P. T. Tseng, B. A. Shere, Y. C. Lai, P. Y. Lin, C. K. Wu and M. L. Stoller (2019). "Open versus robotic partial nephrectomy: Systematic review and meta-analysis of contemporary studies." *Int J Med Robot* 15(1): e1963.
18. Wu, Z., M. Li, B. Liu, C. Cai, H. Ye, C. Lv, Q. Yang, J. Sheng, S. Song, L. Qu, L. Xiao, Y. Sun and L. Wang (2014). "Robotic versus Open Partial Nephrectomy: A Systematic Review and Meta-Analysis." *PLoS One* 9(4): e94878.
19. Würnschimmel, C., G. B. Di Pierro, M. Moschini, P. Grande, P. Baumeister, M. Roth, L. Mordasini and A. Mattei (2020). "Robot-Assisted Laparoscopic Partial Nephrectomy Vs Conventional Laparoscopic Partial Nephrectomy: Functional and Surgical Outcomes of a Prospective Single Surgeon Randomized Study." *Journal of endourology* 34(8): 847-855.
20. Xia, L., X. Wang, T. Xu and T. J. Guzzo (2017). "Systematic Review and Meta-Analysis of Comparative Studies Reporting Perioperative Outcomes of Robot-Assisted Partial Nephrectomy Versus Open Partial Nephrectomy." *J Endourol* 31(9): 893-909.
21. Zhang, X., Z. Shen, S. Zhong, Z. Zhu, X. Wang and T. Xu (2013). "Comparison of peri-operative outcomes of robot-assisted vs laparoscopic partial nephrectomy: a meta-analysis." *BJU Int* 112(8): 1133-1142.
22. Zhang, X., J. Yan, Y. Ren, C. Shen, X. Ying and S. Pan (2014). "Robot-assisted versus laparoscopic partial nephrectomy for localized renal tumors: a meta-analysis." *Int J Clin Exp Med* 7(12): 4770-4779.
23. Zhou, N., X. Zhu and L. Man (2019). "Comparison of the efficacy between robot-assisted and conventional laparoscopic nephron-sparing surgery for early-stage renal cell carcinoma." *International Journal of Clinical and Experimental Medicine* 12(4): 3591-3601.

Hysterectomy for cervical cancer: bibliography July 6, 2020

1. Ascittuo, K. C., G. Kalopotharakos, M. Lofgren, T. Hogberg and C. Borgfeldt (2015). "Robot-assisted surgery in cervical cancer patients reduces the time to normal activities of daily living." *Acta Obstet Gynecol Scand* 94(3): 260-265.
2. Esselen, K. M., A. Vitonis, J. Einarsson, M. G. Muto and S. Cohen (2015). "Health Care Disparities in Hysterectomy for Gynecologic Cancers: Data: From the 2012 National Inpatient Sample." *Obstet Gynecol* 126(5): 1029-1039.
3. Hao, X., S. Han and Y. Wang (2015). "Comparison of conventional laparoscopy and robotic radical hysterectomy for early-stage cervical cancer: A meta-analysis." *J Cancer Res Ther* 11 Suppl(2015): C258-C264.
4. Iavazzo, C. and I. D. Gkegkes (2017). "Cost-benefit analysis of robotic surgery in gynaecological oncology." *Best Pract Res Clin Obstet Gynaecol* 45: 7-18.
5. Jin, Y. M., S. S. Liu, J. Chen, Y. N. Chen and C. C. Ren (2018). "Robotic radical hysterectomy is superior to laparoscopic radical hysterectomy and open radical hysterectomy in the treatment of cervical cancer." *PLoS One* 13(3): e0193033.
6. Liu, Z., X. Li, S. Tian, T. Zhu, Y. Yao and Y. Tao (2017). "Superiority of robotic surgery for cervical cancer in comparison with traditional approaches: A systematic review and meta-analysis." *Int J Surg* 40: 145-154.
7. Luo, C., M. Liu and X. Li (2018). "Efficacy and safety outcomes of robotic radical hysterectomy in Chinese older women with cervical cancer compared with laparoscopic radical hysterectomy." *BMC Women's Health* 18(1): 61.
8. Melamed, A., D. J. Margul, L. Chen, N. L. Keating, M. G. Del Carmen, J. Yang, B. L. L. Seagle, A. Alexander, E. L. Barber, L. W. Rice, J. D. Wright, M. Kocherginsky, S. Shahabi and J. A. Rauh-Hain (2018). "Survival after minimally invasive radical hysterectomy for early-stage cervical cancer." *New England Journal of Medicine* 379(20): 1905-1914.
9. O'Neill, M., P. S. Moran, C. Teljeur, O. E. O'Sullivan, B. A. O'Reilly, M. Hewitt, M. Flattery and M. Ryan (2013). "Robot-assisted hysterectomy compared to open and laparoscopic approaches: systematic review and meta-analysis." *Archives of Gynecology and Obstetrics* 287(5): 907-918.
10. Park, D. A., J. E. Yun, S. W. Kim and S. H. Lee (2017). "Surgical and clinical safety and effectiveness of robot-assisted laparoscopic hysterectomy compared to conventional laparoscopy and laparotomy for cervical cancer: A systematic review and meta-analysis." *Eur J Surg Oncol* 43(6): 994-1002.
11. Piedimonte, S., N. Czuzoj-Shulman, W. Gotlib and H. A. Abenheim (2019). "Robotic Radical Hysterectomy for Cervical Cancer: A Population-Based Study of Adoption and Immediate Postoperative Outcomes in the United States." *J Minim Invasive Gynecol* 26(3): 551-557.
12. Ramirez, P. T., M. Frumovitz, R. Pareja, A. Lopez, M. Vieira, R. Ribeiro, A. Buda, X. Yan, Y. Shuzhong, N. Chetty, D. Isla, M. Tamura, T. Zhu, K. P. Robledo, V. GebSKI, R. Asher, V. Behan, J. L. Nicklin, R. L. Coleman and A. Obermair (2018). "Minimally Invasive versus Abdominal Radical Hysterectomy for Cervical Cancer." *The New England journal of medicine* 379(20): 1895-1904.
13. Reza, M., S. Maeso, J. A. Blasco and E. Andradas (2010). "Meta-analysis of observational studies on the safety and effectiveness of robotic gynaecological surgery." *Br J Surg* 97(12): 1772-1783.
14. Shazly, S. A., M. H. Murad, S. C. Dowdy, B. S. Gostout and A. O. Famuyide (2015). "Robotic radical hysterectomy in early stage cervical cancer: A systematic review and meta-analysis." *Gynecol Oncol* 138(2): 457-471.
15. Wright, J. D., T. J. Herzog, A. I. Neugut, W. M. Burke, Y. S. Lu, S. N. Lewin and D. L. Hershman (2012). "Comparative effectiveness of minimally invasive and abdominal radical hysterectomy for cervical cancer." *Gynecol Oncol* 127(1): 11-17.
16. Yu, J., Y. Wang, Y. Li, X. Li, C. Li and J. Shen (2014). "The safety and effectiveness of Da Vinci surgical system compared with that of open and laparoscopic surgery: A separate meta-analysis of high-quality studies." *Medicine* 93(4): e14171.
17. Zhang, S. S., T. Ding, Z. H. Cui, Y. Lv and R. A. Jiang (2019). "Efficacy of robotic radical hysterectomy for cervical cancer compared with that of open and laparoscopic surgery: A separate meta-analysis of high-quality studies." *Medicine* 98(4): e14171.
18. Zhou, J., B. H. Xiong, L. Ma, Y. Cheng, W. Huang and L. Zhao (2016). "Robotic vs laparoscopic radical hysterectomy for cervical cancer: a meta-analysis." *Int J Med Robot* 12(1): 145-154.
19. Alfonzo, E., E. Wallin, L. Ekdahl, C. Staf, A. F. Radestad, P. Reynisson, K. Stalberg, H. Falconer, J. Persson and P. Dahm-Kahler (2019). "No survival difference between robotic and open radical hysterectomy for women with early-stage cervical cancer: results from a nationwide population-based cohort study." *Eur J Cancer* 116: 169-177.
20. Chen, B., M. Ji, P. Li, P. Liu, W. Zou, Z. Zhao, B. Qu, Z. Li, X. Bin, J. Lang, H. Wang and C. Chen (2020). "Comparison between robot-assisted radical hysterectomy and abdominal radical hysterectomy for cervical cancer: A multicentre retrospective study." *Gynecol Oncol* 157(2): 429-436.
21. Hwang, J. H., B. W. Kim, S. R. Kim and J. H. Kim (2020). "Robotic radical hysterectomy is not superior to laparoscopic radical hysterectomy in perioperative urologic complications: A meta-analysis of 23 studies." *J Minim Invasive Gynecol* 27(1): 38-47. [was #20]
22. Jensen, P. T., T. H. Schnack, L. P. Frøding, S. F. Bjørn, H. Lajer, A. Markauskas, K. M. Jochumsen, K. Fuglsang, J. Dinesen, C. H. Søgaard, E. Søgaard-Andersen, M. M. Jensen, A. Knudsen, L. H. Øster and C. Høgdaal (2020). "Survival after a nationwide adoption of robotic minimally invasive surgery for early-stage cervical cancer - A population-based study." *European Journal of Cancer* 128: 47-56. [was #21]
23. Marra, A. R., M. Puig-Asensio, M. B. Edmond, M. L. Schweizer and D. Bender (2019). "Infectious complications of laparoscopic and robotic hysterectomy: a systematic literature review and meta-analysis." *International journal of gynecological cancer : official journal of the International Gynecological Cancer Society* 29(3): 518-530. [was #22]
24. Matsuo, K., S. Matsuzaki, R. S. Mandelbaum, E. J. Chang, M. Klar, K. Matsushima, B. H. Grubbs, L. D. Roman and J. D. Wright (2020). "Minimally invasive radical hysterectomy for early-stage cervical cancer: Volume-outcome relationship in the early experience period." *Gynecologic Oncology* 158(2): 390-396. [was #23]
25. Narducci, F., E. Bogart, T. Hebert, T. Gauthier, P. Collinet, J. M. Classe, F. Lecuru, A. Delest, S. Motton, V. Conri, C. Ferrer, F. Marchal, G. Ferron, A. Probst, J. Thery, M. C. Le Deley, D. Lefebvre, D. Francon, E. Leblanc and E. Lambaudie (2020). "Severe perioperative morbidity after robot-assisted versus conventional laparoscopy in gynecologic oncology: Results of the randomized ROBOGYN-1004 trial." *Gynecologic Oncology* 158(2): 382-389.
26. Obermair, A., R. Asher, R. Pareja, M. Frumovitz, A. Lopez, R. Moretti-Marques, G. Rendon, R. Ribeiro, A. Tsunoda, V. Behan, A. Buda, M. Q. Bernadini, H. Zhao, M. Vieira, J. Walker, N. M. Spirtos, S. Yao, N. Chetty, T. Zhu, D. Isla, M. Tamura, J. Nicklin, K. P. Robledo, V. GebSKI, R. L. Coleman, G. Salvo and P. T. Ramirez (2020). "Incidence of adverse events in minimally invasive vs open radical hysterectomy in early cervical cancer: results of a randomized controlled trial." *Am J Obstet Gynecol* 222(3): 249 e241-249 e210.

Hysterectomy for endometrial cancer: bibliography July 6, 2020

1. Beck, T. L., M. A. Schiff, B. A. Goff and R. R. Urban (2018). "Robotic, Laparoscopic, or Open Hysterectomy: Surgical Outcomes by Approach in Endometrial Cancer." *J Minim Invasive Gynecol* 25(6): 986-993.
2. Borgfeldt, C., G. Kalapotharakos, K. C. Ascituo, M. Lofgren and T. Hogberg (2016). "A population-based registry study evaluating surgery in newly diagnosed uterine cancer." *Acta Obstet Gynecol Scand* 95(8): 901-911.
3. Bregar, A. J., A. Melamed, E. Diver, J. T. Clemmer, S. Uppal, J. O. Schorge, L. W. Rice, M. G. Del Carmen and J. A. Rauh-Hain (2017). "Minimally Invasive Staging Surgery in Women with Early-Stage Endometrial Cancer: Analysis of the National Cancer Data Base." *Ann Surg Oncol* 24(6): 1677-1687.
4. Chan, J. K., A. B. Gardner, K. Taylor, C. A. Thompson, K. Blansit, X. Yu and D. S. Kapp (2015). "Robotic versus laparoscopic versus open surgery in morbidly obese endometrial cancer patients - a comparative analysis of total charges and complication rates." *Gynecol Oncol* 139(2): 300-305.
5. Chen, S. H., Z. A. Li, R. Huang and H. Q. Xue (2016). "Robot-assisted versus conventional laparoscopic surgery for endometrial cancer staging: A meta-analysis." *Taiwan J Obstet Gynecol* 55(4): 488-494.
6. Dubeshter, B., C. Angel, E. Toy, S. Thomas and J. C. Glantz (2013). "Current role of robotic hysterectomy." *Journal of Gynecologic Surgery* 29(4): 174-178.
7. Esselen, K. M., A. Vitonis, J. Einarsson, M. G. Muto and S. Cohen (2015). "Health Care Disparities in Hysterectomy for Gynecologic Cancers: Data From the 2012 National Inpatient Sample." *Obstet Gynecol* 126(5): 1029-1039.
8. Guy, M. S., J. Sheeder, K. Behbakht, J. D. Wright and S. R. Guntupalli (2016). "Comparative outcomes in older and younger women undergoing laparotomy or robotic surgical staging for endometrial cancer." *Am J Obstet Gynecol* 214(3): 350 e351-350 e310.
9. Iavazzo, C. and I. D. Gkegkes (2017). "Cost-benefit analysis of robotic surgery in gynaecological oncology." *Best Pract Res Clin Obstet Gynaecol* 45: 7-18.
10. Ind, T., A. Laios, M. Hacking and M. Nobbenhuis (2017). "A comparison of operative outcomes between standard and robotic laparoscopic surgery for endometrial cancer: A systematic review and meta-analysis." *Int J Med Robot* 13(4).
11. Jorgensen, S. L., O. Mogensen, C. Wu, K. Lund, M. Iachina, M. Korsholm and P. T. Jensen (2019). "Nationwide Introduction of Minimally Invasive Robotic Surgery for Early-Stage Endometrial Cancer and Its Association With Severe Complications." *JAMA Surg*. In Press.
12. Jørgensen, S. L., O. Mogensen, C. S. Wu, M. Korsholm, K. Lund and P. T. Jensen (2019). "Survival after a nationwide introduction of robotic surgery in women with early-stage endometrial cancer: a population-based prospective cohort study." *European Journal of Cancer* 109: 1-11.
13. Maenpaa, M. M., K. Nieminen, E. I. Tomas, M. Laurila, T. H. Luukkaala and J. U. Maenpaa (2016). "Robotic-assisted vs traditional laparoscopic surgery for endometrial cancer: a randomized controlled trial." *Am J Obstet Gynecol* 215(5): 588 e581-588 e587.
14. O'Neill, M., P. S. Moran, C. Teljeur, O. E. O'Sullivan, B. A. O'Reilly, M. Hewitt, M. Flattery and M. Ryan (2013). "Robot-assisted hysterectomy compared to open and laparoscopic approaches: systematic review and meta-analysis." *Archives of Gynecology and Obstetrics* 287(5): 907-918.
15. Park, D. A., D. H. Lee, S. W. Kim and S. H. Lee (2016). "Comparative safety and effectiveness of robot-assisted laparoscopic hysterectomy versus conventional laparoscopy and laparotomy for endometrial cancer: A systematic review and meta-analysis." *Eur J Surg Oncol* 42(9): 1303-1314.
16. Ran, L., J. Jin, Y. Xu, Y. Bu and F. Song (2014). "Comparison of robotic surgery with laparoscopy and laparotomy for treatment of endometrial cancer: a meta-analysis." *PLoS One* 9(9): e108361.
17. Reza, M., S. Maeso, J. A. Blasco and E. Andradas (2010). "Meta-analysis of observational studies on the safety and effectiveness of robotic gynaecological surgery." *Br J Surg* 97(12): 1772-1783.
18. Safdieh, J., Y. C. Lee, A. Wong, A. Lee, J. P. Weiner, D. Schwartz and D. Schreiber (2017). "A Comparison of Outcomes Between Open Hysterectomy and Robotic-Assisted Hysterectomy for Endometrial Cancer Using the National Cancer Database." *Int J Gynecol Cancer* 27(7): 1508-1516.
19. Salehi, S., E. Åvall-Lundqvist, Y. Brandberg, H. Johansson, C. Suzuki and H. Falconer (2019). "Lymphedema, serious adverse events, and imaging 1 year after comprehensive staging for endometrial cancer: results from the RASHEC trial." *Int J Gynecol Cancer* 29(1): 86-93.
20. Salehi, S., E. Åvall-Lundqvist, B. Legerstam, J. W. Carlson and H. Falconer (2017). "Robot-assisted laparoscopy versus laparotomy for infrarenal paraaortic lymphadenectomy in women with high-risk endometrial cancer: A randomised controlled trial." *European Journal of Cancer* 79: 81-89.
21. Salehi, S., Y. Brandberg, E. Åvall-Lundqvist, C. Suzuki, H. Johansson, B. Legerstam and H. Falconer (2018). "Long-term quality of life after comprehensive surgical staging of high-risk endometrial cancer—results from the RASHEC trial." *Acta Oncologica* 57(12): 1671-1676.
22. Silva, E. S. A., J. P. M. de Carvalho, C. Anton, R. P. Fernandes, E. C. Baracat and J. P. Carvalho (2018). "Introduction of robotic surgery for endometrial cancer into a Brazilian cancer service: a randomized trial evaluating perioperative clinical outcomes and costs." *Clinics (Sao Paulo)* 73(suppl 1): e522s.
23. Somashekhar, S. P., R. C. Jaka and S. S. Zaveri (2014). "Prospective Randomized Study Comparing Robotic-Assisted Hysterectomy and Regional Lymphadenectomy with Traditional Laparotomy for Staging of Endometrial Carcinoma—Initial Indian Experience." *Indian J Surg Oncol* 5(3): 217-223.
24. Wright, J. D., W. M. Burke, A. I. Tergas, J. Y. Hou, Y. Huang, J. C. Hu, G. C. Hillyer, C. V. Ananth, A. I. Neugut and D. L. Hershman (2016). "Comparative Effectiveness of Minimally Invasive Hysterectomy for Endometrial Cancer." *J Clin Oncol* 34(10): 1087-1096.
25. Wright, J. D., W. M. Burke, E. T. Wilde, S. N. Lewin, A. S. Charles, J. H. Kim, N. Goldman, A. I. Neugut, T. J. Herzog and D. L. Hershman (2012). "Comparative effectiveness of robotic versus laparoscopic hysterectomy for endometrial cancer." *J Clin Oncol* 30(8): 783-791.
26. Wright, J. D., L. J. Havrilesky, D. E. Cohn, Y. Huang, J. Rathbun, L. W. Rice, C. L. Brown, R. D. Alvarez and E. M. Ko (2018). "Estimating potential for savings for low risk endometrial cancer using the Endometrial Cancer Alternative Payment Model (ECAP): A companion paper to the Society of Gynecologic Oncology Report on the Endometrial Cancer Alternative Payment Model." *Gynecologic Oncology* 149(2): 241-247.
27. Xie, W., D. Cao, J. Yang, K. Shen and L. Zhao (2016). "Robot-assisted surgery versus conventional laparoscopic surgery for endometrial cancer: a systematic review and meta-analysis." *J Cancer Res Clin Oncol* 142(10): 2173-2183.
28. Yu, X., D. Lum, T. K. Kiet, K. C. Fuh, J. Orr, Jr., R. A. Brooks, S. M. Ueda, L. M. Chen, D. S. Kapp and J. K. Chan (2013). "Utilization of and charges for robotic versus laparoscopic versus open surgery for endometrial cancer." *J Surg Oncol* 107(6): 653-658.
29. Zakhari, A., N. Czuzoj-Shulman, A. R. Spence, W. H. Gotlieb and H. A. Abenhaim (2015). "Laparoscopic and robot-assisted hysterectomy for uterine cancer: a comparison of costs and complications." *Am J Obstet Gynecol* 213(5):665 e661-667.
30. Zakhari, A., N. Czuzoj-Shulman, A. R. Spence, W. H. Gotlieb and H. A. Abenhaim (2016). "Hysterectomy for Uterine Cancer in the Elderly: A Comparison Between Laparoscopic and Robot-Assisted Techniques." *Int J Gynecol Cancer* 26(7): 1222-1227.
31. Casarin, J., C. Song, F. Multinu, S. Cappuccio, E. Liu, K. A. Butler, G. E. Glaser, W. A. Cliby, C. L. Langstraat, F. Ghezzi, A. Z. Fu and A. Mariani (2020). "Implementing robotic surgery for uterine cancer in the United States: Better outcomes without increased costs." *Gynecologic Oncology* 156(2): 451-458.
32. Lundin, E. S., N. B. Wodlin, L. Nilsson and P. Kjolhede (2019). "A prospective randomized assessment of quality of life between open and robotic hysterectomy in early endometrial cancer." *Int J Gynecol Cancer*.
33. Marra, A. R., M. Puig-Asensio, M. B. Edmond, M. L. Schweizer and D. Bender (2019). "Infectious complications of laparoscopic and robotic hysterectomy: a systematic literature review and meta-analysis." *International journal of gynecological cancer : official journal of the International Gynecological Cancer Society* 29(3): 518-530.
34. Narducci, F., E. Bogart, T. Hebert, T. Gauthier, P. Collinet, J. M. Classe, F. Lecuru, A. Delest, S. Motton, V. Conri, C. Ferrer, F. Marchal, G. Ferron, A. Probst, J. Thery, M. C. Le Deley, D. Lefebvre, D. Francon, E. Leblanc and E. Lambaudie (2020). "Severe perioperative morbidity after robot-assisted versus conventional laparoscopy in gynecologic oncology: Results of the randomized ROBOGYN-1004 trial." *Gynecologic Oncology* 158(2): 382-389.
35. Nasioudis, D., M. K. Frey, E. Chapman-Davis, T. A. Caputo and K. Holcomb (2020). "Outcomes of minimally invasive surgery for patients with endometrial carcinoma involving the cervix." *International Journal of Gynecological Cancer* 30(5): 619-625.
36. Rozenholc, A., V. Samouelian, T. Warkus, P. Gauthier, D. Provencher, P. Sauthier, F. Gauthier, P. Drakopoulos and B. Cormier (2019). "Green versus blue: Randomized controlled trial comparing indocyanine green with methylene blue for sentinel lymph node detection in endometrial cancer." *Gynecol Oncol* 153(3): 500-504.
37. Wang, J., X. Li, H. Wu, Y. Zhang and F. Wang (2020). "A Meta-Analysis of Robotic Surgery in Endometrial Cancer: Comparison with Laparoscopy and Laparotomy." *Disease Markers* 2020: 2503753.

LAR/TME/ISR: bibliography July 6, 2020

- Bedrikovetski, S., N. N. Dudi-Venkata, H. M. Kroon, J. W. Moore, R. A. Hunter and T. Sammour (2020). "Outcomes of Minimally Invasive Versus Open Proctectomy for Rectal Cancer: A Propensity-Matched Analysis of Bi-National Colorectal Cancer Audit Data." *Dis Colon Rectum* 63(6): 778-787.
- Bhangu, A. and g. European Society of Coloproctology collaborating (2018). "An international multicentre prospective audit of elective rectal cancer surgery; operative approach versus outcome, including transanal total mesorectal excision (TaTME)." *Colorectal Dis* 20 Suppl 6: 33-46.
- Bolton, W. S., S. J. Chapman, N. Corrigan, J. Croft, F. Collinson, J. M. Brown and D. G. Jayne (2021). "The Incidence of Low Anterior Resection Syndrome as Assessed in an International Randomized Controlled Trial (MRC/NIHR ROLARR)." *Ann Surg* 274(6): e1223-e1229.
- Broholm, M., H. C. Pommergaard and I. Gogenur (2015). "Possible benefits of robot-assisted rectal cancer surgery regarding urological and sexual dysfunction: a systematic review and meta-analysis." *Colorectal Dis* 17(5): 375-381.
- Chapman, B. C., M. Edgcomb, A. Gleisner and J. D. Vogel (2019). "Outcomes in rectal cancer patients undergoing laparoscopic or robotic low anterior resection compared to open: a propensity-matched analysis of the NCDB (2010-2015)." *Surgical Endoscopy* 34(11): 4754-4771.
- Clark, J., K. Shetty, M. H. Sodergren, D. R. James, S. Purkayastha, T. Athanasiou, G. Z. Yang and A. Darzi (2012). "Robot-assisted total mesorectal excision: should it be considered as the technique of choice in the management of rectal cancer?" *J Robot Surg* 6(2): 99-114.
- Corrigan, N., H. Marshall, J. Croft, J. Copeland, D. Jayne and J. Brown (2018). "Exploring and adjusting for potential learning effects in ROLARR: a randomised controlled trial comparing robotic-assisted vs. standard laparoscopic surgery for rectal cancer resection." *Trials* 19(1): 339.
- Cui, Y., C. Li, Z. Xu, Y. Wang, Y. Sun, H. Xu, Z. Li and Y. Sun (2017). "Robot-assisted versus conventional laparoscopic operation in anus-preserving rectal cancer: a meta-analysis." *Ther Clin Risk Manag* 13: 1247-1257.
- Debakay, Y., A. Zaghloul, A. Farag, A. Mahmoud and I. Elattar (2018). "Robotic-Assisted versus Conventional Laparoscopic Approach for Rectal Cancer Surgery, First Egyptian Academic Center Experience, RCT." *Minimally Invasive Surgery* 2018: 5835652.
- Eltair, M., S. Hajibandeh, S. Hajibandeh, A. Nuno, K. H. Abdullah, A. Alkaili-Alyamani, M. I. Aslam, A. Sinha and T. Agarwal (2020). "Meta-analysis and trial sequential analysis of robotic versus laparoscopic total mesorectal excision in management of rectal cancer." *International Journal of Colorectal Disease* 35(8): 1423-1438.
- Gavrilidis, P., J. Wheeler, A. Spinelli, N. de'Angelis, C. Simopoulos and S. Di Saverio (2020). "Robotic vs laparoscopic total mesorectal excision for rectal cancers: has a paradigm change occurred? A systematic review by updated meta-analysis." *Colorectal Dis* 22(11): 1506-1517.
- Gilmore, B., M. A. Adam, K. Rhodin, M. C. Turner, B. Ezekian, C. R. Mantyh and J. Migaly (2020). "Evolution of minimally invasive surgery for rectal cancer: update from the national cancer database." *Surgical Endoscopy* 35(1): 275-290.
- Han, C., P. Yan, W. Jing, M. Li, B. Du, M. Si, J. Yang, K. Yang, H. Cai and T. Guo (2020). "Clinical, pathological, and oncologic outcomes of robotic-assisted versus laparoscopic proctectomy for rectal cancer: A meta-analysis of randomized controlled studies." *Asian Journal of Surgery* 43(9): 880-890.
- Hayes, J. W., E. J. Ryan, P. A. Boland, B. Creavin, M. E. Kelly and D. Beddy (2019). "The prevalence of venous thromboembolism in rectal surgery: a systematic review and meta-analysis." *Int J Colorectal Dis* 34(5): 849-860.
- Hu, K. Y., R. Wu, A. Szabo, T. J. Ridolfi, K. A. Ludwig and C. Y. Peterson (2020). "Laparoscopic Versus Robotic Proctectomy Outcomes: An ACS-NSQIP Analysis." *J Surg Res* 255: 495-501.
- Huang, Y. J., Y. N. Kang, Y. M. Huang, A. T. Wu, W. Wang and P. L. Wei (2019). "Effects of laparoscopic vs robotic-assisted mesorectal excision for rectal cancer: An update systematic review and meta-analysis of randomized controlled trials." *Asian J Surg* 42(6): 657-666.
- Hyde, L. Z., O. Baser, S. Mehendale, D. Guo, M. Shah and R. P. Kiran (2019). "Impact of surgical approach on short-term oncological outcomes and recovery following low anterior resection for rectal cancer." *Colorectal disease : the official journal of the Association of Coloproctology of Great Britain and Ireland* 21(8): 932-942.
- Jayne, D., A. Pigazzi, H. Marshall, J. Croft, N. Corrigan, J. Copeland, P. Quirke, N. West, T. Rautio, N. Thomassen, H. Tilney, M. Gudgeon, P. P. Bianchi, R. Edlin, C. Hulme and J. Brown (2017). "Effect of Robotic-Assisted vs Conventional Laparoscopic Surgery on Risk of Conversion to Open Laparotomy Among Patients Undergoing Resection for Rectal Cancer: The ROLARR Randomized Clinical Trial." *JAMA - Journal of the American Medical Association* 318(16): 1569-1580.
- Jones, K., M. G. Qassem, P. Sains, M. K. Baig and M. S. Sajid (2018). "Robotic total meso-rectal excision for rectal cancer: A systematic review following the publication of the ROLARR trial." *World J Gastrointest Oncol* 10(11): 449-464.
- Kim, M. J., S. C. Park, J. W. Park, H. J. Chang, D. Y. Kim, B. H. Nam, D. K. Sohn and J. H. Oh (2018). "Robot-assisted Versus Laparoscopic Surgery for Rectal Cancer: A Phase II Open Label Prospective Randomized Controlled Trial." *Ann Surg* 267(2): 243-251.
- Kowalewski, K. F., L. Seifert, S. Ali, M. W. Schmidt, S. Seide, C. Haney, C. Tapking, A. Shamiyeh, Y. Kulu, T. Hackert, B. P. Müller-Stich and F. Nickel (2020). "Functional outcomes after laparoscopic versus robotic-assisted rectal resection: a systematic review and meta-analysis." *Surgical endoscopy* 35(1): 81-95.
- Lee, S. H., D. H. Kim and S. W. Lim (2018). "Robotic versus laparoscopic intersphincteric resection for low rectal cancer: a systematic review and meta-analysis." *Int J Colorectal Dis* 33(12): 1741-1753.
- Lee, S. H., S. Lim, J. H. Kim and K. Y. Lee (2015). "Robotic versus conventional laparoscopic surgery for rectal cancer: systematic review and meta-analysis." *Ann Surg Treat Res* 89(4): 190-201.
- Li, L., W. Zhang, Y. Guo, X. Wang, H. Yu, B. Du, X. Yang and Y. Luo (2019). "Robotic Versus Laparoscopic Rectal Surgery for Rectal Cancer: A Meta-Analysis of 7 Randomized Controlled Trials." *Surg Innov* 26(4): 497-504.
- Li, X., T. Wang, L. Yao, L. Hu, P. Jin, T. Guo and K. Yang (2017). "The safety and effectiveness of robot-assisted versus laparoscopic TME in patients with rectal cancer: A meta-analysis and systematic review." *Medicine (Baltimore)* 96(29): e7585.
- Liao, G., Y. B. Li, Z. Zhao, X. Li, H. Deng and G. Li (2016). "Robotic-assisted surgery versus open surgery in the treatment of rectal cancer: the current evidence." *Sci Rep* 6: 26981.
- Lin, S., H. G. Jiang, Z. H. Chen, S. Y. Zhou, X. S. Liu and J. R. Yu (2011). "Meta-analysis of robotic and laparoscopic surgery for treatment of rectal cancer." *World J Gastroenterol* 17(47): 5214-5220.
- Memon, S., A. G. Heriot, D. G. Murphy, M. Bressel and A. C. Lynch (2012). "Robotic versus laparoscopic proctectomy for rectal cancer: a meta-analysis." *Ann Surg Oncol* 19(7): 2095-2101.
- Milone, M., M. Manigrasso, N. Velotti, S. Torino, A. Vozza, G. Sarnelli, G. Aprea, F. Maione, N. Gennarelli, M. Musella and G. D. De Palma (2019). "Completeness of total mesorectum excision of laparoscopic versus robotic surgery: a review with a meta-analysis." *International Journal of Colorectal Disease* 34(6): 983-991.
- Ohtani, H., K. Maeda, S. Nomura, O. Shinto, Y. Mizuyama, H. Nakagawa, H. Nagahara, M. Shibutani, T. Fukuoka, R. Amano, K. Hirakawa and M. Ohira (2018). "Meta-analysis of Robot-assisted Versus Laparoscopic Surgery for Rectal Cancer." *In Vivo* 32(3): 611-623.
- Ouyang, M., T. Liao, Y. Lu, L. Deng, Z. Luo, J. Wu, Y. Ju and X. Yao (2019). "Laparoscopic versus Open Surgery in Lateral Lymph Node Dissection for Advanced Rectal Cancer: A meta-analysis." *Gastroenterology Research and Practice* 2019: 7689082.
- Parascandola, S. A., S. Hota, A. D. Sparks, S. Boulos, K. Cavallo, G. Kim and V. Obias (2021). "Trends in utilization, conversion rates, and outcomes for minimally invasive approaches to non-metastatic rectal cancer: a national cancer database analysis." *Surgical Endoscopy* 35(6): 3154-3165.
- Phan, K., H. R. Kahlaee, S. H. Kim and J. W. T. Toh (2019). "Laparoscopic vs. robotic rectal cancer surgery and the effect on conversion rates: a meta-analysis of randomized controlled trials and propensity-score-matched studies." *Techniques in Coloproctology* 23(3): 221-230.
- Prete, F. P., A. Pezzolla, F. Prete, M. Testini, R. Marzioli, A. Patriti, R. M. Jimenez-Rodriguez, A. Gurrado and G. F. M. Strippoli (2018). "Robotic Versus Laparoscopic Minimally Invasive Surgery for Rectal Cancer: A Systematic Review and Meta-analysis of Randomized Controlled Trials." *Ann Surg* 267(6): 1034-1046.
- Qiu, H., D. Yu, S. Ye, R. Shan, J. Ai and J. Shi (2020). "Long-term oncological outcomes in robotic versus laparoscopic approach for rectal cancer: a systematic review and meta-analysis." *Int J Surg* 80: 225-230.
- Richards, C. R., S. R. Steele, M. B. Lustik, S. M. Gillern, R. B. Lim, J. T. Brady, A. R. Althans and A. T. Schluskel (2019). "Safe surgery in the elderly: A review of outcomes following robotic proctectomy from the Nationwide Inpatient Sample in a cross-sectional study." *Annals of Medicine and Surgery* 44: 39-45.
- Rubinkiewicz, M., J. Witowski, K. Zbroja, K. Rozmus, J. Krzywoń and K. Truskiewicz (2019). "A systematic review and meta-analysis of laparoscopic versus robotic rectal surgery with primary anastomosis." *Polski Przegląd Chirurgiczny/ Polish Journal of Surgery* 92(1): 5-11.
- Simillis, C., N. Lal, S. N. Thoukididou, C. Kontovounisios, J. J. Smith, R. Hompes, M. Adamina and P. P. Tekkis (2019). "Open Versus Laparoscopic Versus Robotic Versus Transanal Mesorectal Excision for Rectal Cancer: A Systematic Review and Network Meta-analysis." *Ann Surg* 270(1): 59-68.
- Somashekhar, S. P., K. R. Ashwin, J. Rajashekhar and S. Zaveri (2015). "Prospective Randomized Study Comparing Robotic-Assisted Surgery with Traditional Laparotomy for Rectal Cancer-Indian Study." *Indian J Surg* 77(Suppl 3): 788-794.
- Sujatha-Bhaskar, S., M. D. Jafari, J. V. Gahagan, C. S. Inaba, C. Y. Koh, S. D. Mills, J. C. Carmichael, M. J. Stamos and A. Pigazzi (2017). "Defining the Role of Minimally Invasive Proctectomy for Locally Advanced Rectal Adenocarcinoma." *Ann Surg* 266(4): 574-581.
- Sun, X. Y., L. Xu, J. Y. Lu and G. N. Zhang (2019). "Robotic versus conventional laparoscopic surgery for rectal cancer: systematic review and meta-analysis." *Minim Invasive Ther Allied Technol* 28(3): 135-142.
- Sun, Y., H. Xu, Z. Li, J. Han, W. Song, J. Wang and Z. Xu (2016). "Robotic versus laparoscopic low anterior resection for rectal cancer: a meta-analysis." *World J Surg Oncol* 14(1): 61.
- Sun, Z., J. Kim, M. A. Adam, D. P. Nussbaum, P. J. Speicher, C. R. Mantyh and J. Migaly (2016). "Minimally Invasive Versus Open Low Anterior Resection: Equivalent Survival in a National Analysis of 14,033 Patients With Rectal Cancer." *Ann Surg* 263(6): 1152-1158.
- Tang, X., Z. Wang, X. Wu, M. Yang and D. Wang (2018). "Robotic versus laparoscopic surgery for rectal cancer in male urogenital function preservation, a meta-analysis." *World Journal of Surgical Oncology* 16(1): 196.
- Taylor, J. P., M. Stem, A. A. Althumairi, S. L. Gearhart, B. Safar, S. H. Fang and J. E. Efron (2020). "Minimally Invasive Proctectomy for Rectal Cancer: A National Perspective on Short-term Outcomes and Morbidity." *World Journal of Surgery* 44(9): 3130-3140.
- Trastulli, S., E. Farinella, R. Cirocchi, D. Cavaliere, N. Avenia, F. Sciannone, N. Gullà, G. Noya and C. Boselli (2012). "Robotic resection compared with laparoscopic rectal resection for cancer: Systematic review and meta-analysis of short-term outcome." *Colorectal Dis* 14(4): e134-e156.
- Wang, G., Z. Wang, Z. Jiang, J. Liu, J. Zhao and J. Li (2017). "Male urinary and sexual function after robotic pelvic autonomic nerve-preserving surgery for rectal cancer." *Int J Med Robot* 13(1).
- Wang, L., Z. Zhang, L. Gong, Y. Zhan, M. Li, S. Li and Y. Xiao (2019). "A Systematic Review and Bayesian Network Meta-Analysis: Short-Term and Long-Term Outcomes of Three Surgery Procedures Following Neoadjuvant Chemoradiotherapy for Rectal Cancer." *J Laparoendosc Adv Surg Tech A* 29(5): 663-670.
- Wee, I. J. Y., L. J. Kuo and J. C. Ngu (2019). "The impact of robotic colorectal surgery in obese patients: a systematic review, meta-analysis, and meta-regression." *Surg Endosc* 33(11): 3558-3566.
- Xiong, B., L. Ma, W. Huang, Q. Zhao, Y. Cheng and J. Liu (2015). "Robotic versus laparoscopic total mesorectal excision for rectal cancer: a meta-analysis of eight studies." *J Gastrointest Surg* 19(3): 516-526.
- Yang, Y., F. Wang, P. Zhang, C. Shi, Y. Zou, H. Qin and Y. Ma (2012). "Robot-assisted versus conventional laparoscopic surgery for colorectal disease, focusing on rectal cancer: a meta-analysis." *Ann Surg Oncol* 19(12): 3727-3736.
- Zheng, B., X. Zhang, X. Wang, L. Ge, M. Wei, L. Bi, X. Deng, Q. Wang, J. Li and Z. Wang (2020). "A comparison of open, laparoscopic and robotic total mesorectal excision: trial sequential analysis and network meta-analysis." *Colorectal Dis* 22(4): 382-391.

Right colectomy: bibliography July 6, 2020

1. Chang, Y. S., J. X. Wang and D. W. Chang (2015). "A meta-analysis of robotic versus laparoscopic colectomy." *J Surg Res* 195(2): 465-474.
2. Davis, B. R., A. C. Yoo, M. Moore and C. Gunnarsson (2014). "Robotic-assisted versus laparoscopic colectomy: cost and clinical outcomes." *Jsls* 18(2): 211-224.
3. Dolejs, S. C., J. A. Waters, E. P. Ceppa and B. L. Zarzaur (2017). "Laparoscopic versus robotic colectomy: a national surgical quality improvement project analysis." *Surg Endosc* 31(6): 2387-2396.
4. Halabi, W. J., C. Y. Kang, M. D. Jafari, V. Q. Nguyen, J. C. Carmichael, S. Mills, M. J. Stamos and A. Pigazzi (2013). "Robotic-assisted Colorectal Surgery in the United States: A Nationwide Analysis of Trends and Outcomes." *World J Surg* 37(12): 2782-2790.
5. Lorenzon, L., F. Bini, G. Balducci, M. Ferri, P. F. Salvi and F. Marinozzi (2016). "Laparoscopic versus robotic-assisted colectomy and rectal resection: a systematic review and meta-analysis." *Int J Colorectal Dis* 31(2): 161-173.
6. Miller, P. E., H. Dao, N. Paluvoi, M. Bailey, D. Margolin, N. Shah and H. D. Vargas (2016). "Comparison of 30-Day Postoperative Outcomes after Laparoscopic vs Robotic Colectomy." *J Am Coll Surg* 223(2): 369-373.
7. Park, J. S., G. S. Choi, S. Y. Park, H. J. Kim and J. P. Ryuk (2012). "Randomized clinical trial of robot-assisted versus standard laparoscopic right colectomy." *British Journal of Surgery* 99(9): 1219-1226.
8. Petruciani, N., D. Sirimarco, G. Nigri, P. Magistri, M. La Torre, P. Aurello, F. Dangelo and G. Ramacciato (2014). "Robotic right colectomy: A worthwhile procedure? Results of a meta-analysis of trials comparing robotic versus laparoscopic right colectomy." *Journal of Minimal Access Surgery* 11(1): 22-28.
9. Rondelli, F., R. Balzarotti, F. Villa, A. Guerra, N. Avenia, E. Mariani and W. Bugiantella (2015). "Is robot-assisted laparoscopic right colectomy more effective than the conventional laparoscopic procedure? A meta-analysis of short-term outcomes." *Int J Surg* 18: 75-82.
10. Solaini, L., F. Bazzocchi, D. Cavaliere, A. Avanzolini, A. Cucchetti and G. Ercolani (2018). "Robotic versus laparoscopic right colectomy: an updated systematic review and meta-analysis." *Surgical Endoscopy and Other Interventional Techniques* 32(3): 1104-1110.
11. Trastulli, S., R. Cirocchi, J. Desiderio, A. Coratti, S. Guarino, C. Renzi, A. Corsi, C. Boselli, A. Santoro, L. Minelli and A. Parisi (2015). "Robotic versus Laparoscopic Approach in Colonic Resections for Cancer and Benign Diseases: Systematic Review and Meta-Analysis." *PLoS One* 10(7): e0134062.
12. Xu, H., J. Li, Y. Sun, Z. Li, Y. Zhen, B. Wang and Z. Xu (2014). "Robotic versus laparoscopic right colectomy: a meta-analysis." *World J Surg Oncol* 12(1): 274.
13. Zarak, A., A. Castillo, K. Kichler, L. de la Cruz, L. Tamariz and S. Kaza (2015). "Robotic versus laparoscopic surgery for colonic disease: a meta-analysis of postoperative variables." *Surg Endosc* 29(6): 1341-1347.
14. Haskins, I. N., T. Ju, M. Skancke, X. Kuang, R. L. Amdur, F. Brody, V. Obias and S. Agarwal (2018). "Right Colon Resection for Colon Cancer: Does Surgical Approach Matter?" *J Laparoendosc Adv Surg Tech A* 28(10): 1202-1206.
15. Khorgami, Z., W. T. Li, T. N. Jackson, C. A. Howard and G. M. Sclabas (2019). "The cost of robotics: an analysis of the added costs of robotic-assisted versus laparoscopic surgery using the National Inpatient Sample." *Surg Endosc* 33(7): 2217-2221.
16. Ma, S., Y. Chen, Y. Chen, T. Guo, X. Yang, Y. Lu, J. Tian and H. Cai (2019). "Short-term outcomes of robotic-assisted right colectomy compared with laparoscopic surgery: A systematic review and meta-analysis." *Asian J Surg* 42(5): 589-598.
17. Park, J. S., H. Kang, S. Y. Park, H. J. Kim, I. T. Woo, I. K. Park and G. S. Choi (2019). "Long-term oncologic after robotic versus laparoscopic right colectomy: a prospective randomized study." *Surg Endosc* 33(9): 2975-2981.
18. Rausa, E., M. E. Kelly, E. Asti, A. Aiolfi, G. Bonitta and L. Bonavina (2019). "Right hemicolectomy: a network meta-analysis comparing open, laparoscopic-assisted, total laparoscopic, and robotic approach." *Surg Endosc* 33(4): 1020-1032.
19. Sweigert, P. J., E. Eguia, A. N. Kothari, K. A. Ban, M. H. Nelson, M. S. Baker and M. A. Singer (2019). "Do prolonged operative times obviate the benefits associated with minimally invasive colectomy?" *Surgery (United States)* 166(3): 336-341.

Lobectomy: bibliography July 6, 2020

1. Agzarian, J., C. Fahim, Y. Shargall, K. Yasufuku, T. K. Waddell and W. C. Hanna (2016). "The Use of Robotic-Assisted Thoracic Surgery for Lung Resection: A Comprehensive Systematic Review." *Semin Thorac Cardiovasc Surg* 28(1): 182-192.
2. Bailey, K. L., N. Merchant, Y. J. Seo, D. Elashoff, P. Benharash and J. Yanagawa (2019). "Short-Term Readmissions After Open, Thoracoscopic, and Robotic Lobectomy for Lung Cancer Based on the Nationwide Readmissions Database." *World J Surg* 43(5): 1377-1384.
3. Emmert, A., C. Straube, J. Buentzel and C. Roever (2017). "Robotic versus thoracoscopic lung resection: A systematic review and meta-analysis." *Medicine (Baltimore)* 96(35): e7633.
4. Glenn, Z. F., M. Zubair, L. Hussain and K. Grannan (2019). "Comparison of pulmonary lobectomies using robotic and video-assisted thoracoscopic approaches: results from 2010-2013 National Inpatient Sample." *J Cardiovasc Surg (Torino)* 60(4): 526-531.
5. Guo, F., D. Ma, S. Li and M. Adamek (2019). "Compare the prognosis of da Vinci robot-Assisted thoracic surgery (RATS) with video-Assisted thoracic surgery (VATS) for non-small cell lung cancer: A Meta-Analysis." *Medicine (United States)* 98(39): e17089.
6. Hendriksen, B. S., C. S. Hollenbeak, M. D. Taylor and M. F. Reed (2019). "Minimally Invasive Lobectomy Modality and Other Predictors of Conversion to Thoracotomy." *Innovations (Phila)* 14(4): 342-352.
7. Hendriksen, B. S., M. F. Reed, M. D. Taylor and C. S. Hollenbeak (2019). "Readmissions After Lobectomy in an Era of Increasing Minimally Invasive Surgery: A Statewide Analysis." *Innovations (Phila)* 14(5): 453-462.
8. Hennon, M. W., L. H. DeGraaff, A. Groman, T. L. Demmy and S. Yendamuri (2020). "The association of nodal upstaging with surgical approach and its impact on long-term survival after resection of non-small-cell lung cancer." *Eur J Cardiothorac Surg* 57(5): 888-895.
9. Hu, J., Y. Chen, J. Dai, X. Zhu, D. Gonzalez-Rivas, G. Jiang, H. Li and P. Zhang (2020). "Perioperative outcomes of robot-assisted vs video-assisted and traditional open thoracic surgery for lung cancer: A systematic review and network meta-analysis." *Int J Med Robot* 16(5): 1-14.
10. Hu, X. and M. Wang (2019). "Efficacy and safety of robot-assisted thoracic surgery (RATS) compare with video-assisted thoracoscopic surgery (VATS) for lung lobectomy in patients with non-small cell lung cancer." *Comb Chem High Throughput Screen* 22(3): 169-178.
11. Huang, J., C. Li, H. Li, F. Lv, L. Jiang, H. Lin, P. Lu, Q. Luo and W. Xu (2019). "Robot-assisted thoracoscopic surgery versus thoracic surgery for c-N2 stage NSCLC: Short-term outcomes of a randomized trial." *Translational Lung Cancer Research* 8(6): 951-958.
12. Kent, M., T. Wang, R. Whyte, T. Curran, R. Flores and S. Gangadharan (2014). "Open, video-assisted thoracic surgery, and robotic lobectomy: review of a national database." *Ann Thorac Surg* 97(1): 236-244.
13. Kim, M. P., D. T. Nguyen, L. M. Meisenbach, E. A. Graviss and E. Y. Chan (2019). "Da Vinci Xi robot decreases the number of thoracotomy cases in pulmonary resection." *J Thorac Dis* 11(1): 145-153.
14. Liang, H., W. Liang, L. Zhao, D. Chen, J. Zhang, Y. Zhang, S. Tang and J. He (2018). "Robotic Versus Video-assisted Lobectomy/Segmentectomy for Lung Cancer: A Meta-analysis." *Ann Surg* 268(2): 254-259.
15. Louie, B. E., J. L. Wilson, S. Kim, R. J. Cerfolio, B. J. Park, A. S. Farivar, E. Vallieres, R. W. Aye, W. R. Burfeind, Jr. and M. I. Block (2016). "Comparison of Video-Assisted Thoracoscopic Surgery and Robotic Approaches for Clinical Stage I and Stage II Non-Small Cell Lung Cancer Using The Society of Thoracic Surgeons Database." *Ann Thorac Surg* 102(3): 917-924.
16. Martin, J. T., E. B. Durbin, L. Chen, T. Gal, A. Mahan, V. Ferraris and J. Zwischenberger (2016). "Nodal Upstaging During Lung Cancer Resection Is Associated With Surgical Approach." *Annals of Thoracic Surgery* 101(1): 238-244.
17. Ng, C. S. H., J. K. Macdonald, S. Gilbert, A. Z. Khan, Y. T. Kim, B. E. Louie, M. Blair Marshall, R. S. Santos, M. Scarci, Y. Shargal and H. C. Fernando (2019). "Optimal Approach to Lobectomy for Non-Small Cell Lung Cancer: Systemic Review and Meta-Analysis." *Innovations: Technology and Techniques in Cardiothoracic and Vascular Surgery* 14(2): 90-116.
18. Nguyen, D. M., I. S. Sarkaria, C. Song, R. M. Reddy, N. Villamizar, L. J. Herrera, L. Shi, E. Liu, D. Rice and D. S. Oh (2020). "Clinical and economic comparative effectiveness of robotic-assisted, video-assisted thoracoscopic, and open lobectomy." *Journal of Thoracic Disease* 12(3): 296-306.
19. Oh, D. S., R. M. Reddy, M. L. Gorrepati, S. Mehendale and M. F. Reed (2017). "Robotic-Assisted, Video-Assisted Thoracoscopic and Open Lobectomy: Propensity-Matched Analysis of Recent Premier Data." *Ann Thorac Surg* 104(5): 1733-1740.
20. Paul, S., J. Jalbert, A. J. Isaacs, N. K. Altorki, O. W. Isom and A. Sedrakyan (2014). "Comparative effectiveness of robotic-assisted vs thoracoscopic lobectomy." *Chest* 146(6): 1505-1512.
21. Rajaram, R., S. Mohanty, D. J. Bentrem, E. S. Pavey, D. D. Odell, A. Bharat, K. Y. Bilimoria and M. M. DeCamp (2017). "Nationwide Assessment of Robotic Lobectomy for Non-Small Cell Lung Cancer." *Ann Thorac Surg* 103(4): 1092-1100.
22. Ray, M. A., N. R. Faris, M. P. Smeltzer, C. Fehnel, C. Houston-Harris, P. Levy, L. Wiggins, V. Sachdev, T. Robbins, D. Spencer and R. U. Osarogiagbon (2018). "Effectiveness of Implemented Interventions on Pathologic Nodal Staging of Non-Small Cell Lung Cancer." *Annals of Thoracic Surgery* 106(1): 228-234.
23. Sesti, J., R. C. Langan, J. Bell, A. Nguyen, A. L. Turner, P. Hilden, K. Leshchuk, M. Dabrowski and S. Paul (2020). "A Comparative Analysis of Long-term Survival of Robotic vs. Thoracoscopic Lobectomy." *Ann Thorac Surg* 110(4): 1139-1146.
24. Subramanian, M. P., J. Liu, W. C. Chapman, M. A. Olsen, Y. Yan, Y. Liu, T. R. Semenkovich, B. F. Meyers, V. Puri and B. D. Kozower (2019). "Utilization Trends, Outcomes, and Cost in Minimally Invasive Lobectomy." *The Annals of thoracic surgery* 108(6): 1648-1655.
25. Swanson, S. J., D. L. Miller, R. J. McKenna, Jr., J. Howington, M. B. Marshall, A. C. Yoo, M. Moore, C. L. Gunnarsson and B. F. Meyers (2014). "Comparing robot-assisted thoracic surgical lobectomy with conventional video-assisted thoracic surgical lobectomy and wedge resection: results from a multihospital database (Premier)." *J Thorac Cardiovasc Surg* 147(3): 929-937.
26. Tang, A., S. Raja, A. C. Bribriesco, D. P. Raymond, M. Sudarshan, S. C. Murthy and U. Ahmad (2020). "Robotic Approach Offers Similar Nodal Upstaging to Open Lobectomy for Clinical Stage I Non-small Cell Lung Cancer." *Ann Thorac Surg* 110(2): 424-433.
27. Veluswamy, R. R., S. A. Whittaker Brown, G. Mhango, K. Sigel, D. G. Nicastrì, C. B. Smith, M. Bonomi, M. D. Galsky, E. Taioli, A. I. Neugut and J. P. Wisnivesky (2020). "Comparative Effectiveness of Robotic-Assisted Surgery for Resectable Lung Cancer in Older Patients." *Chest* 157(5): 1313-1321.
28. Wei, S., M. Chen, N. Chen and L. Liu (2017). "Feasibility and safety of robot-assisted thoracic surgery for lung lobectomy in patients with non-small cell lung cancer: A systematic review and meta-analysis." *World Journal of Surgical Oncology* 15(1): 98.
29. Yu, Z., Q. Xie, L. Guo, X. Chen, C. Ni, W. Luo, W. Li and L. Ma (2017). "Perioperative outcomes of robotic surgery for the treatment of lung cancer compared to a conventional video-assisted thoracoscopic surgery (VATS) technique." *Oncotarget* 8(53): 91076-91084.
30. Zhang, L. and S. Gao (2015). "Robot-assisted thoracic surgery versus open thoracic surgery for lung cancer: A system review and meta-analysis." *International Journal of Clinical and Experimental Medicine* 8(10): 17804-17810.

Glossary

95% CI	95% confidence interval	OR	odds ratio
EBL	estimated blood loss	RAS	robotic-assisted surgery
HTA	health technology assessment	RCT	randomized controlled trial
I²	test statistic for heterogeneity	RD	risk difference
LNY	lymph node yield	SMD	standardized mean difference
LOE	level of evidence	VATS	video-assisted thoracic surgery
LOS	length of hospital stay	WMD	weighted mean difference

Important safety information

Surgical Risks:

Bowel resection and other colorectal procedures (colectomy, sigmoidectomy, low anterior resection, abdominopelvic resection (APR), intersphincteric resection, proctectomy, rectopexy) include: anastomotic leak, anastomotic stricture, colorectal or anorectal dysfunction.

Pulmonary resection (wedge resection, segmentectomy, lobectomy) include: persistent air leak, pneumonia, prolonged mechanical ventilation >48 hours, atrial fibrillation, acute respiratory distress syndrome (ARDS), chylothorax, re-intubation, arrhythmias, bronchopleural fistula, phrenic nerve injury, esophageal injury, difficulty breathing, collapsed lung, pulmonary volvulus, recurrent laryngeal nerve injury leading to vocal cord dysfunction.

Radical prostatectomy include: surrounding nerve damage that can lead to urinary incontinence and/or erectile dysfunction, rectal or bowel injury, urethral stricture, lymphocele, lymphedema; bowel obstruction.

Nephrectomy: renal insufficiency, urine leak, splenic, hepatic or pancreatic laceration, bowel injury, pneumothorax, diaphragmatic injury, urinary fistula, urinoma, renal infarction, lymphocele

Hysterectomy, malignant: urinary tract injury, vaginal cuff problem (separation, adhesions, granulation tissue, infection, cellulitis, hematoma), bladder injury, bowel injury, vaginal tear or laceration, vaginal shortening, voiding dysfunction, urinary incontinence, fistula formation: vesicovaginal, rectovaginal.

Important Safety Information

Serious complications may occur in any surgery, including surgery with a da Vinci system, up to and including death. Examples of serious or life-threatening complications,

which may require prolonged and/or unexpected hospitalization and/or reoperation, include but are not limited to, one or more of the following: injury to tissues/organs, bleeding, infection, and internal scarring that can cause long-lasting dysfunction/pain.

Risks specific to minimally invasive surgery, including surgery with a da Vinci system, include but are

not limited to, one or more of the following: temporary pain/nerve injury associated with positioning; a longer operative time, the need to convert to an open approach, or the need for additional or larger incision sites. Converting the procedure could result in a longer operative time, a longer time under anesthesia, and could lead to increased complications.

Contraindications applicable to the use of conventional endoscopic instruments also apply to the use of all da Vinci instruments.

For important safety information, including surgical risks and considerations, please also refer to www.intuitive.com/safety. For a product's intended use and/or indications for use, risks, full cautions and warnings, please refer to the associated User Manual(s).

Individual outcomes may depend on a number of factors, including but not limited to patient characteristics, disease characteristics, and/or surgeon experience.

Da Vinci Xi®/da Vinci X® system precaution statement

The demonstration of safety and effectiveness for the representative specific procedures did not include evaluation of outcomes related to the treatment of cancer (overall survival, disease-free survival, local recurrence) or treatment of the patient's underlying disease/condition. Device usage in all surgical procedures should be guided by the clinical judgment of an adequately trained surgeon.

Privacy Notice

Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2024 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos are trademarks or registered trademarks of Intuitive Surgical or their respective owner. See www.intuitive.com/trademarks.

INTUITIVE

intuitive.com