

Robotic-Assisted, Video-Assisted Thoracoscopic and Open Lobectomy: Propensity-Matched Analysis of Recent Premier Data

Oh DS, Reddy RM, Gorrepati ML, Mehendale S, Reed MF. Robotic-Assisted, Video-Assisted Thoracoscopic and Open Lobectomy: Propensity-Matched Analysis of Recent Premier Data. *The Annals of Thoracic Surgery*. 2017;104(5):1733-1740. doi:10.1016/j.athoracsur.2017.06.020.

This study's objective was to provide a comparative analysis of perioperative clinical outcomes from elective robotic-assisted lobectomy (RL), VATS lobectomy (VL), and open lobectomy (OL).

The Premier Healthcare Database was analyzed for lobectomies performed from January 1, 2011 to September 30, 2015. International Classification of Diseases, Ninth Revision, Clinical Modification (ICD-9-CM) diagnosis and procedure codes were used to identify surgical approaches (RL, VL, and OL), complications, and conversions to OL.

Propensity score matching (1:1) for patient and hospital characteristics allowed comparison of RL versus OL (n = 2,775 each) and RL versus VL (n = 2,951 each). The following covariates were used for matching: patient characteristics—age, sex, race, Elixhauser comorbidity score, and type of malignancy; and hospital characteristics—payor type, census region, hospital size

(number of beds), type of hospital facility (academic or community), and location of the facility (urban or rural).

Patient populations were statistically similar after matching. Most patients were age 45 to 80 years (mean 67 years) and had lobectomy for primary neoplasm of the lung. Sex distribution and Elixhauser comorbidity scores also were similar between cohorts. Medicare was the predominant payor, and the most prevalent hospital settings were large, urban teaching institutions.

Analyzed data included perioperative outcomes (operating room time, conversion to open, blood transfusion, length of hospital stay, discharge status, complications, and mortality). Complications were assessed as intraoperative, postoperative (surgery through discharge), index hospitalization (admission to discharge), and 30 days (admission to 30 days).


Key takeaways

- Propensity matched analysis shows that RL is associated with a generally lower postoperative complication rate and improved recovery endpoints compared with both OL and VL.
- There continues to be an increase in the number of RL procedures performed. Together, RL and VL now account for more than half of all lobectomies in this US database, indicating a strong trend to more minimally invasive lobectomies.
- The Premier database captures a more heterogeneous, unselected, real-world population with approximately one-third of the patients operated on in community hospitals after propensity matching. In contrast, The Society of Thoracic Surgeons database predominantly favors high-volume academic or teaching institutions.

Data

Adoption of robotic-assisted lobectomy (green line), video-assisted thoracic lobectomy (blue line), and open lobectomy (red line) as percentages of all elective pulmonary lobectomies from January 1, 2011 to September 30, 2015.

(Q = quarter; VATS = video-assisted thoracoscopic surgery.)



Study limitations

- Oncologic data, such as size of the tumor, stage, recurrence, and survival, could not be extracted from the Premier Healthcare Database. There is a possibility of selection bias; patients who underwent OL may have had larger or more central tumors or bulky mediastinal and hilar lymphadenopathy that led to the selection of OL. However, tumor size should presumably not affect the choice between the two minimally invasive approaches.
- As anticipated in any large administrative database, there is potential for coding errors in the data.

PUBLICATION SUMMARY

Data

Summary of Clinical Outcomes and Complications, Matched Analysis (1:1)

Outcome	RL (n = 2,775)	OL (n = 2,775)	p Value (RL vs. OL)	RL (n = 2,951)	VL (n = 2,951)	p Value (RL vs. VL)
OR time, minutes	275.5 ± 94.6	235.3 ± 93.5	<0.0001	275.1 ± 93.9	247.6 ± 86.8	<0.0001
Complications						
Postoperative ^a	961 (34.6)	1,198 (43.2)	<0.0001	1,007 (34.1)	1,109 (37.6)	0.0061
Thirty-day	1,048 (37.8)	1,272 (45.8)	<0.0001	1,101 (37.3)	1,195 (40.5)	0.0130
Index hospital LOS, days	7.0 ± 5.7	8.9 ± 5.9	<0.0001	6.9 ± 5.5	7.3 ± 7.5	0.0060
Median, days	5	7		5	6	
Discharge status						
Health facility	188 (6.8)	284 (10.2)	<0.0001	190 (6.4)	234 (7.9)	0.0302
Home	2,559 (92.2)	2,442 (88.0)	<0.0001	2,734 (92.7)	2,679 (90.9)	0.0108
Conversion to open	179 (6.5)	NA		187 (6.3)	387 (13.1)	<0.0001
Mortality^b						
Index hospital	28 (1.0)	48 (1.7)	0.0282	27 (0.9)	34 (1.2)	0.4400
Thirty-day	36 (1.3)	62 (2.2)	0.0108	35 (1.2)	40 (1.4)	0.6420
Postoperative bleeding ^c	128 (4.6)	290 (10.5)	<0.0001	130 (4.4)	275 (9.3)	<0.0001
Postoperative transfusion ^c	104 (3.8)	151 (5.4)	0.0032	102 (3.5)	109 (3.7)	0.6740

^a Postoperative complications included all complications occurring after surgery through discharge.

^b Mortality rate was measured from admission to discharge ("Index hospital") and admission to 30 days ("30-day").

^c Postoperative bleeding and transfusions that occurred through discharge were included.

Values are mean ± SD or n (%).

LOS: length of stay; NA: not applicable; OL: open lobectomy; OR: operating room

RL: robotic-assisted lobectomy; VATS: video-assisted thoracic surgery; VL: VATS lobectomy

Results

- During the study period, the number of lobectomies performed by OL had an absolute decline of 11.5%, while the numbers of RL and VL cases saw an absolute increase of 10% and 1.5%, respectively.
- RL had lower complication rates both postoperatively and at 30 days when compared to VL and OL.*
- RL was associated with a reduced LOS when compared to VL and OL.
- Patients in the RL group were more likely to be discharged home than to a transitional health care facility, when compared to VL and OL.
- RL had fewer conversions when compared to VL.
- RL had a lower in-hospital and 30-day mortality rate when compared to OL.

*In comparison with OL, additional reviewed clinical studies report a lower or comparable postoperative complication rate, shorter LOS, and comparable in-hospital mortality.¹⁻⁸ In comparison with VL, additional reviewed clinical studies report a comparable, lower, or higher postoperative complication rate,^{1,4,9-19} comparable or shorter LOS^{1-4, 8-19} and comparable conversion rate.⁸⁻¹⁶

PUBLICATION SUMMARY

REFERENCES

1. Kent, M., T. Wang, et al., *Open, Video-Assisted Thoracic Surgery, and Robotic Lobectomy: Review of a National Database*. Ann Thorac Surg, 2014.
2. Adams, R.D., W.D. Bolton, et al., *Initial Multicenter Community Robotic Lobectomy Experience: Comparisons to a National Database*. Ann Thorac Surg, 2014.
3. Deen, S.A., J.L. Wilson, et al., *Defining the Cost of Care for Lobectomy and Segmentectomy: A Comparison of Open, Video-Assisted Thoracoscopic, and Robotic Approaches*. Ann Thorac Surg, 2014.
4. Farivar, A.S., R.J. Cerfolio, et al., *Comparing Robotic Lung Resection With Thoracotomy and Video-Assisted Thoracoscopic Surgery Cases Entered Into The Society of Thoracic Surgeons Database*. Innovations (Phila), 2014. 9(1): p. 10-5.
5. Cerfolio, R.J., A.S. Bryant, et al., *Initial consecutive experience of completely portal robotic pulmonary resection with 4 arms*. Journal of Thoracic and Cardiovascular Surgery, 2011.
6. Veronesi, G.G., D.; Maisonneuve, P.; Melfi, F.; Schmid, R. A.; Borri, A.; Vannucci, F.; Spaggiari, L., *Four-arm robotic lobectomy for the treatment of early-stage lung cancer*. Journal of Thoracic and Cardiovascular Surgery, 2010. 140(1): p. 19-25.
7. Oh, D.S., I. Cho, et al., *Early adoption of robotic pulmonary lobectomy: feasibility and initial outcomes*. Am Surg, 2013. 79(10): p. 1075-80.
8. 510(k) submission for Thoracic Surgery K153276
9. Augustin, F., J. Bodner, et al., *Robotic-assisted minimally invasive vs. thoracoscopic lung lobectomy: comparison of perioperative results in a learning curve setting*. Langenbecks Arch Surg, 2013.
10. Bodner, J., T. Schmid, et al., *Minimally invasive approaches for lung lobectomy - From VATS to robotic and back!* European Surgery - Acta Chirurgica Austriaca, 2011. 43(4): p. 224-228.
11. Jang, H.J.L., H. S.; Park, S. Y.; Zo, J. I., *Comparison of the early robot-assisted lobectomy experience to video-assisted thoracic surgery lobectomy for lung cancer: A single-institution case series matching study*. Innovations: Technology and Techniques in Cardiothoracic and Vascular Surgery, 2011. 6(5): p. 305-310.
12. Lee, B.E., R.J. Korst, et al., *Transitioning from video-assisted thoracic surgical lobectomy to robotics for lung cancer: Are there outcomes advantages?* J Thorac Cardiovasc Surg, 2013.
13. Lee, B.E., M. Shapiro, et al., *Nodal Upstaging in Robotic and Video Assisted Thoracic Surgery Lobectomy for Clinical N0 Lung Cancer*. Annals of Thoracic Surgery, 2015.
14. Louie, B.E.F., A. S.; Aye, R. W.; Vallieres, E., *Early Experience With Robotic Lung Resection Results in Similar Operative Outcomes and Morbidity When Compared With Matched Video-Assisted Thoracoscopic Surgery Cases*. Annals of Thoracic Surgery, 2012.
15. Mahieu, J., P. Rinieri, et al., *Robot-Assisted Thoracoscopic Surgery versus Video-Assisted Thoracoscopic Surgery for Lung Lobectomy: Can a Robotic Approach Improve Short-Term Outcomes and Operative Safety?* Thoracic and Cardiovascular Surgeon, 2015.
16. Paul, S., J. Jalbert, et al., *Comparative Effectiveness of Robotic-Assisted vs. Thoracoscopic Lobectomy*. Chest, 2014.
17. Spillane, J. and P. Brooks, *Developing a robotic program in thoracic surgery at Cape Cod Hospital*. Journal of Robotic Surgery, 2014: p. 1-8.
18. Lee, H.S.J., H. J., *Thoracoscopic mediastinal lymph node dissection for lung cancer*. Semin Thorac Cardiovasc Surg, 2012. 24(2): p. 131-41.
19. Swanson, S.J., D.L. Miller, et al., *Comparing robot-assisted thoracic surgical lobectomy with conventional video-assisted thoracic surgical lobectomy and wedge resection: Results from a multihospital database (Premier)*. Journal of Thoracic and Cardiovascular Surgery, 2013.

FINANCIAL DISCLOSURE

Dr. Oh and Dr. Reddy have received compensation from Intuitive Surgical for consulting and/or educational services. Authors Gorrepati and Mehendale disclose a financial relationship with Intuitive Surgical.

IMPORTANT SAFETY INFORMATION

Serious complications may occur in any surgery, including da Vinci® Surgery, up to and including death. Examples of serious or life-threatening complications, which may require prolonged and/or unexpected hospitalization and/or reoperation, include but are not limited to one or more of the following: injury to tissues/organs, bleeding, infection and internal scarring that can cause long-lasting dysfunction/pain. Individual surgical results may vary.

Risks specific to minimally invasive surgery, including da Vinci® Surgery, include but are not limited to, one or more of the following: temporary pain/nerve injury associated with positioning; a longer operative time, the need to convert to an open approach, or the need for additional or larger incision sites. Converting the procedure could result in a longer operative time, a longer time under anesthesia, and could lead to increased complications. Contraindications applicable to the use of conventional endoscopic instruments also apply to the use of all da Vinci

instruments. You should discuss your surgical experience and review these and all risks with your patients, including the potential for human error and equipment failure. Physicians should review all available information. Clinical studies are available through the National Library of Medicine at www.ncbi.nlm.nih.gov/pubmed.

Be sure to read and understand all information in the applicable user manuals, including full cautions and warnings, before using da Vinci products. Failure to properly follow all instructions may lead to injury and result in improper functioning of the device. Training provided by Intuitive Surgical is limited to the use of its products and does not replace the necessary medical training and experience required to perform surgery. Procedure descriptions are developed with, reviewed and approved by independent surgeons. Other surgical techniques may be documented in publications available at the National Library of Medicine. For Important Safety Information, indications for use, risks, full cautions and warnings, please also refer to www.davincisurgery.com/safety and www.intuitivesurgical.com/safety. Unless otherwise noted, products featured are available for commercial distribution in the U.S. For availability outside the U.S., please check with your local representative or distributor.

PRECAUTION

The demonstration of safety and effectiveness for the specific procedure(s) discussed in this material was based on evaluation of the device as a surgical tool and did not include evaluation of outcomes related to the treatment of cancer (overall survival, disease-free survival, local recurrence) or treatment of the patient's underlying disease/condition. Device usage in all surgical procedures should be guided by the clinical judgment of an adequately trained surgeon.

The friable nature of pulmonary tissue enhances the risk of vascular, bronchiolar or other injury that will be difficult to control when using this device. Published clinical experience as well as clinical studies performed to support this marketing clearance have demonstrated that even surgeons considered expert in laparoscopy/thoracoscopy have substantial learning curves of 10 to 12 cases (Falk, et al., Total endoscopic computer enhanced coronary artery bypass grafting, *Eur J Cardiothorac Surg* 2000; 17: 38-45).

All materials will eventually become obsolete. When referencing printed or digitally replicated materials, please note the revision date that follows the part number (PN). Consult your da Vinci representative or visit the da Vinci Online Community for the latest revision.