
Advancements in Gynecologic Surgery

Presenter's Name

Agenda

- Advancements in GYN Surgery
- What Procedures Do I Offer?
- Potential Patient Benefits & Risks
- Clinical Evidence
- Working Together

My Credentials

- Board certified OB/GYN since ____
 - [current title, fellowship training, hospital]
- Practicing robotic-assisted surgery since ____
- Completed ____ cases:
 - *da Vinci*® *Single-Site*® Hysterectomy ____
 - *da Vinci* multi-port procedures ____

[Please edit with your information and add the *da Vinci* procedures that you perform the most here]

Our Collective Goal: Patients First

- Care for patients' overall, long-term health
 - Fibroids: 1 in 4 women¹
 - Endometriosis: 5+ million²
 - Pelvic prolapse: 1 in 3 women with 1+ children³
 - GYN cancer: 80,000 new cases/year⁴
- Ensure patients are aware of all options for their conditions
 - Lifestyle changes
 - Medications
 - Hormone Therapy
 - Surgical options (e.g. hysterectomy by age 60: 1 in 3 women⁵ - 90% benign⁶)
- Provide the surgical options most appropriate for each patient

1. U.S. Department of Health and Human Services, Office of Women's Health. The Healthy Woman: A Complete Guide for All Ages.


2. Office on Women's Health. <http://www.womenshealth.gov/publications/our-publications/fact-sheet/endometriosis.pdf> 3. International Urogynecologic Association. Pelvic Organ Prolapse: A Guide for Women. 2011 4. CDC. http://www.cdc.gov/cancer/gynecologic/basic_info/index.htm

5. NIH. <http://www.nlm.nih.gov/medlineplus/hysterectomy.html> 6. American College of Surgeons. About Hysterectomy: Surgical removal of the uterus, or womb.


Advancement Towards Less Invasive Surgery



Abdominal
(Open) Surgery




Vaginal Surgery



Robotic-Assisted
Surgery: Multi-port




Laparoscopy



Robotic-Assisted
Surgery: *Single-Site*[®]


5 *Single-Site*[®] Instruments for the *da Vinci*[®] *Si*[™] System bear the CE mark. This device is cleared for commercial distribution in the U.S. for laparoscopic cholecystectomy and hysterectomy and salpingo-oophorectomy for benign conditions.

Incision Comparison




One Large Incision

Open Hysterectomy




Multiple Small Incisions

Traditional Laparoscopic Hysterectomy



da Vinci[®] Hysterectomy:
Multi-port



One Small Incision in the Belly Button

da Vinci[®] *Single-Site*[®] Hysterectomy
or
Traditional Laparoscopic Single Incision Hysterectomy

Why Have I Adopted *da Vinci*[®] Surgery?


- Have been skilled in Vaginal and Lap surgery since

- XX% of my cases were still done Abdominally
 - [please edit reasons why you couldn't perform more MIS previously]
 - [e.g. cases too complex for lap or vaginal surgery]
 - [e.g. hard to access or suture with straight lap instruments]
 - [e.g. “shaky” camera view during lap surgery]
 - [.....]

Why Have I Adopted *da Vinci*® Surgery? (Con't)

- I want to offer minimally invasive surgery to my patients

[please double click the graph below to edit with your own data]




Compared to U.S. Adoption of MIS - Malignant Hyst.

U.S. MALIGNANT HYSTERECTOMY MARKET BY MODALITY

Estimated Adoption of Minimally Invasive Surgery (MIS)

Percentage of all procedures



FDA clearance of *da Vinci*® Surgery GYN, 2005

Rapid decline

Rapid increase

IMPACT OF ROBOTIC-ASSISTED SURGERY:

Prior to the introduction of robotics, most hysterectomies were performed using a large incision (open surgery). The overall rate of minimally invasive surgery (vaginal and laparoscopy) remained relatively unchanged.

Following the introduction of robotics (*da Vinci* Surgery), the rate of open surgery began to rapidly decline, while the rate of minimally invasive surgery (vaginal, laparoscopy and *da Vinci* Surgery) began to rise.

Today, open surgery is used in only about 20% of hysterectomies.

Intuitive Surgical, Inc.
PN 873930 rev F 4/2014


1. Inpatient data: Nationwide Inpatient Sample (NIS), Healthcare Cost and Utilization Project (HCUP), Agency for Healthcare Research and Quality
2. Outpatient data: Solucient® Database - Truven Health Analytics (Formerly Thomson-Reuters) 3. *da Vinci* data: ISI Internal Estimates

Compared to U.S. Adoption of MIS - Benign Hyst.

U.S. BENIGN HYSTERECTOMY MARKET BY MODALITY

Estimated Adoption of Minimally Invasive Surgery (MIS)

Percentage of all procedures



FDA clearance of
da Vinci® Surgery
GYN, 2005

IMPACT OF ROBOTIC-ASSISTED SURGERY:

Prior to the introduction of robotics, most hysterectomies were performed using a large incision (open surgery). The overall rate of minimally invasive surgery (vaginal and laparoscopy) remained relatively unchanged.

Following the introduction of robotics (*da Vinci* Surgery), the rate of open surgery began to rapidly decline, while the rate of minimally invasive surgery (vaginal, laparoscopy and *da Vinci* Surgery) began to rise.

Today, open surgery is used in only about 20% of hysterectomies.

Intuitive Surgical, Inc.
PN 873931 rev F 4/2014

1. Inpatient data: Nationwide Inpatient Sample (NIS), Healthcare Cost and Utilization Project (HCUP), Agency for Healthcare Research and Quality
2. Outpatient data: Solucient® Database - Truven Health Analytics (Formerly Thomson-Reuters) 3. *da Vinci* data: ISI Internal Estimates

da Vinci[®] Surgery: How Does It Work?

da Vinci[®] Surgery: Over 3 Million Procedures Worldwide



High Definition 3D Vision

- Surgeon-controlled
- Stable and immersive view
- Up to 10x zoom

Precision & Dexterity


- Mimics surgeon's hands
- Scales down movements
- With tremor filtration

Intuitive[®] Motion

- Advanced software enables intuitive control (instead of cross-handed)

da Vinci® Single-Site® Surgery

- Curved cannula & semi-rigid instruments → Access from many angles & instrument triangulation
- Designed to limit cannula shift → Minimize potential port-site trauma & post-op pain^{1,2}
- A single incision in the belly button → Virtually scarless^{1,3}



Results, including cosmetic results, may vary.

1. Kroh M, *et al.* First human surgery with a novel single-port robotic system: cholecystectomy using the *da Vinci Single-Site* Platform. *Surg Endosc.* 2011 Nov;25(11):3566-73. Epub 2011 Jun 3rd. 2. Wren SM, Curet MJ. Single-port robotic cholecystectomy results from a first human use clinical study of the new *da Vinci Single-Site* surgical platform. *Arch Surg*, Jun 2011; doi: 10.1001/archsurg.2011.143 3. Cela V, Freschi L, Simi G, Ruggiero M, Tana R, Pluchino N. Robotic single-site hysterectomy: feasibility, learning curve and surgical outcome. *Surg Endosc.* 2013 Jul;27(7):2638-43. doi: 10.1007/s00464-012-2780-8.

da Vinci[®] Surgical Procedures That I Offer

I Offer These *da Vinci*[®] Surgical Procedures

[please edit to show the specific *da Vinci* procedure you offer]

Fibroids /
Abnormal
Uterine Bleeding

Hysterectomy
(*Single-Site*[®] or Multi-port)

Myomectomy
(uterine preservation)



Endometriosis

Endometriosis
Resection
(with or without a
Hysterectomy)



Pelvic Organ
Prolapse

Sacrocolpopexy
(with or without a
Hysterectomy)



GYN Cancer

Hysterectomy
with
Lymphadenectomy
(Multi-port only)



Patient Selection for Referring GYN Patients to My Practice

Patient Selection:

Nearly every GYN patient in need of surgery, with the following conditions:

Conditions:

- Benign Pathology
 - Fibroids
 - Endometriosis
 - Pelvic Prolapse
 - Abnormal Uterine Bleeding
- Cancers
 - Cervical
 - Endometrial (Uterine)
 - Ovarian

Not Limited by:

- Incidence of prior surgeries
- Prior pelvic infection
- Large uteri
- Presence of single or multiple adnexal masses

Patient Case: Hysterectomy [Template]


Patient History	<ul style="list-style-type: none">○ [Sex and Age]○ [Diagnosis]○ [Condition]○ [Symptoms]
Non-Surgical options	<ul style="list-style-type: none">○ [Has your patient tried non-surgical options]○ [Has your patient taken OTC pain medication to relieve pain when needed]
Patient Concerns	<ul style="list-style-type: none">○ [Concern 1]○ [Concern 2]○ [Concern 3]
Post-op	<ul style="list-style-type: none">○ Hospital LOS: _____ hours/days○ Return to normal activities: in _____ days○ Currently: _____ (e.g. disease and complication free)

[Please add/replace content in brackets with your own patient case data]

da Vinci[®]
Hysterectomy

Potential Patient Benefits & Risks

Potential Patient Benefits & Risks *da Vinci*[®] Hysterectomy (Benign)



da Vinci Hysterectomy (Benign) Potential Patient Benefits

vs. Open Surgery

- Reduced complication rate^{1, 2, 3, 9}
- Reduced length of hospital stay^{1, 2, 3, 4, 9}
- Reduced blood loss and less likelihood for transfusion^{1, 3, 4, 9}
- Reduced readmission rate^{4,9}

vs. Traditional Laparoscopic Surgery

- Reduced complication rate^{1, 5, 9}
- Reduced length of hospital stay^{1, 2, 4, 5, 6, 7, 9}
- Reduced blood loss^{1, 2, 4, 7}
- Reduced chance of procedure converting to an abdominal procedure^{2, 5}
- Less likelihood of blood transfusion^{8, 9}

vs. Vaginal Surgery


- Reduced length of hospital stay^{2, 4, 9}
- Reduced blood loss^{2, 4}

Hysterectomy, incl. *da Vinci* Hyst. Potential Patient Risks

- Injury to the ureters (the ureters drain urine from the kidney into the bladder)
- Vaginal cuff problems (scar tissue in vaginal incision, infection, bacterial skin infection, pooling/clotting of blood, incision opens or separates)
- Injury to bladder (organ that holds urine), bowel injury, vaginal shortening, problems urinating (cannot empty bladder, urgent or frequent need to urinate, leaking urine, slow or weak stream)
- Abnormal hole from the vagina into the urinary tract or rectum, vaginal tear or deep cut.
- Uterine tissue may contain unsuspected cancer. The cutting or morcellation of uterine tissue during surgery may spread cancer, and decrease the long-term survival of patients.

Potential Patient Benefits & Risks

da Vinci® *Single-Site*® Benign Hysterectomy



da Vinci Single-Site Benign Hysterectomy

Potential Patient Benefits

Early clinical data suggests:


- Low blood loss^{10,11,12, 13}
- Low complication rate^{10, 11}
- Low likelihood of blood transfusion¹¹
- Low chance of procedure converting to an abdominal procedure^{11,12,13}
- Short hospital stay^{11,12}
- Low post-operative pain¹³

Hysterectomy, incl. *da Vinci* Hyst.

Potential Patient Risks

- Injury to the ureters (the ureters drain urine from the kidney into the bladder)
- Vaginal cuff problem (replaces cervix): scar tissue in vaginal incision, infection, bacterial skin infection, pooling/clotting of blood, incision opens or separates
- Injury to bladder (organ that holds urine), bowel injury, vaginal shortening, problems urinating (cannot empty bladder, urgent or frequent need to urinate, leaking urine, slow or weak stream)
- Abnormal hole from the vagina into the urinary tract or rectum, vaginal tear or deep cut.
- Uterine tissue may contain unsuspected cancer. The cutting or morcellation of uterine tissue during surgery may spread cancer, and decrease the long-term survival of patients .

Potential Patient Benefits & Risks *da Vinci*[®] Myomectomy



da Vinci Myomectomy Potential Patient Benefits

vs. Open Surgery

- Similar rate of complications¹⁴
- Shorter hospital stay^{14,15,16}
- Lower rate of blood transfusions^{14,16}
- Less blood loss^{14,15,16}
- Less chance of post-operative fever¹⁴

vs. Traditional Laparoscopic Surgery


- Similar rate of complications^{14,16,17}
- Similar hospital stay^{14,16,17}
- Similar conversion rate (switch to open surgery)^{14,15,17}
- Similar or less blood loss^{14,15,17}

Myomectomy, incl. *da Vinci* Myomectomy Potential Patient Risks

- Tear or hole in uterus
- Split or bursting of the uterus
- Pre-term (early) birth
- Spontaneous abortion
- Uterine tissue may contain unsuspected cancer. The cutting or morcellation of uterine or fibroid tissue during surgery may spread cancer, and decrease the long-term survival of patients .

Potential Patient Benefits & Risks *da Vinci*[®] Sacrocolpopexy

Sacrocolpopexy Incisions



da Vinci Sacrocolpopexy Potential Patient Benefits

vs. Open Surgery

- Lower rates of complications^{18,19}
- Reduced estimated blood loss^{18,19,20}
- Comparable rates of blood transfusions¹⁸
- Shorter length of hospital stay¹⁸

vs. Traditional Laparoscopic Surgery


- Comparable rates of complications^{18,19,21,22}
- Comparable rates of blood transfusions^{18,21}
- Comparable or shorter length of hospital stay^{18,19,22,23}
- Comparable conversion rate^{18,19,23}
- Comparable or reduced estimated blood loss^{18,19,21}

Sacrocolpopexy, incl. *da Vinci* Sacrocolpopexy Potential Patient Risks

- Mesh erosion/infection caused by mesh moving from vaginal wall into surrounding organs causing the need for another operation
- Injury to rectum/bowel, injury to bladder (organ that holds urine), injury to the ureters (the ureters drain urine from the kidney into the bladder)
- Front wall of the rectum pushes into the back wall of the vagina
- Prolapsed bladder (bladder budes into vagina when supportive tissue weakens)
- Vaginal incision opens or separates, loss of bladder control, pooling of blood between bladder and pubic bone, pooling of blood between the anus and vagina.

Potential Patient Benefits & Risks *da Vinci*® Endometriosis Resection

Endometriosis Resection Incisions



Traditional Lap
Surgery

da Vinci Surgery

da Vinci Endometriosis Resection Potential Patient Benefits


Early Clinical Data Suggest:

- Ability for surgeon to complete difficult dissections (separating of tissue)^{24,25,26}
- Low rate of complications^{24,26,27,28}
- Low blood loss^{26,28,30,31} and low chance for blood transfusion^{26,28}
- Low rate of switching to open surgery (through a large incision)^{24,26,29,30}

Endometriosis Resection, incl. *da Vinci* Endo Resection Potential Patient Risks

- Injury to the bowel, bladder (organ that holds urine) or ureters (the ureters drain urine from the kidney into the bladder)

Potential Patient Benefits & Risks *da Vinci*[®] Hysterectomy (Cancer)



da Vinci Hysterectomy (Simple Total) for Cancer

Potential Patient Benefits

vs. Open Surgery

- Fewer complications ^{32,33,34,35}
- Fewer blood transfusions ^{32,33,34,35}
- Shorter length of stay ^{32,33,34,35,36}
- Reduced estimated blood loss ^{32,33,34,35,36}

vs. Traditional Laparoscopic Surgery

- Comparable or fewer complications ^{32,33,34,35,36}
- Comparable or fewer blood transfusions ^{32,33,34,35,36}
- Comparable or shorter length of stay ^{32,33,34,35,36}
- Comparable operative time ^{32,33,34,35,36}
- Comparable or lower conversion rates ^{32,33,34,35,36}
- Reduced estimated blood loss ^{32,33,34,35,36}

Hysterectomy, incl. *da Vinci* Hyst. Potential Patient Risks

- Injury to the ureters (the ureters drain urine from the kidney into the bladder)
- Vaginal cuff problem (replaces cervix): scar tissue in vaginal incision, infection, bacterial skin infection, pooling/clotting of blood, incision opens or separates
- Injury to bladder (organ that holds urine), bowel injury, vaginal shortening, problems urinating (cannot empty bladder, urgent or frequent need to urinate, leaking urine, slow or weak stream)
- Abnormal hole from the vagina into the urinary tract or rectum, vaginal tear or deep cut

Clinical Evidence



Multicenter analysis comparing robotic, open, laparoscopic, and vaginal hysterectomies performed by high-volume surgeons for benign indications

Peter C. Lim, John T. Crane, Eric J. English, Richard W. Farnam, Devin M. Garza, Marc L. Winter, Jerry L. Rozeboom

Robotic patients had:

- Larger uteri
- Fewer complications
- Shorter hospital stay
- Similar conversions, readmissions & reoperation rate as other MIS approaches

	Robotic-Assisted (n = 2300)	Abdominal (n = 9745)	Vaginal (n = 8121)	Laparoscopic (n = 11952)
Age	49.3 ± 11.5	46.7 ± 10.7*	48.7 ± 13.3*	43.9 ± 9.4*
Large uterus (>250g)	366 (15.9%)	368 (3.8%)*	589 (7.3%)*	1,671 (14%)*
Intraoperative complications	17 (0.7%)	174 (1.8%)*	142 (1.8%)*	142 (1.2%)
Postoperative complications	131 (6.3%)	2,047 (21%)*	1,314 (16.2%)*	1,953 (16.3%)*
Conversion to open surgery	2 (0.1%)	NA	1 (0.0%)	11 (0.1%)
Inpatient length of stay (days)	1.37 ± 1.1	3.0 ± 1.6*	1.9 ± 1.0*	1.7 ± 1.2*
Hospital readmission (related to index surgery)	28 (1.3%)	340 (3.5%)*	130 (1.6%)	186 (1.6%)
Reoperation rate	17 (0.8%)	187 (1.9%)*	84 (1.0%)	118 (1.0%)

In this study, two sided $p < 0.05$ was considered significant. * Indicates statistically significant difference ($p < 0.05$) versus robotic-assisted hysterectomy.

Study limitations: Retrospective data review—missing data is a common, inherent limitation of retrospective data collection; The Premier database relies on ICD-9-CM diagnostic and procedure codes; there is a potential for miscoding.; If patients were readmitted to non-Premier hospitals, patient readmission data could have been lost, which raises the potential for under-reporting. The length of stay for outpatients was not available in the Premier database.

Financial disclosure: This study was funded by Intuitive Surgical for independent research and editorial support. Dr. Lim, Dr. Crane, Dr. English, Dr. Farnam, Dr. Garza, Dr. Winter, and Dr. Rozeboom have received compensation from Intuitive Surgical for consulting and/or educational services.

A Comparison of Quality Outcome Measures in Patients Having a Hysterectomy for Benign Disease: Robotic vs. Non-robotic Approaches

Martin A. Martino, MD*, Elizabeth A. Berger, DO, Jeffrey T. McFetridge, MS, Jocelyn Shubella, BS, Gabrielle Gosciniak, BA, Taylor Wejkszner, BA, Gregory F. Kainz, DO, Jeremy Patriarco, BS, M. Bijoy Thomas, MD, and Richard Boulay, MD



ARTICLE IN PRESS

First published online 28 October 2013

	Robotic-assisted (N=601)	Laparoscopic (N=427)	Abdominal (Open) (N=1,194)	Vaginal (N=332)	<i>p</i> -Value (non-robotic vs. robotic-assisted)
Mean Estimated Blood Loss (mL)	108.2	315.08	318.8	340.8	<0.05
Mean Length of Stay (minutes)	1570.3	3038.5	3440.5	3789.2	<0.05
Total Readmissions <30 days (Total No.)	1.00% (6)	2.58% (11)	3.52% (42)	2.41% (8)	≤0.03
Total Readmission Cost (adjusted for inflation to 2012 \$)	\$32,946	\$50,290	\$328,230	\$51,264	<0.05

Compared to non-robotic approaches:

- Less blood loss
- Shorter hospital stay
- Lower rate of <30-day readmission
- Significant readmission related cost savings

Study limitations: Retrospective data review, single institution experience, inability to account for patients who were readmitted to outside hospitals (though the authors state that this finding was expected to be similar for all 4 cohorts)

Financial disclosure: "Dr. Martino has received travel reimbursement from Intuitive Surgical for educational research."

Outcomes of Abdominal and Minimally Invasive Sacrocolpopexy: A Retrospective Cohort Study



Patrick A. Nosti, MD,* Uduak Umoh Andy, MD,† Sarah Kane, MD,‡ Dena E. White, MD,§ Heidi S. Harvie, MD, MBA, MSCE,† Lior Lowenstein, MD, MS,|| and Robert E. Gutman, MD*

ASC=abdominal sacrocolpopexy, MISC=minimally invasive sacrocolpopexy, LSC=laparoscopic sacrocolpopexy, RSC= robotic sacrocolpopexy

Intraoperative & Postoperative Complications	ASC (n=589)	MISC (n=535)	P-Value	LSC (n=273)	RSC (n=262)	P-Value
All Complications	20%	12.7%	<0.01	18%	7%	<0.01
<i>Cystotomy</i>	4%	2%	<0.01	2.5%	1.5%	0.7
<i>DVT/PE</i>	1.5%	1%	0.3	3%	0.0%	<0.01
<i>Ileus/SBO</i>	5%	2%	<0.01	1.8%	1.5%	1
<i>Conversion to Open</i>	-	-	-	4%	0.4%	<0.01
Anatomical Failure						
Anatomical Failure	ASC	MISC	P-Value	LSC	RSC	P-Value
At or beyond hymen	15.1%	7.4%	<0.001	6.5%	8.4%	0.49
Stage 2 or higher	25.3%	14.2%	<0.001	11.3%	17.2%	0.069

The MISC group had shorter hospitalization, less blood loss, but longer operative times compared with ASC group

LSC was associated with more complications compared with RSC. Patients who underwent LSC were more likely to have their procedure converted to open

Study limitations: Retrospective study; study investigators used a composite complication score which may not account for all adverse events; definition of failure did not include subjective data which were inconsistently available

Financial disclosure: None

Robotic-Assisted, Laparoscopic, and Abdominal Myomectomy: A Comparison of Surgical Outcomes

Ehab E. Barakat, MD, Mohamed A. Bedaiwy, MD, Stephen Zimberg, MD, Benjamin Nutter, Mohsen Nosseir, MD, and Tommaso Falcone, MD

FEBRUARY 2011



	Abdominal (N=393)	Lap (N=93)	Robotic (N=89)	P-Value
Surgical Time (min)	126	155	181	Open vs. Robotic=.003 Lap vs. Robotic=.083
Myoma Weight (g)	263.00	96.65	223.00	Open vs. Robotic=.36 Lap vs. Robotic=.021
Estimated Blood Loss (mL)	200	150	100	Open vs. Robotic<.001 Lap vs. Robotic=.818
Hemoglobin Drop (g/dL)	2.00	1.55	1.30	Open vs. Robotic<.001 Lap vs. Robotic=.431
Length of Hospital Stay (days)	3	1	1	Open vs. Robotic<.001 Lap vs. Robotic=.506

Compared to Lap:

- Bigger myomas can be removed

Compared to open surgery:

- Less blood loss
- Less hemoglobin drop
- Shorter hospital stay

Study limitations: Retrospective, Lack of long-term outcomes
Financial disclosure: None.

Retrospective analysis of robot-assisted versus standard laparoscopy in the treatment of pelvic pain indicative of endometriosis

John F. Dulemba · Cyndi Pelzel · Helen B. Hubert

J Robotic Surg (2013) 7:163–169

DOI 10.1007/s11701-012-0361-4



	Robot-assisted	Laparoscopy	p Value
Operative time (minutes)	77.4	72.3	0.23
Patients w/ confirmed endometriosis	80.0%	56.8%	<0.001
Appendectomy	23.3%	30.0%	0.32
Appendix positive for endometriosis	28.6%	3.3%	0.02

Histological confirmation of endometriosis was significantly higher for the robotic cohort, suggesting robot-assisted techniques may provide more accurate visualization and, thus, excision of existing endometriosis

Study limitations: Retrospective, single surgeon experience, lack of validated and longer-term outcome measures

Financial disclosure: Dr. Dulemba has received compensation from Intuitive Surgical for consulting and/or educational services.

Dulemba, John F., Cyndi Pelzel, and Helen B. Hubert. "Retrospective Analysis of Robot-assisted versus Standard Laparoscopy in the Treatment of Pelvic Pain Indicative of Endometriosis." Journal of Robotic Surgery 7.2 (2013): 163-69. Print.

A comparative study of 3 surgical methods for hysterectomy with staging for endometrial cancer: robotic assistance, laparoscopy, laparotomy

John F. Boggess, MD; Paola A. Gehrig, MD; Leigh Cantrell, MD; Aaron Shafer, MD; Mildred Ridgway, MD; Elizabeth N. Skinner, MD; Wesley C. Fowler, MD

OCTOBER 2008



	Open (n=138)	Lap (n=81)	Robotic (n=103)	p-Value
Mean age (years)	64.0	62.0	61.9	Open vs. dV=.06 Lap vs. dV=.95
Mean BMI (kg/m ²)	34.7	29.0	32.9	Open vs. dV=.17 Lap vs. dV=.0008
Mean EBL (ml)	266.0	145.8	74.5	Open vs. dV<.0001 Lap vs. dV<.0001
Mean operative time (skin-to-skin) (min)	146.5	213.4	191.2	Open vs. dV<.0001 Lap vs. dV<.0001
Mean total lymph nodes (n)	14.9	23.1	32.9	Open vs. dV<.0001 Lap vs. dV<.0001
Mean hospital stay (days)	4.4	1.2	1.0	Open vs. dV<.0001 Lap vs. dV=.001
Complications	29.7%	13.6%	5.8%	Open vs. dV<.0001 Lap vs. dV<.0001
Conversions to Open	---	4.9%	2.9%	Lap vs. dV=.7

Compared to open surgery:

- Fewer complications


Compared to open and lap surgery:

- Less blood loss
- Higher lymph node yield
- Shorter hospital stay

Study limitations discussed by the authors: "it was not randomized; because of our relatively recent incorporation of robotic technology, we could not examine long-term oncologic results."

Financial disclosure: None were included in the publication.

What Can We Do Together?



My Contact Information:

[Surgeon's Name]

[Address]

[Phone]

[Email]

[Contact Name for managing referrals]

[Website]

Important Safety Information

Serious complications may occur in any surgery, including *da Vinci*[®] Surgery, up to and including death. Examples of serious and life-threatening complications, which may require prolonged and/or unexpected hospitalization and/or reoperation, include, but are not limited to, one or more of the following:

- Injury to tissues and/or organs
- Bleeding
- Infection
- Internal scarring that can cause long-lasting dysfunction or pain

Patients should consider that risks of any surgery also include, but are not limited to, the following:

- Potential for human error
- Potential for equipment failure
- Potential for anesthesia complications

Individual surgical results may vary.

Risks specific to minimally invasive surgery, including *da Vinci* Surgery, include, but are not

Important Safety Information

limited to, the following:

- Temporary pain and/or nerve injury associated with positioning;
- Temporary pain and/or discomfort from the use of air or gas in the procedure;
- A longer operative time and time under anesthesia;
- The need to convert the procedure to an open surgery;
- Converting the procedure could result in a longer operative time, a longer time under anesthesia, and/or the need for additional or larger incisions and/or increased complications.

Medical Advice & Surgeon Training

Patients should talk to their doctor to decide if *da Vinci* Surgery is right for them. Other options may be available and appropriate. Only a doctor can determine whether *da Vinci* Surgery is appropriate for a patient's situation. Patients and doctors should review all available information on both non-surgical and surgical options in order to make an informed decision.

Only surgeons who have received specific training in the use of the *da Vinci* Surgical System should use the system. Training provided by Intuitive Surgical is limited to the use of its products and does not replace the necessary medical training and experience required to

Important Safety Information

perform surgery.

***da Vinci* System Description**

There are several models of the *da Vinci* Surgical System. The *da Vinci* Surgical Systems are designed to help doctors perform minimally invasive surgery. The *da Vinci* Surgical System is not programmed to perform surgery on its own. Instead, the surgery is performed entirely by a doctor, who controls the system. *da Vinci* Systems offer doctors high-definition 3D vision, a magnified view, and robotic and computer assistance. They use specialized instrumentation, including a miniaturized surgical camera and wristed instruments (i.e., scissors, scalpels and forceps) – that are designed to help with precise dissection and reconstruction deep inside the body.

When is the *da Vinci* System Used?

One or more of the *da Vinci* Surgical System models are commercially available for use in the following specialty areas.

- Urologic surgery
- General laparoscopic surgery
- Gynecologic surgery
- Transoral robotic surgery restricted to benign (non-cancer) and malignant tumors (cancer) classified as T1 and T2 (early stage cancer) and for benign base of tongue

Important Safety Information

- resection procedures. * The safety and effectiveness of this device for use in the options of obstructive sleep apnea have not been established.
- Thoracic surgery
- Heart surgery

**Not cleared for use with the da Vinci Xi Surgical System.*

When the *da Vinci* System Is Not Used?

Patients who are not candidates for non-robotic minimally invasive surgery are also not candidates for *da Vinci* Surgery.

More information about the *da Vinci* System and Locating a Doctor

If you have questions about the *da Vinci* System or *da Vinci* procedures, consult with a qualified surgeon. Surgeons experienced with the *da Vinci* System can be found using the [Surgeon Locator](http://www.davincisurgery.com) at www.davincisurgery.com. Intuitive Surgical provides surgeons training on the use of the *da Vinci* System but does not certify, credential or qualify the surgeons listed in the Surgeon Locator.

When is *Single-Site* Used and What are the Risks?

da Vinci Surgery with *Single-Site*[®] Instruments is cleared for use in gallbladder removal, and for

Important Safety Information

hysterectomy and ovary removal for benign conditions. Patients who are not candidates for non-robotic minimally invasive surgery are also not candidates for *da Vinci* Surgery, including *da Vinci* Surgery with *Single-Site*[®] Instruments. There may be an increased risk of incision-site hernia with single-incision surgery, including Single-Site surgery with *da Vinci*.

Precaution

The demonstration of safety and effectiveness for the specific procedure(s) discussed in this material was based on evaluation of the device as a surgical tool and did not include evaluation of outcomes related to the treatment of cancer (overall survival, disease-free survival, local recurrence) or treatment of the patient's underlying disease/condition. Device usage in all surgical procedures should be guided by the clinical judgment of an adequately trained surgeon.

Spontaneous opinions expressed during live presentations by individual doctors or patients belong to those individuals. These opinions are not necessarily shared by Intuitive Surgical, Inc.

Intraoperative video, including video labeled as showing full-length procedures, may have been edited for content, for length or to meet file-size limitations.

Thank You!

[name]
[practice name]
[phone]
[email]
[website]

References

References:

1. Ho C, Tsakonas E, Tran K, Cimon K, Severn M, Mierzwinski-Urban M, Corcos J, Pautler S. "Robot-Assisted Surgery Compared with Open Surgery and Laparoscopic Surgery: Clinical Effectiveness and Economic Analyses." Ottawa (ON): Canadian Agency for Drugs and Technologies in Health; 2011 Sep.
2. Landeen, Laurie B., MD, MBA, Maria C. Bell, MD, MPH, Helen B. Hubert, MPH, PhD, Larissa Y. Bennis, MD, Siri S. Knutsten-Larsen, MD, and Usha Seshari-Kreaden, MSc. "Clinical and Cost Comparisons for Hysterectomy via Abdominal, Standard Laparoscopic, Vaginal and Robot-assisted Approaches." *South Dakota Medicine* 64.6 (2011): 197-209. Print.
3. Geppert B, Lönnerfors C, Persson J. "Robot-assisted laparoscopic hysterectomy in obese and morbidly obese women: surgical technique and comparison with open surgery." *Acta Obstet Gynecol Scand.* 90.11 (2011): 1210-1217. doi: 10.1111/j.1600-0412.2011.01253.x. Epub.
4. Martino, Martin A., MD, Elizabeth A. Berger, DO, Jeffrey T. McFetridge, MD, Jocelyn Shubella, BS, Gabrielle Gosciniak, BA, Taylor Wejkszner, BA, Gregory F. Kainz, DO, Jeremy Patriarco, BS, M. B. Thomas, MD, and Richard Boulay, MD. "A Comparison of Quality Outcome Measures in Patients Having a Hysterectomy for Benign Disease: Robotic vs. Non-robotic Approaches." *Journal of Minimally Invasive Gynecology* 21.3 (2014): 389-93. Web.
5. Scandola, Michele, Lorenzo Grespan, Marco Vicentini, and Paolo Fiorini. "Robot-Assisted Laparoscopic Hysterectomy vs Traditional Laparoscopic Hysterectomy: Five Metaanalyses." *Journal of Minimally Invasive Gynecology* 18.6 (2011): 705-15. Print.
6. Wright, Jason D., Cande V. Ananth, Sharyn N. Lewin, William M. Burke, Yu-Shiang Lu, Alfred I. Neugut, Thomas J. Herzog, and Dawn L. Hershman. "Robotically Assisted vs Laparoscopic Hysterectomy Among Women With Benign Gynecologic Disease." *Jama* 309.7 (2013): 689-98. Print.
7. Orady, Mona, Alexander Hrynewych, A. Karim Nawfal, and Ganesa Wegienka. "Comparison of Robotic-Assisted Hysterectomy to Other Minimally Invasive Approaches." *JSL, Journal of the Society of Laparoendoscopic Surgeons* 16.4 (2012): 542-48. Print.
8. Rosero, Eric B., Kimberly A. Kho, Girish P. Joshi, Martin Giesecke, and Joseph I. Schaffer. "Comparison of Robotic and Laparoscopic Hysterectomy for Benign Gynecologic Disease." *Obstetrics & Gynecology* 122.4 (2013): 778-86. Print.
9. Lim, Peter C., John T. Crane, Eric J. English, Richard W. Farnam, Devin M. Garza, Marc L. Winter, and Jerry L. Rozeboom. "Multicenter analysis comparing robotic, open, laparoscopic, and vaginal hysterectomies performed by high-volume surgeons for benign indications." *International Journal of Gynecology & Obstetrics* 133.3 (2016): 359-364. Print.
10. Cela V, Freschi L, Simi G, Ruggiero M, Tana R, Pluchino N. "Robotic single-site hysterectomy: feasibility, learning curve and surgical outcome." *Surg Endosc.* 2013 Jul;27(7):2638-43. doi: 10.1007/s00464-012-2780-8. Epub 2013 Feb 8.
11. Akdemir A, Zeybek B, Ozgurel B, Oztekin MK, Sendag F. "Learning curve analysis of intracorporeal cuff suturing during robotic single-site total hysterectomy." *J Minim Invasive Gynecol.* 2015 Mar-Apr;22(3):384-9. doi: 10.1016/j.jmig.2014.06.006. Epub 2014 Jun 19.
12. Scheib SA, Fader AN. "Gynecologic robotic laparoendoscopic single-site surgery: prospective analysis of feasibility, safety, and technique." *Am J Obstet Gynecol.* 2015 Feb;212(2):179.e1-8. doi: 10.1016/j.ajog.2014.07.057. Epub 2014 Aug 1.
13. Bogliolo S, Mereu L, Cassani C, Gardella B, Zanellini F, Dominoni M, Babilonti L, Delpezzo C, Tateo S, Spinillo A. "Robotic single-site hysterectomy: two institutions' preliminary experience." *Int J Med Robot.* 2014 Sep 18. doi: 10.1002/rcs.1613. [Epub ahead of print]
14. Pundir, Jyotsna, Vishal Pundir, Rajalaxmi Walavalkar, Kireki Omanwa, Gillian Lancaster, and Salma Kayani. "Robotic-Assisted Laparoscopic vs Abdominal and Laparoscopic Myomectomy: Systematic Review and Meta-Analysis." *Journal of Minimally Invasive Gynecology* 20.3 (2013): 335-45. Print.
15. Reza, M., S. Maeso, J. A. Blasco, and E. Andradas. "Meta-analysis of Observational Studies on the Safety and Effectiveness of Robotic Gynaecological Surgery." *British Journal of Surgery* 97.12 (2010): 1772-783. Print.
16. Govern, Joseph M., C. J. Rosemeyer, James F. Barter, and Albert J. Steren. "Comparison of Robotic, Laparoscopic, and Abdominal Myomectomy in a Community Hospital." *JSL, Journal of the Society of Laparoendoscopic Surgeons* 17.1 (2013): 116-20. Print.

References, cont.

References:

17. Pluchino, Nicola, Piero Litta, Letizia Freschi, Marinella Russo, Giovanna Simi, Anna N. Santoro, Stefano Angioni, Angiolo Gadducci, and Vito Cela. "Comparison of the Initial Surgical Experience with Robotic and Laparoscopic Myomectomy." *The International Journal of Medical Robotics and Computer Assisted Surgery* (2013): N/a. Web.
18. Serati, et al (2014). Robot-assisted Sacrocolpopexy for Pelvic Organ Prolapse: A Systematic Review and Meta-Analysis of Comparative Studies. *European Urology*. 66:202-318.
19. Nosti, Patrick A., Uduak Umoh Andy, Sarah Kane, Dena E. White, Heidi S. Harvie, Lior Lowenstein, and Robert E. Gutman. "Outcomes of Abdominal and Minimally Invasive Sacrocolpopexy." *Female Pelvic Medicine & Reconstructive Surgery* 20.1 (2014): 33-37. Print.
20. Siddiqui, Nazema Y., Elizabeth J. Geller, and Anthony G. Visco. "Symptomatic and Anatomic 1-year Outcomes after Robotic and Abdominal Sacrocolpopexy." *American Journal of Obstetrics and Gynecology* 206.5 (2012): 435.e1-35.e5. Print.
21. Anger, et al (2014). Robotic Compared with Laparoscopic Sacrocolpopexy. A Randomized Controlled Trial. *Obstetrics & Gynecology*. 123(1)5-12.
22. Flack, et al (2015). National Trends in the Performance of Robot-Assisted Sacrocolpopexy. *J Endourology* Jul;29(7):777-83. doi: 10.1089/end.2014.0710. Epub 2015 Mar 10.
23. Paraiso, et al (2011). Laparoscopic Compared with Robotic Sacrocolpopexy for Vaginal Prolapse. A Randomized Controlled Trial. *Obstetrics & Gynecology*. 118(5)1005-1013.
24. Collinet P, Leguevaque P, Neme RM, Cela V, Barton-Smith P, Hébert T, Hanssens S, Nishi H, Nisolle M. "Robot-assisted laparoscopy for deep infiltrating endometriosis: international multicentric retrospective study." *Surgical Endoscopy* 28.8 (2014):2474-2479.
25. Nezhat, Camran, Anna M. Modest, and Louise P. King. "The Role of the Robot in Treating Urinary Tract Endometriosis." *Current Opinion in Obstetrics and Gynecology* 25.4 (2013): 308-11.
26. Siesto, Gabriele, Nicoletta Ieda, Riccardo Rosati, and Domenico Vitobello. "Robotic Surgery for Deep Endometriosis: A Paradigm Shift." *The International Journal of Medical Robotics and Computer Assisted Surgery* 10 (2013): 140-46.
27. Bedaiwy, Mohamed A., Mohamed Abdel Y. Rahman, Mark Chapman, Heidi Frasure, Sangeeta Mahajan, Vivian E. Von Gruenigen, William Hurd, and Kristine Zanotti. "Robotic-Assisted Hysterectomy for the Management of Severe Endometriosis: A Retrospective Review of Short-Term Surgical Outcomes." *JSLs, Journal of the Society of Laparoendoscopic Surgeons* 17.1 (2013): 95-99.
28. Ercoli, A., M. D'asta, A. Fagotti, F. Fanfani, F. Romano, G. Baldazzi, M. G. Salerno, and G. Scambia. "Robotic Treatment of Colorectal Endometriosis: Technique, Feasibility and Short-term Results." *Human Reproduction* 27.3 (2012): 722-26.
29. Dulemba, John F., Cyndi Pelzel, and Helen B. Hubert. "Retrospective Analysis of Robot-assisted versus Standard Laparoscopy in the Treatment of Pelvic Pain Indicative of Endometriosis." *Journal of Robotic Surgery* 7.2 (2013): 163-69.
30. Nezhat, C. L., M.; Kotikela, S.; Veeraswamy, A.; Saadat, L.; Hajhosseini, B. (2010). "Robotic versus standard laparoscopy for the treatment of endometriosis." *Fertility and Sterility*.
31. Nezhat, CR; Stevens, A; Balassiano, E; and Rose Soliemannjad. "Robotic-Assisted Laparoscopy vs Conventional Laparoscopy for the Treatment of Advanced Stage Endometriosis." *JMIG* 22.1 (2014): 40-44.
32. O'Neill, M., et al. (2013). Robot-assisted hysterectomy compared to open and laparoscopic approaches: systematic review and meta-analysis. *Archives of Gynecology and Obstetrics*. 287: 907-918.
33. O'Sullivan, S. (2011). HIQA Ireland Health technology assessment of robot-assisted surgery in selected surgical procedures.
34. Ran, L., et al. (2014). Comparison of robotic surgery with laparoscopy and laparotomy for treatment of endometrial cancer: a meta-analysis. *PLoS ONE*. 9: e108361.
35. Reza, M., et al. (2010). Meta-analysis of observational studies on the safety and effectiveness of robotic gynaecological surgery. *British Journal of Surgery*.
36. Gaia, G., et al. (2010). Robotic-assisted hysterectomy for endometrial cancer compared with traditional laparoscopic and laparotomy approaches: a systematic review. *Obstetrics and Gynecology*. 116: 1422-1431.

My Results: [Procedure]

(optional slide)

Mean	Open Surgery	Traditional Laparoscopy	<i>da Vinci</i> [®] Surgery	<i>p</i> -Value
N	xx	yy	zz	p
BMI	xx	yy	zz	p
EBL or Transfusion (mL)	xx	yy	zz	p
Complications (%)	xx	yy	zz	p
Conversion Rate	xx	yy	zz	p
Length of Hospital Stay (days)	xx	yy	zz	p
Other 1	xx	yy	zz	p
Other 2	xx	yy	zz	p

[Please share your own data]