

INTUITIVE



Intuitive da Vinci SP Instrument and Accessory Catalogue Europe

February 2024

For more information, call Intuitive.

+800 0821 2020 (toll-free) / +41 21 821 2020
or visit www.intuitive.com

For the name of your local distributor,
call Intuitive or visit www.intuitive.com.

Table of contents

6 mm instruments

- 4 Monopolar cautery instruments
- 4 Bipolar cautery instruments
- 5 Instruments

Accessories

- 6 Reusable accessories
- 6 Disposable accessories

Vision equipment

- 6 Endoscope

Training equipment

- 7 Training instruments

Cleaning and sterilization

- 7 Cleaning and sterilization

Documentation

- 7 Documentation

Specifications




- 8 Instrument specifications
- 9 Common carton dimensions and weights

Contact information



- 10 Technical support and customer service

6 mm instruments

Monopolar cautery instruments

	Product description	Part number
	Monopolar curved scissors 25 uses	430004
	Shown with monopolar curved scissors tip (sold separately - 1 procedure) ¹	430035
	Monopolar cautery instrument 25 uses	430007
	Shown with cautery hook tip (sold separately - 1 procedure) ¹	400156
	Monopolar cautery instrument 25 uses	430007
	Shown with cautery spatula tip (sold separately - 1 procedure) ¹	400160





Bipolar cautery instruments

	Product description	Part number
	Maryland bipolar forceps 25 uses	430010
	Fenestrated bipolar forceps 25 uses	430011

¹ The cautery spatula tip, cautery hook tip, and monopolar curved scissors tip can be ordered separately. Cautery hook tip (box of 10) part number 400156. Cautery spatula tip (box of 10) part number 400160. Monopolar curved scissors tip (box of 10) part number 430035.

6 mm
instruments

Instruments

	Product description	Part number
	Needle driver 25 uses	430006
	Cadiere forceps 25 uses	430009
	Round tooth retractor 25 uses	430002
	Medium-large clip applier² 150 uses <small>(Shown with Weck® Hem-o-lok® polymer ligating clips medium-large.)</small> Validated for use only with Weck® Hem-o-lok® polymer ligating clips medium-large.	430005
	Weck product code:	544230









² Weck® Hem-o-lok® polymer ligating clips medium-large (Weck product code 544230) are the only compatible clips used with this instrument. To order, contact Teleflex Medical (www.teleflex.com) or your local distributor.

Accessories

Reusable accessories

	Product description	Qty/ box	Uses	Part number
	Cannula, 25 mm x 100 mm	1	NA	430036
	Obturator, 25 mm x 100 mm	1	NA	430016
	Bipolar cautery cord	1	20	470384
	EnergyShield monopolar cautery cord	1	NA	430068

Disposable accessories

	Product description	Qty/ box	Uses	Part number
	Monopolar curved scissors tip (box of 10)	10	1 procedure	430035
	Cautery hook tip (box of 10)	10	1 procedure	400156
	Cautery spatula tip (box of 10)	10	1 procedure	400160
	Instrument sheath (box of 30)	30	1 procedure	430012
	Camera sheath (box of 10)	10	1 procedure	430020
	Instrument arm drape	1	1 procedure	430013
	Access Port kit, small incision (2.7 cm to 4 cm incision size (box of 6) includes: • Access port • Short entry guide • Wound protector	6	1 procedure	430073
	Access Port kit, large incision (2.7 cm to 7 cm incision size (box of 6) includes: • Access port • Short entry guide • Wound protector	6	1 procedure	430075
	Entry guide kit (box of 6) includes: • Entry guide cannula insert • Entry guide cannula seal	6	1 procedure	430022

Vision equipment

Vision instruments

	Product description	Qty/ box	Uses	Part number
	SP Endoscope	1	33	430077





Training equipment

Training instruments

Product description	Qty/ box	Uses	Part number
Round tooth retractor	1	64	430002-T
Monopolar curved scissors	1	64	430004-T
Needle driver	1	64	430006-T
Monopolar cautery	1	64	430007-T
Maryland bipolar forceps	1	64	430010-T
Fenestrated bipolar forceps	1	64	430011-T

Cleaning and sterilization

Cleaning and sterilization

	Product description	Qty/ box	Uses	Part number
	SP instrument reprocessing tray	1	NA	430850
	SP endoscope reprocessing tray	1	NA	430851
	SP reprocessing connector (1 device)	1	NA	430852
	SP reprocessing connector (2 devices)	1	NA	430853
	Cleaning and sterilization kit includes: <ul style="list-style-type: none"> • Flush port adaptor • Brush with nylon bristles • Pipe-to-hose adapter • Leak test adapter 	1	NA	430048

Documentation


Documentation

Product description	Part number
Reprocessing documentation kit <ul style="list-style-type: none"> • Reprocessing instructions for instruments, accessories and endoscope 	400281
System user manual	554432
Instrument and accessory user manual	554438
Wall chart kit	381476

Documentation is listed in English only. For other languages contact Intuitive customer service. Documentation may not be available in all languages.

Instrument specifications

Part number	Uses	Product description	Jaw open angle	Jaw length	Maximum elbow articulation	Maximum working length	Maximum instrument length
430002	25	Round tooth retractor	40°	2.0 cm (0.80")	4.5 cm	27 cm	64.5 cm (25.4")
430004	25	Monopolar curved scissors	40°	1.3 cm (0.50")	3.0 cm	27 cm	64.5 cm (25.4")
430005	150	Medium-large clip applier	NA	1.5 cm (0.57")	3.0 cm	27 cm	64.5 cm (25.4")
430006	25	Needle driver	35°	1.0 cm (0.39")	3.0 cm	27 cm	64.5 cm (25.4")
430007	25	Monopolar cautery instrument	NA	NA	3.0 cm	27 cm	64.5 cm (25.4")
430009	25	Cadiere forceps	40°	1.8 cm (0.70")	3.0 cm	27 cm	64.5 cm (25.4")
430010	25	Maryland bipolar forceps	40°	1.8 cm (0.71")	3.0 cm	27 cm	64.5 cm (25.4")
430011	25	Fenestrated bipolar forceps	40°	2.0 cm (0.80")	3.0 cm	27 cm	64.5 cm (25.4")



Common carton dimensions and weights

Item	Item code	Outside dimensions	Average weight
Instruments	Most	73.66 x 6.83 x 18.10 cm 29" x 2.69" x 7.13"	0.54 kg 1.18 lbs
Monopolar curved scissors tip	430035	28.8925 x 11.43 x 8.59 cm 11.375" x 4.50" x 3.38"	0.34 kg 0.75 lbs
Access Port kit (small and large incision)	430073 430075	24.1 x 11.7 x 43.9 cm 9.5" x 4.6" x 17.3"	0.9 kg 1.9 lbs
Entry guide kit	430022	25.55 x 13.18 x 7.14 cm 10.063" x 5.188" x 2.81"	2.05 kg 4.52 lbs
Sheaths	430020	65.10 x 8.26 x 13.34 cm 25.63" x 3.25" x 5.25"	0.61 kg 1.34 lbs
	430012		0.72 kg 1.58 lbs
Cannula, 25 mm x 100 mm	430036	36.35 x 13.97 x 12.34 cm 14.31" x 5.50" x 4.86"	0.48 kg 1.06 lbs
Obturator, 25 mm x 100 mm	430016	36.35 x 13.97 x 12.34 cm 14.31" x 5.50" x 4.86"	0.47 kg 1.04 lbs
Shielded monopolar cord	430068	24.77 x 16.19 x 3.49 cm 9.75" x 6.38" x 1.38"	0.28 kg 0.62 lbs
Bipolar cautery cord	470384	24.77 x 16.19 x 3.49 cm 9.75" x 6.38" x 1.38"	0.26 kg 0.57 lbs
Drape	430013	39.70 x 38.58 x 18.64 cm 15.63" x 15.19" x 7.34"	2.95 kg 6.50 lbs
SP Endoscope	430077	81.6102 x 29.21 x 9.3726 cm 32.13" x 11.5" x 3.69"	2.27 kg 5.00 lbs

Technical Support and Customer Service

Contact us by phone

+800 0821 2020 (international toll-free) or use a local/local toll-free number:

Austria 0126 752 25 0800 298 890	Finland 0800 415 729 0800 416 321	Ireland 019 131 227 1800 851 011	Switzerland 021 821 2020 0800 821 200
Belgium 02 888 17 49 0800 716 96	France 05 37 10 02 26 0800 905 720	Norway 21 98 41 47	The Netherlands 020 8081 348 0800 0249 547
Denmark 89 872 907 808 31158	Germany 0761 488 8962 0800 1825 068	Sweden 08 580 987 81 020 127 661	United Kingdom 01865 521 820 0800 066 8902

The options during your call

Option 1: Emergency Technical Support

Support in local language:

Business days,
Mo. – Fr. 07:00 – 17:00 CET

Support in English: 365/24/7

Option 2: Non-urgent Technical Question

Business days,
Mo. – Fr. 07:00 – 17:00 CET

Option 3: Customer Service

- Purchase order inquiries
- Complaint
- Product return authorization
- Product exchanges

Business days,
Mo. – Fr. 08:00 – 18:00 CET

You can access services also 24/7
on the Intuitive Customer Portal.

For more information please visit:
my.intuitive.com

Option 4: Digital Products Support

- User account management
- Login and password issues
- Troubleshooting etc.

Business days,
Mo. – Fr. 08:00 – 18:00 CET

Contact us by e-mail

Austria:
Support.AT@intusurg.com

Belgium:
Support.BE@intusurg.com

Denmark:
Support.DK@intusurg.com

Finland:
Support.FI@intusurg.com

France:
Support.FR@intusurg.com

Germany:
Support.DE@intusurg.com

Ireland:
Support.IE@intusurg.com

Luxembourg:
Support.LU@intusurg.com

Norway:
Support.NO@intusurg.com

Sweden:
Support.SE@intusurg.com

Switzerland:
Support.CH@intusurg.com

The Netherlands:
Support.NL@intusurg.com

United Kingdom:
Support.UK@intusurg.com

Digital Products Support:
EU.DigitalSupport@intusurg.com

European Regional Offices

European Headquarters (Switzerland)

Intuitive Surgical Sàrl
1, chemin des Müriers
1170 Aubonne, Suisse
Tel: +41 21 821 2000
Fax: +41 21 821 2001
VAT #: CHE-113.407.162

Germany

Intuitive Surgical Deutschland GmbH
Am Flughafen 6
79108 Freiburg, DE
MWST #: DE200038727
Tel: +49 761 8878 7700
Amtsgericht Freiburg, HRB 720585
Geschäftsführer: Dirk Barten,
Jamie Samath, Jean-Yves Raimon-
Dacunha-Castelle, Gary Loeb

France

Intuitive Surgical SAS
Cité de la Photonique,
bâtiment Gienah
11 avenue de Canteranne
33600 Pessac, FR
Tel: +33 5 86 07 01 01
TVA #: FR 64 417 889 177
SIREN 417 889 177
Capital Social 37000€
RCS Bordeaux Code APE 4646Z

United Kingdom

Intuitive Surgical Ltd
The Schrodinger Building
Oxford Science Park
Heatley Rd, Oxford, OX4 4GE, UK
Tel: +44 (0) 1865 411561
VAT #: GB943 2213 49
Company registration
No 6727762

Global Headquarters

USA

Intuitive Surgical, Inc.
1020 Kifer Road, Building 108
Sunnyvale, CA 94086-5304
Tel: +1.408.523.2100
Fax: +1.408.523.1390

For more information, call Intuitive.

+800 0821 2020 (toll-free) / +41 21 821 2020
or visit www.intuitive.com

For the name of your local distributor,
call Intuitive or visit www.intuitive.com.

For intended use and/or indications for use, risks,
cautions, and warnings and full prescribing in-
formation, refer to the associated user manu-
al(s) or visit <https://manuals.intuitivesurgical.com/market>.

Intuitive's Privacy Notice is available at
www.intuitive.com/privacy.

© 2024 Intuitive Surgical Operations, Inc.
All rights reserved. Product and brand names/
logos are trademarks or registered trademarks
of Intuitive Surgical or their respective owner.