

EndoWrist Stapler 30

For the da Vinci Xi and
da Vinci X Surgical Systems

The Stapler 30 Curved-Tip provides enhanced placement and visualization around tenuous vasculature. Combined with the da Vinci Xi and da Vinci X Surgical System's crystal clear 3D HD vision, access, and surgeon control and autonomy, the EndoWrist Stapler 30 naturally expands the da Vinci stapling portfolio.

Surgeon Control & Autonomy

Positioning, gripping, clamping, and firing of the EndoWrist Stapler are completely controlled through the da Vinci Surgeon Console.

Fully Wristed Articulation

The EndoWrist Stapler has the greatest side-to-side articulation currently available, with a range of 108° total side-to-side. And unlike other staplers currently available, the EndoWrist Stapler articulates through a range of 54° total up-and-down, providing a fully-wristed cone of articulation.

Smartclamp Feedback



SmartClamp feedback reduces the guesswork associated with conventional stapling. Prior to firing, SmartClamp detects if the Stapler jaws are adequately closed on the tissue for the installed reload color.



Instruments

Stapler 30 Instruments





Stapler 30 shares the same cannula accessories and sheath as the da Vinci Xi and da Vinci X Stapler 45.

		Part Number
	Stapler 30 instrument 1 unit/50 firings	470430
	Stapler 30 Curved-Tip instrument 1 unit/50 firings	470530

Reloads

Stapler 30 Reloads

All reloads include a new knife.

	Product description	Qty/box	Uses	Item code
	Gray reload 6-row of staples	12	1 firing	48630M
	White reload 6-row of staples	12	1 firing	48630W
	Blue reload 6-row of staples	12	1 firing	48630B
	Green reload 6-row of staples	12	1 firing	48630G

To contact a representative or receive additional information:
Visit www.intuitive.com or call Intuitive Customer Service:
+800.0821.2020 (toll free) or +41.21.821.2020

Disclaimer

The EndoWrist Stapler 30 Instrument and Stapler 30 Reloads are intended to be used with the da Vinci Xi and da Vinci X Surgical Systems for resection, transection and/or creation of anastomoses in General, Thoracic, Gynecologic, and Urologic surgery. The device is indicated for adult and pediatric use. The device can be used with staple line or tissue buttressing material.

The EndoWrist Stapler 30 Instrument and Stapler 30 Reloads are class IIa and IIb medical devices CE marked (CE 2460) under the European Medical Devices Directive (93/42/EEC), manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.

© 2021 Intuitive Surgical, Inc. All rights reserved. Product names are trademarks or registered trademarks of their respective holders.